



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Flamingo Curo

Company:	Ship To:
Account Number:	Date:
Quote: PO Number:	Address:
Clinician / Sales Rep:	City: Province:
Marked For:	Postal Code: Phone:
Associated Quote/Order#:	Email/Fax:
	Effective: 4/13/2026

Base Frame and Mobility Base Options

Flamingo Curo Seat and Back



Seat and Back, Size 1 (882001)	\$945
Seat and Back, Size 2 (882002)	\$1015
Seat and Back, Size 3 (882003)	\$1113

Seat and Back Upholstery

Upholstery must match seat size selection



Seat Base Upholstery Inlay

Seat Base Upholstery Inlay, Blue, Size 1 (882071)	\$209
Seat Base Upholstery Inlay, Blue, Size 2 (882072)	\$227
Seat Base Upholstery Inlay, Blue, Size 3 (882073)	\$245



Back Upholstery Inlay

Back Upholstery Inlay, Blue, Size 1 (882075)	\$215
Back Upholstery Inlay, Blue, Size 2 (882076)	\$224
Back Upholstery Inlay, Blue, Size 3 (882077)	\$245

Flamingo Curo Base Frame

Frames must match seat size selection

Flamingo Curo Compact Frame



Flamingo Curo Compact Frame with 3" Casters, Size 1 (882011)	\$1393
Flamingo Curo Compact Frame with 3" Casters, Size 2 (882012)	\$1442
Flamingo Curo Compact Frame with 3" Casters, Size 3 (882013)	\$1505



Potty Frame

Potty Frame, Size 1 (882050)	\$455
------------------------------	-------

Available only with Size 1 Seat

Flamingo Curo High-Low Frame




Flamingo Curo High-Low Frame with 3" Casters, Size 1 (882031)	\$4410
Flamingo Curo High-Low Frame with 3" Casters, Size 2 (882032)	\$4550
Flamingo Curo High-Low Frame with 3" Casters, Size 3 (882033)	\$4690



Hardware for Toilet Mounting

Direct Toilet Mounting, Size 1 (882061)	\$378
Direct Toilet Mounting, Size 2 (882062)	\$399
Direct Toilet Mounting, Size 3 (882063)	\$413



Suction Cups

Suction Cups (880150)	\$245
-----------------------	-------

Requires Commode Pan / Potty Ring Holders



Casters

4" Casters (880100-S)	\$499
-----------------------	-------

Footrest

Footrest options must match seat size selection



Foot Support with Quick Release

Removeable Foot Support, Size 1 (882151)	\$507
Removeable Foot Support, Size 2 (882152)	\$538
Removeable Foot Support, Size 3 (882153)	\$567



Calf Support

Calf Support, Size 1 (882135)	\$213
Calf Support, Size 2 (882136)	\$217
Calf Support, Size 3 (882137)	\$223

Pelvic and Trunk Positioning

Pelvic Positioning



Seat Infill

Must match seat size selection

Seat Infill for Seat, Blue PU, Size 1 (882121)	\$210
Seat Infill for Seat, Blue PU, Size 2 (882122)	\$210
Seat Infill for Seat, Blue PU, Size 3 (882123)	\$210



Splashguard

Splashguard, Small (882145)	\$182
Splashguard, Large (882147)	\$182



Hip Guide

Requires Armrests or Hip Guide Attachment Rod

Hip Guide, Small (882111)	\$245
Hip Guide, Large (882112)	\$252



Hip Belt

Neoprene Hip Belt, Size 1 (85423-1)	\$150
Neoprene Hip Belt, Size 2 (85423-2)	\$150
Neoprene Hip Belt, Size 3 (85423-3)	\$160



Hip Guide Attachment Rod (882116)	\$42
-----------------------------------	------

Not compatible with Armrests

Trunk Positioning



Shoulder Support with PU

Not compatible with Side Supports. Must match seat size selection.

Shoulder Support Inlay, Blue PU, Size 1 (882081)	\$315
Shoulder Support Inlay, Blue PU, Size 2 (882082)	\$340
Shoulder Support Inlay, Blue PU, Size 3 (882083)	\$371



Lower Back Support

Must match seat size selection

Lower Back Support, Size 1 (882131)	\$231
Lower Back Support, Size 2 (882132)	\$231
Lower Back Support, Size 3 (882133)	\$231



Fixed Side Supports

Fixed Side Supports, Small (882090)	\$252
Fixed Side Supports, Medium (882091)	\$273
Fixed Side Supports, Large (882092)	\$301



Swing-Away Side Supports

Swing-Away Side Supports, Small (882095)	\$455
Swing-Away Side Supports, Medium (882096)	\$483
Swing-Away Side Supports, Large (882097)	\$511



Cross Vest

Cross Vest, Net, Size 0 (89625-040)	\$300
Cross Vest, Net, Size 1 (89625-140)	\$300
Cross Vest, Net, Size 2 (89625-240)	\$300
Cross Vest, Net, Size 3 (89625-340)	\$300
Cross Vest, Net, Size 4 (89625-440)	\$300



Vest

Vest, Net, Size 0 (89111-040)	\$307
Vest, Net, Size 3/4 (89111-07540)	\$307
Vest, Net, Size 1 (89111-140)	\$307
Vest, Net, Size 2 (89111-240)	\$307
Vest, Net, Size 3 (89111-340)	\$307
Vest, Net, Size 4 (89111-440)	\$307

Headrests and Hardware



Head Support, Small, Curved, Standard Hardware (9931056-61)	\$651
Head Support, Medium, Curved, Standard Hardware (9932056-61)	\$679
Head Support, Small, Curved, Extended Hardware (9931066-61)	\$665
Head Support, Medium, Curved, Extended Hardware (9932066-61)	\$693

Armrest Options



Armrests

Armrests, Small (882113)	\$329
Armrests, Large (882114)	\$357



Tray

Requires Armrests to mount. Must match seat size selection.

Tray, Blue, Size 1 (882141)	\$350
Tray, Blue, Size 2 (882142)	\$385
Tray, Blue, Size 3 (882143)	\$413

Accessories

Handles and Push Braces



Push Brace

Push Brace (880110)	\$277
---------------------	-------



Push Handle

Requires High-Low Frame

Helper Handle, Size 1 (880110-HL)	\$110
Helper Handle, Sizes 2 and 3 (880110-HLL)	\$110

Potty Options



Commode Pan

Requires Commode Pan / Potty Ring Holders

Commode Pan, Size 1 (870800)	\$28
Commode Pan, Sizes 2 and 3 (870801)	\$63



Potty Ring

Requires Commode Pan / Potty Ring Holders

Potty Ring, Size 1 (870805)	\$63
Potty Ring, Sizes 2 and 3 (870807)	\$77



Commode Pan Lid

Commode Pan Lid, Size 1 (870800-L)	\$118
Commode Pan Lid, Sizes 2 and 3 (870802-L)	\$118



Extended Potty Ring

Requires Commode Pan / Potty Ring Holders




Potty Ring, Extended, Sizes 2 and 3 (870807-E)	\$77
--	------

















Commode Pan / Potty Ring Holders









Commode Pan / Potty Ring Holders, Size 1 (880125-1)	\$42
Commode Pan / Potty Ring Holders, Size 2 (880125-2)	\$49
Commode Pan / Potty Ring Holders, Size 3 (880125-3)	\$56

Measurement Guide

Seat size 1, 2 & 3			
mm (inch) kg (lb)	 1	 2	 3
 30° Backrest recline angle	(-20°)-(30°)	(-20°)-(30°)	(-20°)-(30°)
 Overall backrest height**	425-495 mm (16¾-12½")	470-550 mm (18½-21¾")	515-615 mm (20¼-24¼")
 Effective backrest height**	375/410/445 mm (14¾/16¼/17½")	405/445/485 mm (16/17½/19")	450/500/550 mm (17¾/19¾/21¾")
 Effective back width (with side supports)	140-240 mm (5½-9½")	170-270 mm (6¾-10¾")	220-320 mm (8¾-12½")
 Effective back width (without side supports)	300 mm (11¾")	335 mm (13¼")	390 mm (15¼")
 Effective seat depth**	240-300 mm (9½-11¾")	285-345 mm (11¼-13½")	335-425 mm (13¼-16¾")
 Overall length	360-375 mm (14¼-14¾")	395-430 mm (15½-17")	435-495 mm (17¼-19½")
 Effective seat width (without upper leg supports)	265 mm (10½")	300 mm (11¾")	355 mm (14")
 Effective seat width (with upper leg supports)	185/225/265 mm (7¼/8¾/10½")	220/260/300 mm (8¾/10¼/11¾")	275/315/355 mm (10¾/12½/14")
 Effective foot support to seat distance	220-285 mm (8¾-11¼")	260-345 mm (10¼-13½")	300-400 mm (11¾-15¾")
 Mass device*	5.2 kg (11.5 lb)	6.1 kg (13.4 lb)	7.5 kg (16.5 lb)
 Maximum user mass	<35 kg (<77 lb)	<55 kg (<121 lb)	<70 kg (<154 lb)
 Max rated load	42.2 kg (93 lb)	63.1 kg (139.1 lb)	79.5 kg (175.3 lb)

*seat only (without accessories) **including PUR

Seat & Compact frame				Seat & High-low frame			
mm (inch) kg (lb)	 1	 2	 3	 1	 2	 3	
 Overall height*	905-1000 mm (35¾-39¾")	950-1055 mm (37½-41½")	995-1120 mm (39¾-44")	895-1075 mm (35¼-42¼")	930-1205 mm (36½-47½")	975-1270 mm (38½-50")	
 Overall height, to pushbrace*	850-1020 mm (33½-40¾")	890-1070 mm (35-42¼")	930-1130 mm (36½-44½")	900-1180 mm (35½-46½")	945-1235 mm (37¼-48½")	990-1350 mm (39-53¼")	
 Overall length with leg rest*	935 mm (36¾")	955 mm (37½")	1110 mm (43¾")	930 mm (36½")	1065 mm (42")	1200 mm (47¼")	
 Overall width*	470 mm (18½")	510 mm (20")	560 mm (22")	485 mm (19")	580 mm (22¾")	580 mm (22¾")	
 Mass device*	12.4 kg (27.3 lb)	13.4 kg (29.5 lb)	14.9 kg (32.8 lb)	14.8 kg (32.6 lb)	17 kg (37.5 lb)	18.1 kg (39.9 lb)	
 Overall minimum turning diameter*	850 mm (33½")	875 mm (34½")	905 mm (35¾")	785 mm (31")	955 mm (37½")	955 mm (37½")	
 Seat surface height at front edge ***	480-580 mm (19-22¾")	480-580 mm (19-22¾")	480-580 mm (19-22¾")	530-740 mm (20¾-29¾")	530-740 mm (20¾-29¾")	530-740 mm (20¾-29¾")	
 Seat plane angle (tilt) **	(-10°)-(30°)	(-10°)-(30°)	(-10°)-(30°)	(-5°)-(20°)	(-5°)-(20°)	(-5°)-(20°)	

Seat & Potty frame							
	mm (inch) kg (lb)						
	 Overall height	 Overall length	 Overall length with leg rest***	 Overall width	 Mass device	 Effective foot support to seat distance***	 Seat surface height at front edge***
 1	630-750 mm (24¾-29½")	555 mm (21¾")	870 mm (34¾")	430 mm (17")	7.1 kg (16 lb)	220-285 mm (8¾-11¼")	290-330 mm (11½-13")

*seat and High-low frame with 75 mm wheels **(+/-5°) ***including PUR