












Company:	Ship To:		
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	Province:
Clinician / Sales Rep:	Postal Code:	Phone:	
Marked For:	Email/Fax:		
Associated Quote/Order#:		Effective: 1/10/2024	

Crocodile

	0	1	2	3
 Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
 Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼"*	21¾ - 33"	29½ - 41¼"
 Max. user height	41¼"	45¼"	59"	70¾"
 Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
 Length	25¼"	28¼"	33"	34¼"
 Overall width	22¾"	24¾"	27¼"	27¼"
 Folded length	25½"	28½"	33¾"	39¼"
 Folded width	22¾"	24½"	27¼"	27¼"
 Folded height	11¾"	12¾"	10½"	11½"
 Minimum turning radius	19"	19¾"	22"	45"
 Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

*Mounted with the Starter Kit accessories. Further measurements at our website.

- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



Crocodile Base

Crocodile, White, Size 0	\$895
Crocodile, Yellow, Size 1	\$987
Crocodile, Orange, Size 2	\$987
Crocodile, Metallic Black, Size 2	\$987
Crocodile, Anthracite Grey, Size 3	\$1628

Mounting Kits



Crocodile Starter Kit	\$651
<i>For size 1</i>	

Anti-Tippers (Pair)



Anti-Tip Device, Swing-away

Anti-tip Swing-away Pair for Sizes 0-2	\$214
Anti-tip Swing-away Pair for Size 3	\$214



Anti-Tip Device, Fixed

Anti-tip Fixed Pair for Sizes 0-2	\$102
Anti-tip Fixed Pair for Size 3	\$102

Arm Supports



Forearm Support Pair for Sizes 0-1

Forearm Support Pair for Sizes 0-1	\$576
Forearm Support Pair for Size 2	\$576
Forearm Support Pair for Size 3	\$576



Padded Forearm Positioner Strap

For use with Forearm Support with Handgrip. For quantity 2 select both options below.

Padded Forearm Positioner Strap	\$91
Padded Forearm Positioner Strap	\$91



Handle Knobs

Handle Knobs, Pair	\$68
--------------------	------

Brake

Hand brake eliminates the function of anti-reverse wheel stop. If ordering hand brake with anti-reverse wheel stop function, a custom order required



Hand Brake

Hand Brake Kit for Sizes 1-2	\$320
Hand Brake Kit for Size 3	\$320



Deceleration Hand Brake Kit

Deceleration Hand Brake Kit for Sizes 1-2	\$320
Deceleration Hand Brake Kit for Size 3	\$320

Holder for Hand Brake

Required with use of Hand brake and Forearm support. Allows for mounting of handbrake to Forearm Support.

Small for Size 1	\$487
Medium for Size 2	\$549
Medium for Size 2	\$590

Back and Pelvic Supports



Surrounding Supports

Surrounding Support _____ \$654
Available only with size 0



Back Supports

Back Support - Adjustable Padded _____ \$674
For Size 3



Pelvic Supports

Pelvic Support With Laterals for Sizes 0-2 _____ \$350
Pelvic Support With Laterals for Sizes 1-2 _____ \$350

Abdominal Support Belt

Abdominal Support Belt _____ \$164
Used with Back Support above, for size 3



Lateral Strap

Lateral Strap _____ \$68
For sizes 1 and 2



Cross Bar

Cross Bar for Mounting Back/Pelvic Supports for Size 3 _____ \$129
Required with all Back and Pelvic Supports for size 3

Handles and Push Braces



Attendant Handle

Attendant Handle _____ \$185
Not available with size 3



Grip Bar For Crocodile

Grip Bar for Size 0 _____ \$173
Grip Bar for Size 1 _____ \$198

Hip Support Pair



Hip Supports, Pair _____ \$313
Not available with size 3

Bags and Baskets



Storage Bag for Sizes 0-2 _____ \$95

Storage Bag, for Size 3 _____ \$280

Seating



Lightweight Seat Sling

Requires Pelvic Support with Laterals and Lateral Strap or Surrounding Support

Lightweight Sling Seat for Size 0-1 _____ \$179
Lightweight Sling Seat for Size 2 _____ \$179



Seat Sling

Seat Sling for Size 1 _____ \$179
Seat Sling for Size 2 _____ \$179

Flip Down Seat

White, Size 0 _____ \$226
Yellow, Size 1 _____ \$226
Not compatible with Crocodile Starter Kit
Orange, Size 2 _____ \$226
Metallic Black, Size 2 _____ \$226
Anthracite Grey, Size 3 _____ \$226