

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:












Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025

## Crocodile

		0	1	2	3
	Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
	Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼" *	21¾ - 33"	29½ - 41¼"
	Max. user height	41¼"	45¼"	59"	70¾"
	Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
	Length	25¼"	28¼"	33"	34¼"
	Overall width	22¾"	24¾"	27¼"	27¼"
	Folded length	25½"	28½"	33¾"	39¼"
	Folded width	22¾"	24½"	27¼"	27¼"
	Folded height	11¾"	12¼"	10½"	11½"
	Minimum turning radius	19"	19¾"	22"	45"
	Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

\*Mounted with the Starter Kit accessories. Further measurements at our website.

- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



## Crocodile Frame

### Crocodile Base

#### Frame

Crocodile, White, Size 0 (86800)	\$913
Crocodile, Yellow, Size 1 (86801)	\$1007
Crocodile, Orange, Size 2 (86802)	\$1007
Crocodile, Metallic Black, Size 2 (86852)	\$1007
Crocodile, Anthracite Grey, Size 3 (86803)	\$1661

### Mounting Kits



Crocodile Starter Kit, for Crocodile size 1 (86801-0)	\$665
<i>Used for Size 1</i>	

### Anti-Tippers (Pair)



#### Anti-Tip Device, Swing-away

Anti-tip Swing-away, Pair, for Size 0-2 (86815)	\$219
Anti-tip Swing-away, Pair, for Size 3 (86815-3S)	\$219



#### Anti-Tip Device, Fixed

Anti-tip Fixed, Pair, for Size 1-2 (86815-F)	\$105
Anti-tip Fixed, Pair, for Size 3 (86815-3F)	\$105

## Arm and Forearm Options

### Arm Supports



#### Forearm Support with Handgrip

Forearm Support Pair - for Sizes 0-1 (86858-1)	\$588
Forearm Support Pair - for Size 2 (86859-1)	\$588
Forearm Support Pair - for Size 3 (86860-1)	\$588



#### Padded Forearm Positioner (Each)

*For use with Forearm Support with Handgrip. For qty 2 selection both options below.*

Padded Forearm Positioner Strap (869370)	\$93
Padded Forearm Positioner Strap (869370)	\$93



#### Handle Knobs

Handle Knobs, Pair (86861)	\$70
----------------------------	------

## Handles and Grip Bars

### Grip Bar



#### Grip Bar For Crocodile

Grip Bar, for Size 0 (86879-0)	\$177
Grip Bar, for Size 1 (86879-1)	\$202
Grip Bar, for Size 2 (86879-2)	\$255

### Handles and Push Braces



#### Attendant Guide Bar

Attendant Handle (86866-1)	\$289
----------------------------	-------

#### Angle Adjustable Grip Bar

Angle Adjustable Grip Bar, Size 0 (86878-0)	\$329
Angle Adjustable Grip Bar, Size 1 (86878-1)	\$385
Angle Adjustable Grip Bar, Size 2 (86878-2)	\$427

## Brakes

### Brake

*Hand brake eliminates the function of anti-reverse wheel stop. If ordering hand brake with anti-reverse wheel stop function, a custom order required.*



#### Hand Brake

Hand Brake Kit, for Size 1-2 (86834) \$327  
Hand Brake Kit, for Size 3 (86873) \$327



#### Deceleration Hand Brake Kit

Deceleration Hand Brake Kit, Sizes 1-2 (86834-DEC) \$327  
Deceleration Hand Brake Kit, Size 3 (86873-DEC) \$327

### Holder for Hand Brake

*Required with use of Hand brake and Forearm support. Allows for mounting of handbrake to Forearm support.*

#### Holder for Hand Brake

Small, for Size 1 (86874-1) \$497  
Medium, for Size 2 (86874-2) \$560  
Large, for Size 3 (86874-3) \$602

## Back and Pelvic and Hip Supports

### Back and Pelvic Supports



#### Surrounding Supports

Surrounding Support for Size 0 (868610) \$668



#### Back Support

Back support - Adjustable Padded for Size 3 (86886) \$688



#### Pelvic Support With Laterals

Pelvic Support With Laterals, for Sizes 0-2 (86808) \$357  
Pelvic Support With Laterals, For Sizes 1-2 (86809) \$357

#### Abdominal Support Belt

*For size 3 back support only*

Abdominal Support Belt, for Size 3 (85490-2) \$168



#### Lateral Strap

Lateral Strap (905011) \$70



#### Mounting Hardware

*Required for use of pelvic and back support on size 3*

Cross Bar for Mounting Back/Pelvic Supports, for Size 3 (86880) \$132

### Hip Support Pair



#### Hip Support Pair

Hip Supports, Pair (86819) \$320

## Seating Options

Solid Flip Down Seat


	<b>Solid Flip Down Seat</b>	
	White, Size 0 (86820-01)	\$231
	Yellow, Size 1 (86820)	\$231
	<i>Not compatible with Crocodile starter kit</i>	
	Orange, Size 2 (86821)	\$231
	Metallic Black, Size 2 (86823)	\$231
	Anthracite Grey, Size 3 (86863)	\$231

Seat Slings

	Seat Sling, for Size 1 (86824)	\$183		<i>Requires pelvic support with lateral and lateral strap or surrounding support</i>		
	Seat Sling, for Size 2 (86825)	\$183				
					Lightweight Sling Seat, for Size 0-1 (86824)	\$183
					Lightweight Sling Seat, for Size 2 (86825-1)	\$183

## Accessories

Bags and Baskets

	<b>Bags and Baskets</b>	
	Storage Bag for Sizes 0-2 (86839)	\$97
	Storage Bag, for Size 3 (86843)	\$286