



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

R82 Crocodile

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

[Submit Form](#)

Effective: 1/10/2024

Crocodile

		0	1	2	3
	Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
	Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼"*	21¾ - 33"	29½ - 41¼"
	Max. user height	41¼"	45¼"	59"	70¾"
	Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
	Length	25¼"	28¼"	33"	34¼"
	Overall width	22¾"	24¾"	27¼"	27¼"
	Folded length	25½"	28½"	33¾"	39¼"
	Folded width	22¾"	24½"	27¼"	27¼"
	Folded height	11¾"	12¾"	10½"	11½"
	Minimum turning radius	19"	19¾"	22"	45"
	Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

*Mounted with the Starter Kit accessories. Further measurements at our website.

- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1197.1 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

R82 Crocodile

Crocodile Base

Crocodile, White, Size 0	\$718
Crocodile, Yellow, Size 1	\$792
Crocodile, Orange, Size 2	\$792
Crocodile, Metallic Black, Size 2	\$792
Crocodile, Anthracite Grey, Size 3	\$1306

Mounting Kits



Crocodile Starter Kit	\$502
<i>For size 1</i>	

Anti-Tippers (Pair)



Anti-Tip Device, Swing-away

Anti-tip Swing-away Pair for Sizes 0-2	\$171
Anti-tip Swing-away Pair for Size 3	\$171



Anti-Tip Device, Fixed

Anti-tip Fixed Pair for Sizes 0-2	\$82
Anti-tip Fixed Pair for Size 3	\$82

Arm Supports



Forearm Support Pair for Sizes 0-1

Forearm Support Pair for Sizes 0-1	\$462
Forearm Support Pair for Size 2	\$462
Forearm Support Pair for Size 3	\$462



Padded Forearm Positioner Strap

For use with Forearm Support with Handgrip. For quantity 2 select both options below.

Padded Forearm Positioner Strap	\$73
Padded Forearm Positioner Strap	\$73



Handle Knobs

Handle Knobs, Pair	\$54
--------------------	------

Brake

Hand brake eliminates the function of anti-reverse wheel stop. If ordering hand brake with anti-reverse wheel stop function, a custom order required



Hand Brake

Hand Brake Kit for Sizes 1-2	\$256
Hand Brake Kit for Size 3	\$256



Deceleration Hand Brake Kit

Deceleration Hand Brake Kit for Sizes 1-2	\$256
Deceleration Hand Brake Kit for Size 3	\$256

Holder for Hand Brake

Required with use of Hand brake and Forearm support. Allows for mounting of handbrake to Forearm Support.

Small for Size 1	\$359
Medium for Size 2	\$405
Medium for Size 2	\$436



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

R82 Crocodile

Back and Pelvic Supports



Surrounding Supports

Surrounding Support _____ \$524
Available only with size 0



Back Supports

Back Support - Adjustable Padded _____ \$541
For Size 3



Pelvic Supports

Pelvic Support With Laterals for Sizes 0-2 _____ \$281
Pelvic Support With Laterals for Sizes 1-2 _____ \$281

Abdominal Support Belt

Abdominal Support Belt _____ \$131
Used with Back Support above, for size 3



Lateral Strap

Lateral Strap _____ \$54
For sizes 1 and 2



Cross Bar

Cross Bar for Mounting Back/Pelvic Supports for Size 3 _____ \$104
Required with all Back and Pelvic Supports for size 3

Handles and Push Braces



Attendant Handle

Attendant Handle _____ \$149
Not available with size 3



Grip Bar For Crocodile

Grip Bar for Size 0 _____ \$138
Grip Bar for Size 1 _____ \$159

Hip Support Pair



Hip Supports, Pair _____ \$251
Not available with size 3

Bags and Baskets



Storage Bag for Sizes 0-2 _____ \$76

Storage Bag, for Size 3 _____ \$225

Seating



Lightweight Seat Sling

Requires Pelvic Support with Laterals and Lateral Strap or Surrounding Support

Lightweight Sling Seat, For Size 0-1 _____ \$144
Lightweight Sling Seat, For Size 2 _____ \$144



Seat Sling

Seat Sling for Size 1 _____ \$144
Seat Sling for Size 2 _____ \$144

Flip Down Seat

White, Size 0 _____ \$182
Yellow, Size 1 _____ \$182
Not compatible with Crocodile Starter Kit
Orange, Size 2 _____ \$182
Metallic Black, Size 2 _____ \$182
Anthracite Grey, Size 3 _____ \$182