











Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



# Caribou

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote: PO Number: \_\_\_\_\_ City: \_\_\_\_\_ Province: \_\_\_\_\_  
 Clinician / Sales Rep: \_\_\_\_\_ Postal Code: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_  
 Associated Quote/Order#: \_\_\_\_\_ Effective: 3/3/2025

## Caribou

mm (inch) kg (lb)				
	 1	 2	 3	 4
 Overall width	510 mm (20")	570 mm (22½")	650 mm (25½")	750 mm (29½")
 Length	750 mm (29½")	920 mm (36¼")	1110 mm (43¾")	1110 mm (43¾")
 Overall height, vertical	460–570 mm (18–22½")	530–640 mm (20¾–25¼")	680–790 mm (26¾–31")	880–990 mm (34¾–39")
 Overall height, horizontal	700 mm (27½")	700 mm (27½")	800 mm (31½")	800 mm (31½")
 Angling	0°–87°	0°–87°	0°–87°	0°–87°
 Foot plate angle	(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)
 Total mass	20.5 kg (45.1 lb)	24 kg (52.9 lb)	30.5 kg (67.2 lb)	35 kg (77.1 lb)
 Maximum rated load	62.5 kg (137.7 lb)	86 kg (189.6 lb)	112.5 kg (248 lb)	137 kg (302 lb)
 Total mass, with power	22.5 kg (49.6 lb)	25.5 kg (56.2 lb)	32.5 kg (71.7 lb)	37 kg (81.6 lb)
 Maximum rated load, with power	64.5 kg (142.1 lb)	87.5 kg (192.9 lb)	114.5 kg (252.4 lb)	139 kg (306.4 lb)
 Maximum user mass	40 kg (88.1 lb)	60 kg (132.2 lb)	80 kg (176.3 lb)	100 kg (220.4 lb)
 User height	850–1120 mm (33½–44")	980–1290 mm (38½–50¾")	1250–1560 mm (49¾–61½")	1530–1840 mm (60¾–72½")
 Static stability downhill	max. 5°	max. 5°	max. 5°	max. 5°
 Static stability uphill	max. 5°	max. 5°	max. 5°	max. 5°
 Sideways stability sideways	max. 5°	max. 5°	max. 5°	max. 5°

## Caribou Base

### Caribou Base

*Required selection*

Caribou Base, Size 1 (860101-FP)	\$1982
Caribou Base, Size 2 (860102-FP)	\$2045
Caribou Base, Size 3 (860113-FP)	\$2279
Caribou Base, Size 4 (860114-FP)	\$2688

### Caribou Base with Power Actuator

*Power Actuator recommended for Size 3 and 4.*

Caribou Base with Power Actuator, Size 3 (860103-ELFP)	\$4017
Caribou Base with Power Actuator, Size 4 (860104-ELFP)	\$4385



### Caribou Mast with Torso - Lower Body Pad

*Required Selection*

Lower Body Pad, Size 1 (X860101-P)	\$1508
Lower Body Pad, Size 2 (X860101-P)	\$1553
Lower Body Pad, Size 3 (X860113-P)	\$1734
Lower Body Pad, Size 4 (X860114-P)	\$2044


### Multi Adjustable Footrests

*Required Selection*

Multi Adjustable Footrest, Size 1 (X860101-F)	\$524
Multi Adjustable Footrest, Size 2 (X860102-F)	\$564
Multi Adjustable Footrest, Size 3 (X860113-F)	\$597
Multi Adjustable Footrest, Size 4 (X860114-F)	\$612

## Chest Pad


*Chest Pad is a required selection.*



Chest Pad, Size 1, A: 6", B: 7", C: 7.75" (861101-1)	N/C
Chest Pad, Size 2, A: 7.5", B: 9", C: 9.75" (861101-2)	N/C
Chest Pad, Size 3, A: 8.75", B: 10.5", C: 12.5" (861101-3)	N/C
Chest Pad, Size 4, A: 9.75", B: 12.25", C: 15.75" (861101-4)	N/C

## Extra Cushion


*Choose the same size as the Caribou Base.*



### Extra Cushion

Extra Cushion, Size 1, A: 9.75", B: 12.25", C: 15.75" (861101-4) (Each)	\$314
Extra Cushion, Size 2: W: 11.75" x H: 4" (861106-2) (Each)	\$314
Extra Cushion, Size 3, W: 13.75" x H: 4" (861106-3) (Each)	\$341
Extra Cushion, Size 4, W: 15.75" x H: 4" (861106-4) (Each)	\$341

## Supports, Straps and Wedges



### Wedges

Wedges for Side Slots, Size 1 (861107-1) (Pair)	\$23
Wedges for Side Slots, Size 2 (861107-2) (Pair)	\$23
Wedges for Side Slots, Size 3 (861107-3) (Pair)	\$23



### Back Supports


Back Supports, Size 1 and 2 (861115-2) (Each)	\$402
Back Supports, Size 3 and 4 (861115-3) (Each)	\$408



### Upper Arm Supports

*Can not be used with shoulder supports. Can not be used with lateral supports.*

Upper Arm Support, Small (861109-10)	\$533
Upper Arm Support, Large (861109-20)	\$544



### Shoulder Support (Pair)

*Pad size must match Chest Pad size.*

Shoulder Support, Size 1, A: 3", B: 7.5", C: 2" (861126-1) (Pair)	\$348
Shoulder Support, Size 2, A: 3.25", B: 9.75", C: 2" (861126-2) (Pair)	\$348
Shoulder Support, Size 3, A: 4", B: 12.5", C: 2" (861126-3) (Pair)	\$348
Shoulder Support, Size 4, A: 4.25", B: 15.25", C: 2" (861126-4) (Pair)	\$348

## Trunk Positioning- Vests

Measuring Shoulder to Pelvis Diagonally. See pictural under Belts and Straps.

### Cross Vest



Cross Vest, XS, S-P: 7.75" - 12.5" (85720-XS) (Each)	\$259
Cross Vest, Small, S-P: 10.75" - 12.5" (85720-S) (Each)	\$259
Cross Vest, Medium, S-P: 12.5" - 16" (85720-M) (Each)	\$259
Cross Vest, Large, S-P: 16" - 18.5" (85720-L) (Each)	\$259
Cross Vest, XL, S-P: 20.5" - 23.5" (85720-XL) (Each)	\$259




### Torso Vest (Ea)

Torso Vest, XS, S-P: 7.75"-10.75" (85705-XS)	\$273
Torso Vest, Small, S-P: 10.75"-12.5" (85705-S)	\$273
Torso Vest, Medium, S-P: 12.5"-16" (85705-M)	\$273
Torso Vest, Large, S-P: 16"-18.5" (85705-L)	\$280
Torso Vest, XL, S-P: 18.5"-20.5" (85705-XL)	\$273
Torso Vest, XXL, S-P: 20.5"-23.5" (85705-XXL)	\$273

## Vest Mounting Hardware

Required Selection for All Vests



### Fix Locks

Fix-Locks, Black, Bag of 6 Pcs (862109-35)	\$58
--	------

## Belts and Straps



### Chest Strap (Ea)

If multiple straps are required, please add to notes section below or contact customer service.




Chest Strap and Pelvic Strap, Small, W: 16.5" - 24.75" (862101-2)	\$164
Chest Strap and Pelvic Strap, Medium, W: 19" - 30.25" (862101-3)	\$164
Chest Strap and Pelvic Strap, Large, W: 22.75" - 37.75" (862101-4)	\$164



### Cross Straps

Cross Strap, Size 1, 15" - 19.25" (862102-1) (Pair)	\$381
Cross Strap, Size 2, 17" - 21.5" (862102-2) (Pair)	\$381
Cross Strap, Size 3, 20.75" - 27.5" (862102-3) (Pair)	\$402
Cross Strap, Size 4, 24.75" - 33.75" (Pair)	\$402

## Side Supports



### Side Supports, Fixed (Pair)

Supports, Fixed, Size 1 and 2, A: 4.75" x B: 4.25" (861102-10) (Pair)	\$341
Side Supports, Fixed, Type 5, Size 3 and 4, A: 6" x 6.25" (861102-20) (Pair)	\$408

### Adjustable Plate for Fixed Side Supports(Pair)

Adjustable plate used with fixed side support to increase width

Adjustable Plate for Fixed Supports (861102-PL) (Pair)	\$191
--	-------

### Straight Mounting Bracket

To be used on the extra cushion.

Straight Mounting Bracket for Extra Pad (861106-Kit)	\$148
--	-------

### Removeable Lateral Supports (Pair)



Removeable Lateral Supports, Size 1, A: 5.5" B: 5.5" C: 3.75" (861102-1) (Pair)	\$387
Removeable Lateral Supports, Size 2, A: 6.5" B: 7" C: 5" (861102-2) (Pair)	\$387
Removeable Lateral Supports, Size 3, A: 7.75" B: 7" C: 5" (861102-3) (Pair)	\$387

### Lateral Support Strap (Each)


For use with the Removeable Lateral Supports



Lateral Support Strap, Small: 10-12.5" (862103-1)(Each)	\$70
Lateral Support Strap, Medium: 13-16" (862103-2)(Each)	\$70
Lateral Support Strap, Large: 13.5-9.75" (862103-3) (Each)	\$70
Lateral Support Strap, X-Large: 17.75-23.5" (862103-4) (Each)	\$70

## Anatomic Head Support

If different style headrest desired, please include in notes section.




### Anatomic Head Support

Anatomic Head Support, Size Small (5099001)	\$348
Anatomic Head Support, Size Medium (5099002)	\$348
Anatomic Head Support, Size Large (5099003)	\$348

## Anatomic Head Support Cover

If Anatomic Head Support is selected, you must select a cover. Does not come standard.

Head Support Cover, Small (9923-83)	N/C
Head Support Cover, Medium (9923-83)	N/C
Head Support Cover, Large (9967-83)	N/C



### Headrest Mounting Hardware

Required with All Head Supports

Headrest Mounting Hardware, Short (98619-25)	\$314
Headrest Mounting Hardware, Long (98619-26)	\$314


Other headrest options available, please see Head supports order form. Please verify Compatibility.

## Head Pad



Head Pad, Size 1 and 2 (861116-1)	\$688
Head Pad, Size 3 and 4 (861116-2)	\$688

## Pommel



### Pommel

Pommel, Supine, A: 2", B: 2.75", C: 6", D: 3.5" (861103-1)	\$185
Pommel, Prone, A: 6", B: 2.75", C: 6", D: 6.75" (861103-2)	\$185




### Straps for Pommels

Straps for Pommel, Size 1, L: 9.84" - 11.81" (862104-1)	\$197
Straps for Pommel, Size 2, L: 10.63" - 18.5" (862104-2)	\$197

## Knee Supports

### Knee Support, Fixed



Knee Support, Standard Pad, Fixed Bracket, Size 1 and 2 (861144-2)	\$525
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 1 and 2 (861144-20)	\$729
Knee Support, Standard Pad, Fixed Bracket, Size 3 and 4 (861144-3)	\$525
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 3 and 4 (861144-30)	\$729

### Knee Support, Swing-Away




Knee Support, Standard Pad, Swing Away Bracket, Size 1 and 2 (861154-2)	\$905
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 1 and 2 (861154-20)	\$1109
Knee Support, Standard Pad, Swing Away Bracket, Size 3 and 4 (861154-3)	\$905
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 3 and 4 (861154-30)	\$1109

## Footplates and Straps



### Adjustable Sandals

Adjustable Sandals, Small (869490-RT)	\$226
Adjustable Sandals, Medium (869491-RT)	\$226
Adjustable Sandals, Large (869492-RT)	\$226



### Heel Support for Parted Foot Support


Heel Support for Parted Foot Support, Size 1 and 2 (861251) (Pair)	\$77
Heel Support for Parted Foot Support, Size 3 and 4 (861252) (Pair)	\$77

### Additional Straps for Adjustable Sandals

Additional Straps, Small/Medium (869497)	
Additional Straps, Large (869498)	\$92


### Upholstered Ankle Straps

Upholstered Ankle Straps, XS (85745-XS)	\$251
Upholstered Ankle Straps, Medium (85745-M)	\$251
Upholstered Ankle Straps, XL (85745-XL)	\$251



### Foot Straps

Foot Straps, Plastic Lock (85096) (Pair)	\$77
Foot Straps, Velcro (85102) (Pair)	\$70



### Neoprene Sandals

Neoprene Sandals, Small (85696-1) (Pair)	\$177
Neoprene Sandals, Medium (85696-2) (Pair)	\$177

## Tray Accessories



### Tray, Plastic, with Bowls and Cover

- Tray, Plastic, with Bowls and Cover, Size 1 and 2 (862130-2) (Each) \$763
- Tray, Plastic, with Bowls and Cover, Size 3 and 4 (862130-3) (Each) \$763



### Tray Armrests

- Armrest for Plastic Tray (861131) \$219

### Push Handle

- Attendant Push Handle (861108-U) \$191

### Mounting Kit

*For supine use only*

- Low Mount Kit for Plastic Tray (862131-KIT) \$280

### Padded Tray Tops

- Square Padded Tray Top (862135-1) \$564
- Small Rounded Tray Top \$564
- Large Rounded Tray Top \$564



### Toybar for Plastic Tray

- Toybar for Plastic Tray (861135) \$327