



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Caribou

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:



















Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/21/2025

Caribou

mm (inch)	kg (lb)				
					
 Overall width		510 mm (20")	570 mm (22½")	650 mm (25½")	750 mm (29½")
 Length		750 mm (29½")	920 mm (36¼")	1110 mm (43¾")	1110 mm (43¾")
 Overall height, vertical		460–570 mm (18–22½")	530–640 mm (20¾–25¼")	680–790 mm (26¾–31")	880–990 mm (34¾–39")
 Overall height, horizontal		700 mm (27½")	700 mm (27½")	800 mm (31½")	800 mm (31½")
 Angling		0°–87°	0°–87°	0°–87°	0°–87°
 Foot plate angle		(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)
 Total mass		20.5 kg (45.1 lb)	24 kg (52.9 lb)	30.5 kg (67.2 lb)	35 kg (77.1 lb)
 Maximum rated load		62.5 kg (137.7 lb)	86 kg (189.6 lb)	112.5 kg (248 lb)	137 kg (302 lb)
 Total mass, with power		22.5 kg (49.6 lb)	25.5 kg (56.2 lb)	32.5 kg (71.7 lb)	37 kg (81.6 lb)
 Maximum rated load, with power		64.5 kg (142.1 lb)	87.5 kg (192.9 lb)	114.5 kg (252.4 lb)	139 kg (306.4 lb)
 Maximum user mass		40 kg (88.1 lb)	60 kg (132.2 lb)	80 kg (176.3 lb)	100 kg (220.4 lb)
 User height		850–1120 mm (33½–44")	980–1290 mm (38½–50¾")	1250–1560 mm (49¾–61½")	1530–1840 mm (60¾–72½")
 Static stability downhill		max. 5°	max. 5°	max. 5°	max. 5°
 Static stability uphill		max. 5°	max. 5°	max. 5°	max. 5°
 Sidewards stability sideways		max. 5°	max. 5°	max. 5°	max. 5°

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1231.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Caribou

Caribou Base

Caribou Base

Required selection

Caribou Base, Size 1 (860101-FP)	\$1416
Caribou Base, Size 2 (860102-FP)	\$1461
Caribou Base, Size 3 (860113-FP)	\$1628
Caribou Base, Size 4 (860114-FP)	\$1920

Caribou Base with Power Actuator

Power Actuator recommended for Size 3 and 4.

Caribou Base with Power Actuator, Size 3 (860103-ELFP)	\$3448
Caribou Base with Power Actuator, Size 4 (860104-ELFP)	\$3765



Caribou Mast with Torso - Lower Body Pad

Required Selection

Lower Body Pad, Size 1 (X860101-P)	\$1160
Lower Body Pad, Size 2 (X860101-P)	\$1195
Lower Body Pad, Size 3 (X860113-P)	\$1334
Lower Body Pad, Size 4 (X860114-P)	\$1572


Multi Adjustable Footrests

Required Selection

Multi Adjustable Footrest, Size 1 (X860101-F)	\$403
Multi Adjustable Footrest, Size 2 (X860102-F)	\$434
Multi Adjustable Footrest, Size 3 (X860113-F)	\$459
Multi Adjustable Footrest, Size 4 (X860114-F)	\$485

Chest Pad


Chest Pad is a required selection.



Chest Pad, Size 1 A: 9.75", B: 12.25", C: 15.75" (861101-4)	\$395
Chest Pad, Size 2 A: 9.75", B: 12.25", C: 15.75" (861101-4)	\$425
Chest Pad, Size 3, A: 8.75", B: 10.5", C: 12.5" (861101-3)	\$450
Chest Pad, Size 4, A: 9.75", B: 12.25", C: 15.75" (861101-4)	\$475

Extra Cushion


Choose the same size as the Caribou Base.



Extra Cushion

Extra Cushion, Size 1, A: 9.75", B: 12.25", C: 15.75" (861101-4) (Each)	\$251
Extra Cushion, Size 2: W: 11.75" x H: 4" (861106-2) (Each)	\$251
Extra Cushion, Size 3, W: 13.75" x H: 4" (861106-3) (Each)	\$273
Extra Cushion, Size 4, W: 15.75" x H: 4" (861106-4) (Each)	\$273

Supports, Straps and Wedges



Wedges

Wedges for Side Slots, Size 1 (861107-1) (Pair)	\$17
Wedges for Side Slots, Size 2 (861107-2) (Pair)	\$17
Wedges for Side Slots, Size 3 (861107-3) (Pair)	\$17



Back Supports


Back Supports, Size 1 and 2 (861115-2) (Each)	\$322
Back Supports, Size 3 and 4 (861115-3) (Each)	\$328



Upper Arm Supports

Can not be used with shoulder supports. Can not be used with lateral supports.

Upper Arm Support, Small (861109-10)	\$410
Upper Arm Support, Large (861109-20)	\$437



Shoulder Support (Pair)

Pad size must match Chest Pad size.

Shoulder Support, Size 1, A: 3", B: 7.5", C: 2" (861126-1) (Pair)	\$279
Shoulder Support, Size 2, A: 3.25", B: 9.75", C: 2" (861126-2) (Pair)	\$279
Shoulder Support, Size 3, A: 4", B: 12.5", C: 2" (861126-3) (Pair)	\$279
Shoulder Support, Size 4, A: 4.25", B: 15.25", C: 2" (861126-4) (Pair)	\$279



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Caribou

Trunk Positioning- Vests

Measuring Shoulder to Pelvis Diagonally. See pictural under Belts and Straps.

Cross Vest



Cross Vest, XS, S-P: 7.75" - 12.5" (85720-XS) (Each)	\$208
Cross Vest, Small, S-P: 10.75" - 12.5" (85720-S) (Each)	\$208
Cross Vest, Medium, S-P: 12.5" - 16" (85720-M) (Each)	\$208
Cross Vest, Large, S-P: 16" - 18.5" (85720-L) (Each)	\$208
Cross Vest, XL, S-P: 20.5" - 23.5" (85720-XL) (Each)	\$208




Torso Vest (Ea)

Torso Vest, XS, S-P: 7.75"-10.75" (85705-XS)	\$219
Torso Vest, Small, S-P: 10.75"-12.5" (85705-S)	\$219
Torso Vest, Medium, S-P: 12.5"-16" (85705-M)	\$219
Torso Vest, Large, S-P: 16"-18.5" (85705-L)	\$219
Torso Vest, XL, S-P: 18.5"-20.5" (85705-XL)	\$219
Torso Vest, XXL, S-P: 20.5"-23.5" (85705-XXL)	\$219

Vest Mounting Hardware

Required Selection for All Vests



Fix Locks

Fix-Locks, Black, Bag of 6 Pcs (862109-35)	\$46
--	------

Belts and Straps



Chest Strap (Ea)

If multiple straps are required, please add to notes section below or contact customer service.




Chest Strap and Pelvic Strap, Small, W: 16.5" - 24.75" (862101-2)	\$131
Chest Strap and Pelvic Strap, Medium, W: 19" - 30.25" (862101-3)	\$131
Chest Strap and Pelvic Strap, Large, W: 22.75" - 37.75" (862101-4)	\$131



Cross Straps

Cross Strap, Size 1, 15" - 19.25" (862102-1) (Pair)	\$305
Cross Strap, Size 2, 17" - 21.5" (862102-2) (Pair)	\$305
Cross Strap, Size 3, 20.75" - 27.5" (Pair)	\$322
Cross Strap, Size 4, 24.75" - 33.75" (Pair)	\$322

Side Supports



Side Supports, Fixed (Pair)

Supports, Fixed, Size 1 and 2, A: 4.75" x B: 4.25" (861102-10) (Pair)	\$273
Side Supports, Fixed, Type 5, Size 3 and 4, A: 6" x 6.25" (861102-20) (Pair)	\$328



Removeable Lateral Supports (Pair)

Removeable Lateral Supports, Size 1, A: 5.5" B: 5.5" C: 3.75" (861102-1) (Pair)	\$312
Removeable Lateral Supports, Size 2, A: 6.5" B: 7" C: 5" (861102-2) (Pair)	\$312
Removeable Lateral Supports, Size 3, A: 7.75" B: 7" C: 5" (861102-3) (Pair)	\$312

Adjustable Plate for Fixed Side Supports(Pair)

Adjustable plate used with fixed side support to increase width

Adjustable Plate for Fixed Supports (861102-PL) (Pair)	\$155
--	-------

Straight Mounting Bracket

To be used on the extra cushion.

Straight Mounting Bracket for Extra Pad (861106-Kit)	\$97
--	------

Lateral Support Strap (Each)

For use with the Removeable Lateral Supports



Lateral Support Strap, Small: 10-12.5" (862103-1)(Each)	\$56
Lateral Support Strap, Medium: 13-16" (862103-2)(Each)	\$56
Lateral Support Strap, Large: 13.5-9.75" (862103-3) (Each)	\$56
Lateral Support Strap, X-Large: 17.75-23.5" (862103-4) (Each)	\$56




Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Caribou

Anatomic Head Support

If different style headrest desired, please include in notes section.




Anatomic Head Support

Anatomic Head Support, Size Small (5099001)	\$279
Anatomic Head Support, Size Medium (5099002)	\$279
Anatomic Head Support, Size Large (5099003)	\$279

Anatomic Head Support Cover

If Anatomic Head Support is selected, you must select a cover. Does not come standard.

Head Support Cover, Small (9923-83)	N/C
Head Support Cover, Medium (9923-83)	N/C
Head Support Cover, Large (9967-83)	N/C



Headrest Mounting Hardware

Required with All Head Supports

Headrest Mounting Hardware, Short (98619-25)	\$251
Headrest Mounting Hardware, Long (98619-26)	\$251


Other headrest options available, please see Head supports order form. Please verify Compatibility.

Head Pad



Head Pad, Size 1 and 2 (861116-1)	\$552
Head Pad, Size 3 and 4 (861116-2)	\$552

Pommel



Pommel

Pommel, Supine, A: 2", B: 2.75", C: 6", D: 3.5" (861103-1)	\$147
Pommel, Prone, A: 6", B: 2.75", C: 6", D: 6.75" (861103-2)	\$147




Straps for Pommels

Straps for Pommel, Size 1, L: 9.84" - 11.81" (862104-1)	\$160
Straps for Pommel, Size 2, L: 10.63" - 18.5" (862104-2)	\$160

Knee Supports

Knee Support, Fixed



Knee Support, Standard Pad, Fixed Bracket, Size 1 and 2 (861144-2)	\$421
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 1 and 2 (861144-20)	\$584
Knee Support, Standard Pad, Fixed Bracket, Size 3 and 4 (861144-3)	\$421
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 3 and 4 (861144-30)	\$584

Knee Support, Swing-Away



Knee Support, Standard Pad, Swing Away Bracket, Size 1 and 2 (861154-2)	\$727
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 1 and 2 (861154-20)	\$905
Knee Support, Standard Pad, Swing Away Bracket, Size 3 and 4 (861154-3)	\$727
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 3 and 4 (861154-30)	\$905

Footplates and Straps



Adjustable Sandals

Adjustable Sandals, Small (869490-RT)	\$180
Adjustable Sandals, Medium (869491-RT)	\$180
Adjustable Sandals, Large (869492-RT)	\$180

Heel Support for Parted Foot Support




Heel Support for Parted Foot Support, Size 1 and 2 (861251) (Pair)	\$61
Heel Support for Parted Foot Support, Size 3 and 4 (861252) (Pair)	\$61

Additional Straps for Adjustable Sandals

Additional Straps, Small/Medium (869497)	
Additional Straps, Large (869498)	\$62


Upholstered Ankle Straps

Upholstered Ankle Straps, XS (85745-XS)	\$202
Upholstered Ankle Straps, Medium (85745-M)	\$202
Upholstered Ankle Straps, XL (85745-XL)	\$202



Foot Straps

Foot Straps, Plastic Lock (85096) (Pair)	\$40
Foot Straps, Velcro (85102) (Pair)	\$56



Neoprene Sandals

Neoprene Sandals, Small (85696-1) (Pair)	\$141
Neoprene Sandals, Medium (85696-2) (Pair)	\$141



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Caribou

Tray Accessories



Tray, Plastic, with Bowls and Cover

Tray, Plastic, with Bowls and Cover, Size 1 and 2 (862130-2) (Each)	\$611
Tray, Plastic, with Bowls and Cover, Size 3 and 4 (862130-3) (Each)	\$611



Tray Armrests

Armrest for Plastic Tray (861131)	\$175
-----------------------------------	-------

Push Handle

Attendant Push Handle (861108-U)	\$155
----------------------------------	-------

Mounting Kit

For supine use only

Low Mount Kit for Plastic Tray (862131-KIT)	\$224
---	-------

Padded Tray Tops

Square Padded Tray Top (862135-1)	\$403
Small Rounded Tray Top	\$403
Large Rounded Tray Top	\$403



Toybar for Plastic Tray

Toybar for Plastic Tray (861135)	\$262
----------------------------------	-------