



# Chilla

Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:









Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 3/12/2025

## Chilla Base Frame

Frame size 1 & 2					
		mm (inch), kg (lb)			
		Size 1 		Size 2 	
		Normal	Suspension	Normal	Suspension
	Overall height, to pushbrace	1200 mm (47¼ inch)	1200 mm (47¼ inch)	1200 mm (47¼ inch)	1200 mm (47¼ inch)
	Overall width	600 mm (23½ inch)	620 mm (24½ inch)	660 mm (26 inch)	680 mm (26¾ inch)
	Maximum load	43 kg (94 lb)	43 kg (94 lb)	60 kg (132 lb)	60 kg (132 lb)
	Mass device**	9.2 kg (20 lb)	9.7 kg (21 lb)	9.5 kg (20 lb)	10 kg (22 lb)
	Minimum user mass in transportation	22 kg (48 lb)	22 kg (48 lb)	22 kg (48 lb)	22 kg (48 lb)
	Maximum user mass in transportation	35 kg (77 lb)	35 kg (77 lb)	50 kg (110 lb)	50 kg (110 lb)

\*380 mm (15 inch) without seat and rear wheels


















\*\*with solid wheels

\*\*\*frame without rear wheels



# Chilla

Phone (877) 227-0097 | Fax (704) 882-0751  
 sales.us@etac.com

Seat size 1 & 2			
		mm (inch), kg (lb)	
		Size 1 	Size 2 
	Backrest recline angle	80-160°	80-160°
			
	Backrest height with head support plate	505-665 mm (20-26¼ inch)	620-810 mm (24½-32 inch)
	Backrest height without head support plate	365-525 mm (14½-20¾ inch)	420-610 mm (16½-24 inch)
	Back width, top of backrest (with side supports)	150-225 mm (6-8¾ inch)	170-260 mm (6¾-10¼ inch)
	Back width, top of backrest (without side supports)	315 mm (12½ inch)	365 mm (14½ inch)
	Effective seat depth	180-340 mm (7-13½ inch)	240-420 mm (9½-16½ inch)
	Effective seat width (without hip supports)	330 mm (13 inch)	380 mm (15 inch)
	Effective seat width minimum/maximum (with hip supports)	180-260 mm (7-10¼ inch)	230-310 mm (9-12¼ inch)
	Foot support to seat	130-260 mm (5-10¼ inch)	185-420 mm (7¼-16½ inch)
	Mass device*	5.3 kg (11.7 lb)	6.3 kg (13.8 lb)
	Maximum user mass	35 kg (77 lb)	50 kg (110 lb)
	Maximum user mass in transportation	35 kg (77 lb)	50 kg (110 lb)
	Minimum user mass in transportation	22 kg (48 lb)	22 kg (48 lb)
	Seat recline angle	(-5°)-(40°)	(-5°)-(40°)

\*seat only (without accessories)



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com



# Chilla

## Base Frame

*Frame, Seat and Back size, must match*



Chilla Size 1 Frame Solid Adjustable Suspension (E1234) (881111110)	\$3682
Chilla Size 2 Frame Solid Adjustable Suspension (E1234) (881121110)	\$3682

## Transit

Transit Fittings (881090)	\$314
---------------------------	-------

## Storage Basket

Under Storage Basket (881092)	N/C
-------------------------------	-----



## Seat and Back

*Seat and Back size must match frame size. Seat and Back are both required selections.*

### Seat and Back Size 1

Size 1 Seat (E2292) (X881012)	\$609
Size 1 Back (E2291) (XB881012)	\$763

### Seat and Back Size 2

Size 2 Seat (E2292) (X881022)	\$637
Size 2 Back (E2291) (XB881022)	\$791

## Reclining Back

*Selection required*

Reclining Back (E1225) (X26000)	\$609
---------------------------------	-------

## Seat and Back Upholstery

### Back Cushion

#### Back Cushion

Back Cushion Size 1 Pink (881111-10)	N/C
Back Cushion Size 1 Blue (881111-40)	N/C
Back Cushion Size 1 Red (881111-50)	N/C
Back Cushion Size 1 Grey (881111-60)	N/C
Back Cushion Size 2 Blue (881121-40)	N/C
Back Cushion Size 2 Red (881121-50)	N/C
Back Cushion Size 2 Grey (881121-60)	N/C

#### Mesh Back Cushion

Mesh Back Cushion Size 1 Blue (881112-45)	\$42
Mesh Back Cushion Size 1 Red (881112-55)	\$42
Mesh Back Cushion Size 1 Grey (881112-65)	\$42
Mesh Back Cushion Size 2 Blue (881122-45)	\$42
Mesh Back Cushion Size 2 Red (881122-55)	\$42
Mesh Back Cushion Size 2 Grey (881122-65)	\$42

### Seat Cushion

#### Seat Cushion

Seat Cushion Size 1 Pink (881211-10)	N/C
Seat Cushion Size 1 Blue (881211-40)	N/C
Seat Cushion Size 1 Red (881211-50)	N/C
Seat Cushion Size 1 Grey (881211-60)	N/C
Seat Cushion Size 2 Blue (881221-40)	N/C
Seat Cushion Size 2 Red (881221-50)	N/C
Seat Cushion Size 2 Grey (881221-60)	N/C

#### Mesh Seat Cushion

Mesh Seat Cushion Size 1 Blue (881212-45)	\$42
Mesh Seat Cushion Size 1 Red (881212-55)	\$42
Mesh Seat Cushion Size 1 Grey (881212-65)	\$42
Mesh Seat Cushion Size 2 Blue (881222-45)	\$42
Mesh Seat Cushion Size 2 Red (881222-55)	\$42
Mesh Seat Cushion Size 2 Grey (881222-65)	\$42

### Infant Insert Cushion

*Available only with Size 1*

Infant Cushion Pink (K0108) (881001-10)	\$525
Infant Cushion Blue (K0108) (881001-20)	\$525
Infant Cushion Black (K0108) (881001-80)	\$525

### Back and Side Cover

*Selection required*

Cover for Back Sides Size 1 (K0108) (881111)	\$231
Cover for Back Sides Size 2 (K0108) (881112)	\$287



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com



# Chilla

## Footrests

*Mounting Rod and Footplate are required for complete footrest*

### Footrest Mounting Rod

Footrest Mounting Rod Short 7.5" (E0990) (881050-1)	\$154
Footrest Mounting Rod Long 11" (E0990) (881050-2)	\$154

### Angle and Height Adjustable Footplate

Angle and Height Adj Footplate - Short Bracket (K0040) (881055)	\$147
Angle and Height Adj Footplate - Long Bracket (K0040) (881056)	\$147

## Headrests

### Headrest Pads

*Headrest pads in this section do not require mounting hardware*

Contoured Head Cushion Size 1 (E0955) (874082-1)	\$300
Contoured Head Cushion Size 2 (E0955) (874082-2)	\$300
Sculpted Head Support Size 1 (E0955) (874084-1)	\$308
Sculpted Head Support Size 2 (E0955) (874084-2)	\$308
Two-Parted Head Support (E0955) (874083)	\$665
Anatomic Head Support Small (E0955) (874085-1)	\$371
Anatomic Head Support Medium (E0955) (874085-2)	\$385
Curved Head Support Small (E0955) (874086)	\$301



### Head Supports and Hardware

*Head supports in this section require Multi-Adjustable Headrest Hardware and Bracket for Head Support Mount*

Anatomic Head Support Small (E0955) (9923012-80)	\$371
Anatomic Head Support Medium (E0955) (9945012-80)	\$385
Curved Head Support Small (E0955) (9921012-80)	\$301

### Bracket for Head Support

Bracket for Head Support Mount (K0108) (881098)	\$77
---	------

### Multi-Adjustable Headrest Hardware

Multi-Adjustable Headrest Hardware (E1028) (X99XX012-80)	\$312
--	-------

## Postural Supports

### Hip Supports Pad and Hardware

Hip Support Pad Size 1 (E0956) (x8740250)	\$273
Hip Support Pad Size 2 (E0956) (x8740251)	\$315

### Hip Support Hardware

*Required when Hip Support Pad is selected*

Adjustable Removeable Hardware Size 1 (E1028) (881023-1)	\$245
Adjustable Removeable Hardware Size 2 (E1028) (881023-2)	\$245
Adjustable Removeable Rotational Hardware Size 1 (E1028) (881033-1)	\$630
Adjustable Removeable Rotational Hardware Size 2 (E1028) (881033-2)	\$630

### Pommel Hardware

*Required when Pommel is selected*

Pommel Hardware (E1028) (X881044)	\$245
-----------------------------------	-------

### Anti Thrust Wedge

Anti Thrust Wedges (K0108) (8720405)	\$301
--------------------------------------	-------

### Pommel

Pommel Size 1 (E0957) (881044-1)	\$203
Pommel Size 2 (E0957) (881044-2)	\$203

### Trunk Supports



Side Supports Size 1 (E0956) (881030-1)	\$308
Side Supports Size 2 (E0956) (881030-2)	\$308
Side Supports Size 3 (E0956) (881030-3)	\$308

### Outer Knee Supports and Hardware

Outer Knee Pad Size 1 (E0953) (881043-1)	\$273
Outer Knee Pad Size 2 (E0953) (881043-2)	\$301

### Adjustable Outer Knee Hardware

*Required when Outer Knee Pad is selected*

Outer Knee Support Adjustable Removable Hardware (E1028) (X881043-X)	\$245
--	-------



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

## Vests and Belts

### Torso Vest

Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$325
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$325
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$325
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$325

### Cross Vest

Cross Vest X-Small - S-P: 7.75" - 12.5" (Each) (85720-XS)(E0960)	\$308
Cross Vest Small - S-P: 10.25" - 13.5" (Each) (85720-S) (E0960)	\$308
Cross Vest Medium - S-P: 13.5" - 17.75" (Each) (85720-M) (E0960)	\$308
Cross Vest Large - S-P: 17.75" - 20.75" (Each) (85720-L) (E0960)	\$308
Cross Vest X-Large - S-P: 20.75" - 24.5" (Each) (85720-XL) (E0960)	\$308

### 2 Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$157
---	-------

### 4 Point Padded Hip Belt

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$218
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$218
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M)(E0960)	\$218
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$218

### Pelvic Harness

Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0) (K0108)	\$398
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1) (K0108)	\$398
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2) (K0108)	\$398

## Foot Positioning

### Foot Straps

Foot Straps, Plastic Lock (Pair) (K0040) (85096)	\$59
Foot Straps, Velcro (Pair) (K0040) (85102)	\$83

### Neoprene Sandals

Neoprene Sandals Small (85696-1) (Pair)	\$193
Neoprene Sandals Medium (85696-2) (Pair)	\$193

### Ankle Straps

Ankle Straps X-Small - C: 6.25"-9" (K0108) (85745-XS) (K0108)	\$300
Ankle Straps Medium - C: 9"-12.25" (K0108) (85745-M)	\$300
Ankle Straps X-Large - C: 11.5"-15" (K0108) (85745-XL)	\$300

## Accessories

### Canopy



For more information on canopy options, please see etac.com

Canopy, Size 1 (K0108) (881080-10)	N/C
Canopy, Size 2 (K0108) (881080-20)	N/C
Premium Canopy, Size 1 (K0108) (881080-30)	\$553
Premium Canopy, Size 2 (K0108) (881080-40)	\$553

### Additional Accessories

Bracket for Communication Aid (K0108) (881099)	\$235
Rain Cover, Size 1 (881093-1)	\$175
<i>Requires Canopy</i>	
Rain Cover, Size 2 (881093-2)	\$175
<i>Requires Canopy</i>	
Hand Brake Kit (881020)	\$525
Backpack Size 1 (88043184-1)	\$175
Backpack Size 2 (88043184-2)	\$175
Vent Tray (E1029) (881095)	\$462
<i>Not available with Under Storage Basket</i>	
O2 Holder (E2208) (881096)	\$245
IV Pole (K0105) (881097)	\$665

### Tray / Grip Bar and Hardware



Trays and Grip Bars use the same mounting hardware. Both may be selected; however, only one may be used at a time. Mounting hardware is required for both options.

Tray (E0950) (881071)	\$156
-----------------------	-------



Grip Bar Size 1 (K0108) (881070-1)	\$133
Grip Bar Size 2 (K0108) (881070-2)	\$133

### Tray and Grip Bar Hardware

Mounting Hardware for Tray or Grip Bar (E0128) (881060)	\$246.40
---	----------