

# R82 Cricket

## USD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


Effective: 1.10.24 - Rev 9

- Cricket
- Cushion
- Arm rests
- Head support
- Tray
- Grip bar
- Folding top
- Transport fittings



- A (Back Height) \_\_\_\_\_
- B (Arm Height) \_\_\_\_\_
- C (Lower Leg Length) \_\_\_\_\_
- D (Seat Depth) \_\_\_\_\_
- E (Chest Width) \_\_\_\_\_
- F (Seat Width) \_\_\_\_\_
- G (Height from Seat to Top of Head) \_\_\_\_\_
- H (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



#### Size 1 Inches

#### Size 2 Inches

Effective seat width (A)	13½"	15½"
Effective seat depth (B)	8½"-12¼"	12¼"-17"
Backrest height (C1)	19¾"	22½"
Seat surface height at front edge (E)	19¾"-21"	21¼"-22¼"
Footrest to seat distance (D)	7"-14¼"	8¾"-14½"
Leg to seat surface angle	360°	360°
Backrest width, outside (F1)	13½"	15½"
Backrest width, inside (F)	10¾"	13"
Overall width (G)	20¾"	22¾"
Length (H)	29½"	29½"
Overall height (I)	39¼"	41¼"
Folded width	14"	14"
Folded height	19"	19"
Folding length	30"	30"
Backrest angle	110°-125°	110°-125°
Minimum turning radius	27½"	30"
Total mass	27.3 lbs.	29½ lbs.
Max load	77 lbs.	110 lbs.
Max load at transportation	59½ lbs.	94½ lbs.
Required width of right angled corridor	32"	32"
Static stability (sideways, downhill, uphill)	10°	10°

### Notes:

## Cricket Standard Features




- Back Recline
- Adjustable Seat Depth
- Height and Angle Adjustable Foot Support
- Armrest, Fixed (Pr.)

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Cricket sz 1	87043100	\$2,044	E1236
<input type="checkbox"/> Cricket sz 2	87043200	\$2,044	E1236

Description	Size	Part No.	USD Price	HCPCS
-------------	------	----------	-----------	-------


**Required, must select a color**

- Blue = 1       Red = 2



- |                                  |   |          |     |
|----------------------------------|---|----------|-----|
| <input type="checkbox"/> Cushion | 1 | 8704311x | \$0 |
| <input type="checkbox"/> Cushion | 2 | 8704321x | \$0 |

◆ Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item



- |   |  |            |      |       |
|---|--|------------|------|-------|
| <input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆ |  | 87043130-1 | \$77 | K0108 |
| <input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆ |  | 87043130-2 | \$77 | K0108 |
- ◆ Requires Seat Side Parts (87043130)

- |  |          |          |      |
|--|----------|----------|------|
| <input type="checkbox"/> Seat Side Parts ◆ | One Size | 87043130 | \$55 |
|--|----------|----------|------|
- ◆ Required with Seat Width Reduction Kit

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



- |   |  |          |          |       |
|---|--|----------|----------|-------|
| <input type="checkbox"/> Transit Option |  | One Size | 87043188 | \$214 |
|---|--|----------|----------|-------|


Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



- |   |  |   |            |     |
|---|--|---|------------|-----|
| <input type="checkbox"/> Armrest, Fixed (Pr.) ◆ |  | 1 | 87043120-1 | \$0 |
| <input type="checkbox"/> Armrest, Fixed (Pr.) ◆ |  | 2 | 87043120-2 | \$0 |

◆ Required with Grip Bar, In Armrests




Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



- |   |               |    |          |       |       |
|---|---------------|----|----------|-------|-------|
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 7¼–10¾"  | XS | 85705-XS | \$232 | E0960 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 10¾–12½" | S  | 85705-S  | \$232 | E0960 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 12½–16"  | M  | 85705-M  | \$232 | E0961 |
| <input type="checkbox"/> Torso vest (Ea.) | S-P: 16–18½"  | L  | 85705-L  | \$232 | E0962 |



- |                                     |  |   |         |       |
|-------------------------------------|--|---|---------|-------|
| <input type="checkbox"/> Torso Belt |  | 1 | 85480-1 | \$148 |
| <input type="checkbox"/> Torso Belt |  | 2 | 85480-2 | \$148 |

Accessories	Description: Pelvic Positioning	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$112	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 6¾-12½"	XS	85725-XS	\$156	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 8¼-13¾"	S	85725-S	\$156	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 9¾-16¼"	M	85725-M	\$156	E0979
	<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦ ♦ <i>Required for use of Vests</i>		One Size	862109-35	\$45	

Description: Foot Positioning	Measurements	Size	Part No.	USD Price	HCPCS
-------------------------------	--------------	------	----------	-----------	-------




<input type="checkbox"/> Heel Straps (Pr.)		One Size	87043185	\$57	E0951
--	--	----------	----------	------	-------

Other types of vests and belts available, please see Vest and Belts Order Form

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Headrest Cushion, Black (Ea.)		One Size	87043150	\$214	E0955
--	--	----------	----------	-------	-------







<input type="checkbox"/> Backrest Height Extension (Ea.)		1	87043160-1	\$58	K0108
<input type="checkbox"/> Backrest Height Extension (Ea.)		2	87043160-2	\$58	K0108

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



<input type="checkbox"/> Grip Bar, In Armrests ♦		1	87043181-1	\$184	
<input type="checkbox"/> Grip Bar, In Armrests ♦		2	87043181-2	\$184	

♦ *Requires Armrest*

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		1	87043170	\$244	
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		2	87043270	\$244	
	<input type="checkbox"/> Underseat Storage Basket ♦		1	87043184-10	\$75	
	<input type="checkbox"/> Underseat Storage Basket ♦ ♦ <b>22 lbs. weight limit</b>		2	87043184-20	\$75	
	<input type="checkbox"/> Tray ♦		1	87043180-1	\$361	E0950
	<input type="checkbox"/> Tray ♦ ♦ <b>Must select same size as armrest</b>		2	87043180-2	\$361	E0950
	<input type="checkbox"/> Storage Strap		One Size	85100	\$33	

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Cricket complies with ISO7176-19, for safe transportation in motor vehicles. The R82 Cricket was tested with the Cricket back rest height extension and Y-hip belt and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

**Prices in USD are subject to change without notice**

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**

OF5100-ENG-REV.09

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

  
Creating Possibilities