



Account #: _____
 Contact: _____
 Bill To: _____
 Address: _____
 City/State/Zip: _____
 Phone: _____
 Email: _____


Purchase Order No.: _____
 Mark For: _____
 Ship To: _____
 Address: _____
 City/State/Zip: _____
 Phone: _____

Date: _____

Canada

Effective: 01.01.2022 - Rev 4

Measurements



No.	Type	Measurement
9	Total height	
10	Shoulder height	
11	Chest height / Armpit height	
15	Top of pelvic height	
18	Chest circumference	
19	Pelvic circumference	

		0	1	2	3
	Width between handles	4 - 9¾"	7½ - 13¾"	9½ - 16½"	13 - 19"
	Height from floor to handles	11¾ - 24"	16½ - 28¾" 15 - 21¼"*	21¾ - 33"	29½ - 41¼"
	Max. user height	41¼"	45¼"	59"	70¾"
	Max. user mass	66 lbs	66 lbs	99 lbs	176 lbs
	Length	25¼"	28¼"	33"	34¼"
	Overall width	22¾"	24¾"	27¼"	27¼"
	Folded length	25½"	28½"	33¾"	39¼"
	Folded width	22¾"	24½"	27¼"	27¼"
	Folded height	11¾"	12¼"	10½"	11½"
	Minimum turning radius	19"	19¾"	22"	45"
	Mass device	13.4 lbs	15.7 lbs	15.7 lbs	26.7 lbs

*Mounted with the Starter Kit accessories. Further measurements at our website.

R82 Crocodile

AADL CAD Order Form



Features, incl. in standard price

- Width and height adjustable handles
- Easy to fold
- Lightweight durable aluminium frame – Sizes 0, 1 and 2
- Stable steel and aluminium frame – Size 3
- Upholstery and cover for top frame
- Solid front and rear wheels
- Swivel front wheels with directional stabilizer
- Tension adjustable rear wheels
- On/off anti-reverse wheel stop
- Padded Rear Stabilizer Bar w/ Tool Kit



Description	Measurements	Size	Article number	CAD
-------------	--------------	------	----------------	-----

Main products

Crocodile Color and Size

	Crocodile, white Standard		Size 0	86800	872.90 <input type="checkbox"/>
	Crocodile, yellow Standard		Size 1	86801	963.00 <input type="checkbox"/>
	Crocodile, orange Standard		Size 2	86802	963.00 <input type="checkbox"/>
	Crocodile, metallic black Standard		Size 2	86852	963.00 <input type="checkbox"/>
	Crocodile, anthracite grey Standard		Size 3	86803	1,586.00 <input type="checkbox"/>

Description	Measurements	Size	Article number	CAD
-------------	--------------	------	----------------	-----

Accessories

Mounting kits

	Crocodile Starter Kit <i>For Crocodile size 1.</i>		One Size	86801-0	569.00 <input type="checkbox"/>
--	--	--	----------	---------	---------------------------------








Anti-tippers

	Anti-tip - Swing-away Anti-tip Swing-away Pair, for size 0-2		Size 1	86815	209.00 <input type="checkbox"/>
	Anti-tip Swing-away, Pair, for size 3		Size 2	86815-3S	209.00 <input type="checkbox"/>
	Anti-tip - Fixed Anti-tip fixed, Pair - for size 1-2		Size 1	86815-F	101.00 <input type="checkbox"/>
	Anti-tip fixed, Pair - for size 3		Size 2	86815-3F	107.00 <input type="checkbox"/>



Description	Measurements	Size	Article number	CAD
Arm supports				
	Forearm support w/handgrip			
	Forearm support, Pair - for size 0-1	L:5¾ x W:3"	Size 1	86858-1 481.00 <input type="checkbox"/>
	Forearm support, Pair - for size 2	L:7 x W:4"	Size 2	86859-1 481.00 <input type="checkbox"/>
	Forearm support, Pair - for size 3	L:7 x W:4"	Size 3	86860-1 481.00 <input type="checkbox"/>
	Padded forearm positioner strap*			
	Padded forearm positioner strap, EA		One Size	869370 89.00 <input type="checkbox"/>
	Padded forearm positioner strap, EA		One Size	869370 89.00 <input type="checkbox"/>
	<i>*For use with Forearm support w/ handgrip. For quantity 2 select both options above.</i>			
	Handle knobs - Pair		One Size	86861 66.00 <input type="checkbox"/>
Brakes				
	Hand Brake Kit*			
	Hand Brake Kit, for size 1-2		Size 1	86834 312.00 <input type="checkbox"/>
	Hand Brake Kit, for size 3		Size 2	86873 312.00 <input type="checkbox"/>
	<i>*Hand brake eliminates the function of anti-reverse wheel stop</i>			
	Decelerating Hand Brake Kit*			
	Decelerating Hand Brake Kit, for sizes 1-2		Size 1	86834-DEC 312.00 <input type="checkbox"/>
	Decelerating Hand Brake Kit, for size 3		Size 2	86873-DEC 312.00 <input type="checkbox"/>
	<i>*Hand brake eliminates the function of anti-reverse wheel stop</i>			
	<i>*If ordering hand brake with anti-reverse wheel stop function, a custom order required</i>			
Back and pelvic supports				
	Surrounding Support*		One Size	868610 637.22 <input type="checkbox"/>
	<i>*Only available with size 0</i>			
	Pelvic Support w/ Laterals			
	Pelvic Support w/ Laterals, for sizes 0-2	3½" x 4¾" - 9"	Small	86808 341.00 <input type="checkbox"/>
	Pelvic Support w/ Laterals, for sizes 1-2	3½" x 7" - 11½"	Large	86809 341.00 <input type="checkbox"/>
	Lateral strap*		One Size	905011 66.00 <input type="checkbox"/>
	<i>*For sizes 1 and 2. Back and side support size small and large</i>			
	Back support - Adjustable Padded*	9½" - 11½" x 7"		86886 657.00 <input type="checkbox"/>
	<i>*For Size 3</i>			
	Abdominal Support Belt*	13½" x 4"	-	85490-2 113.00 <input type="checkbox"/>
	<i>*For Back Support #86886, for size 3</i>			
	Cross bar for mounting back/pelvic supports, for Size 3		-	86880 126.00 <input type="checkbox"/>
	<i>*Required with all Back and Pelvic Supports for size 3.</i>			



Description	Measurements	Size	Article number	CAD
Handles and push braces				
 Attendant Handle* <i>*Not available with size 3.</i>		One Size	86866-1	180.82 <input type="checkbox"/>
 Grip Bar Grip Bar, for size 0 Grip Bar, for size 1		Size 0 Size 1	86879-0 86879-1	168.35 <input type="checkbox"/> 193.29 <input type="checkbox"/>
Postural supports				
 Hip Supports, Pair* <i>*Not available with size 3</i>		One Size	86819	305.52 <input type="checkbox"/>
Seating				
Flip Down Seat White Yellow* Orange Metallic Black Anthracite Grey <i>*Not compatible with 86801-0, Crocodile Starter Kit</i>		Size 0 Size 1 Size 2 Size 2 Size 3	86820-01 86820 86821 86823 86863	220.72 <input type="checkbox"/> 221.00 <input type="checkbox"/> 221.00 <input type="checkbox"/> 221.00 <input type="checkbox"/> 221.00 <input type="checkbox"/>
 Lightweight Sling Seat* Lightweight Sling Seat, for size 0-1 Lightweight Sling Seat, for size 2 <i>*Requires Pelvic Support w/laterals 8680X and strap 905011 or 868610 for Size 0.</i>	23½" (Waist) 27¼" (Waist)	Size 1 Size 2	86824-1 86825-1	176.00 <input type="checkbox"/> 176.00 <input type="checkbox"/>
 Sling seat Sling Seat, for size 1 Sling Seat, for size 2	23½" (Waist) 27¼" (Waist)	Size 1 Size 2	86824 86825	176.00 <input type="checkbox"/> 176.00 <input type="checkbox"/>
Bags and baskets				
 Bag Storage Bag, for sizes 0-2		One Size	86839	93.00 <input type="checkbox"/>
 Bag Storage Bag, for size 3		One Size	86843	273.00 <input type="checkbox"/>

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice

For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5202-ENG-REV.04

etac[®]
Creating Possibilities