



Phone (877) 227-0097 | Fax (704) 882-0751



# Flamingo

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:

Marked For:







Email/Fax:







Associated Quote/Order#:

Effective: 1/20/2025

## Flamingo

mm (inch) kg (lb)

	1	2	3	4
 Overall width, seat	350 mm (13¾")	380 mm (15")	430 mm (17")	480 mm (19")
 Length, seat	260 mm (10¼")	315 mm (12½")	360 mm (14¼")	425 mm (16¾")
 Overall height, seat	400 mm (15¾")	430 mm (17")	480 mm (19")	510 mm (20")
 Effective seat width	240 mm (9½")	280 mm (11")	320 mm (12½")	340 mm (13½")
 Effective seat depth	250 mm (9¾")	290 mm (11½")	350 mm (13¾")	400 mm (15¾")
 Back width, top of backrest	200 mm (7¾")	220 mm (8¾")	260 mm (10¼")	280 mm (11")
 Backrest height	350 mm (13¾")	380 mm (15")	420 mm (16½")	470 mm (18½")
 Folded height, seat	240 mm (9½")	260 mm (10¼")	270 mm (10¾")	270 mm (10¾")
 Folded length, seat	370 mm (15")	400 mm (15¾")	430 mm (17")	460 mm (18")
 Backrest angle	30°	30°	30°	30°
 Mass device, seat	2.5 kg (5.5 lb)	3 kg (6.5 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max rated load	39.5 kg (87 lb)	55 kg (121.2 lb)	66 kg (145.5)	76.5 kg (168.6 lb)

kg (lb)	1	2	3	4
 Mass device, frame	11,5 kg (25.3 lb)	11,5 kg (25.3 lb)	12,5 kg (27.5 lb)	12,5 kg (27.5 lb)
 Mass device, seat	2,5 kg (5.5 lb)	3 kg (6.6 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Mass device	14 kg (30.8 lb)	14,5 kg (31.9 lb)	16,5 kg (36.3 lb)	17 kg (37.4 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max. load, operating force	35 kg (77.1 lb)	35 kg (77.1 lb)	55 kg (121.2 lb)	55 kg (121.2 lb)
 Max rated load	51 kg (112.4 lb)	66,5 kg (146.6 lb)	78,5 kg (173 lb)	89 kg (196 lb)

## Seat and Back



### Seat and Back

Seat and Back, Blue, Size 1 (880101-61)	\$962
Seat and Back, Blue, Size 2 (880102-61)	\$1028
Seat and Back, Blue, Size 3 (880103-61)	\$1095
Seat and Back, Blue, Size 4 (880104-61)	\$1161

## Flamingo Mobility Base

Footrest only included with the High-low frame.

### Flamingo Frame



Flamingo Frame W/ 3" Casters, Size 1 (880201-NL)	\$1128
Flamingo Frame W/ 3" Casters, Size 2 (880202-NL)	\$1188
Flamingo Frame W/ 3" Casters, Size 3 (880203-NL)	\$1293
Flamingo Frame W/ 3" Casters, Size 4 (880204-NL)	\$1388

### High-Low Mobility Frame



High-Low Bath Frame with 4" Casters, Size 1 (880411)	\$3007
High-Low Bath Frame with 4" Casters, Size 2 (880412)	\$3007
High-Low Bath Frame with 4" Casters, Size 3 (880413)	\$3007
High-Low Bath Frame with 4" Casters, Size 4 (880414)	\$3089



### Potty Frame

Potty Frame, Size 1 (880160-1)	\$509
Potty Frame, Size 2 (880160-2)	\$509

### Fitting for Direct Toilet Mounting

Direct Toilet Mount, Size 1 (880141)	\$830
Direct Toilet Mount, Size 2 (880142)	\$830
Direct Toilet Mount, Size 3 (880143)	\$830
Direct Toilet Mount, Size 4 (880144)	\$830

## Footrest for Mobility Base

*If wanting a footrest for Flamingo frame or Direct Toilet Mounting, please quote the footrest below.*



Flamingo Foot Support Pair, Size 1 (880211)	\$273
Flamingo Foot Support Pair, Size 2 (880212)	\$273
Flamingo Foot Support Pair, Size 3 (880213)	\$273
Flamingo Foot Support Pair, Size 4 (880214)	\$273

## Calf and Footrest Supports

Calf Support (880199)	\$89
-----------------------	------



## Footrest Support

Footrests Straps (Pair) (85096)	\$77
---------------------------------	------

## Pelvic Positioning

### Upholstery, Seat w/wide sides, Blue



Upholstery, Seat w/wide sides, Blue, Size 2 (880166-61)	\$142
Upholstery, Seat w/wide sides, Blue, Size 3 (880167-61)	\$150
Upholstery, Seat w/wide sides, Blue, Size 4 (880168-61)	\$158



### Splashguard, Low Profile

Splashguard, Low Profile, Size 1, Blue (880195-61)	\$106
Splashguard, Low Profile, Size 2,3,4, Blue (880196-61)	\$106



### Seat Infill for Seat

Seat Infill for Seat, Blue, Size 1 (880185-61)	\$75
Seat Infill for Seat, Blue, Size 2 (880186-61)	\$75
Seat Infill for Seat, Blue, Size 3 (880187-61)	\$82
Seat Infill for Seat, Blue, Size 4 (880188-61)	\$82



### Splashguard, High Profile

Splashguard, High Profile, Size 1, Blue (880191-61)	\$112
Splashguard, High Profile, Size 2,3,4, Blue (880192-61)	\$112

## Suction Cups for Seat


Suction Cups for Seat (880150)	\$185
--------------------------------	-------

## Seat Belt w/ Plastic Lock

Seat Belt w/ Plastic Lock (85100)	\$41
-----------------------------------	------

## Trunk Positioning

### Sides for Back W/ Upholstery, Blue (Each)



Sides for Back W/ Upholstery, Blue (Each), Size 1 (880111-61)	\$177
Sides for Back W/ Upholstery, Blue (Each), Size 2 (880112-61)	\$177
Sides for Back W/ Upholstery, Blue (Each), Size 3 (880113-61)	\$177
Sides for Back W/ Upholstery, Blue (Each), Size 4 (880114-61)	\$191




### Vest, Water-Proof Mesh (Each)

Vest, Water-Proof Mesh (Each), Size 3/4 (89111-07540)	\$273
Vest, Water-Proof Mesh (Each), Size 1 (89111-140)	\$273
Vest, Water-Proof Mesh (Each), Size 2 (89111-240)	\$273

## Cross vest, Water-proof mesh (net)

Cross vest, Water-proof mesh, Size 0 (89625-040)	\$278
Cross vest, Water-proof mesh, Size 1 (89625-140)	\$278
Cross vest, Water-proof mesh, Size 2 (89625-240)	\$278
Cross vest, Water-proof mesh, Size 3 (89625-340)	\$278
Cross vest, Water-proof mesh, Size 4 (89625-440)	\$278



## Back Reinforcement Kit

Back Reinforcement Kit 3 and 4 (880155)	\$348
---	-------



## Back Extension

Trunk Positioning, Back Extension, Blue Size 1 and 2 (880118-61)	\$232
Trunk Positioning, Back Extension, Blue, Size 3 and 4 (880119-61)	\$232

## Tilt Support Rod

Tilt Support Rod, Size 3 and 4	\$199
--------------------------------	-------

## Mounting Rod for Back Extension L:8"

Mounting Rod for Back Extension L:8" (880120)	\$64
---	------



Phone (877) 227-0097 | Fax (704) 882-0751



# Flamingo


## Fixed Side Supports



Fixed Side Supports, Gray Cover, Set, Size 1, 3.5" x 4" (88401-10) \$319  
Fixed Side Supports, Gray Cover, Set, Size 2, 4.75" x 5.5" (88401-20) \$319

## Bar, Straight, Stainless

*Qty 2 must be quoted*



Bar, Straight, Stainless L: 4" (82569) N/C  
Bar, Straight, Stainless L: 6" (82648) N/C

## Swing-Away Brackets

*Must quote swing away brackets, bar and cushion for a complete Swing away lateral.*

Swing-Away Bracket, Right, Each (83300-R) \$308  
Swing-Away Bracket, Right, Offset (83500-R) \$308

Swing-Away Bracket, Left, Each (83300-L) \$308  
Swing-Away Bracket, Left, Offset (83500-L) \$308

## Cushion, Curved

*Qty 2 must be quoted*




Cushion, Curved, Blue, Each, L: 5.75" x H: 2.5" (6075409) \$142  
Cushion, Curved, Blue, Each, L: 6.5" x H: 2.5" (6075408) \$142  
Cushion, Curved, Blue, Each, L: 7.5" x H: 2.75" (6075407) \$142

## Push Handles



Push Brace (One Size) (880110) \$129



### Helper Lift Handle

*For use with High-low Base*

Helper Lift Handle, Size 1 and 2 (880110-HL) \$117  
Helper Lift Handle, Size 3 and 4 (880110-HLL) \$117

## Protection Pad Covers

Protection Pad Covers, Gray (89151) \$126

## Armrests



Swing-Away Armrest/Tray, Blue, Size 1 (880131-61) \$303  
Swing-Away Armrest/Tray, Blue, Size 2 (880132-61) \$303  
Swing-Away Armrest/Tray, Blue, Size 3 (880133-61) \$303  
Swing-Away Armrest/Tray, Blue, Size 4 (880134-61) \$303




## Separate Armrests

Separate Armrests, Set, 9" (880152) \$240

## Head Support



Head Support No. 1, Blue PU (6071000) \$280  
Head Support No. 2, Blue PU (6072000) \$280  
Head Support No. 3, Blue PU (6073000) \$280  
Head Support No. 4, Blue PU (6074000) \$280



### Turtle Bar System

*Turtle bar required for headrest*

Turtle Bar System, Short, 5" (98619-24) \$314

## Commode Feature



### Commode pan w/holder

Commode pan w/holder, Size 1 (880121) \$164  
Commode pan w/holder, Size 2 (880122) \$164  
Commode pan w/holder, Size 3 (880123) \$164  
Commode pan w/holder, Size 4 (880124) \$164

## Potty Ring



Potty Ring W/ Holders, Size 1 (880121-PR) \$177  
Potty Ring W/ Holders, Size 2 (880122-PR) \$177  
Potty Ring W/ Holders, Size 3 (880123-PR) \$177  
Potty Ring W/ Holders, Size 4 (880124-PR) \$177