



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:


Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025

Flamingo







mm (inch) kg (lb)

	 1	 2	 3	 4
 Overall width, seat	350 mm (13¾")	380 mm (15")	430 mm (17")	480 mm (19")
 Length, seat	260 mm (10¼")	315 mm (12½")	360 mm (14¼")	425 mm (16¾")
 Overall height, seat	400 mm (15¾")	430 mm (17")	480 mm (19")	510 mm (20")
 Effective seat width	240 mm (9½")	280 mm (11")	320 mm (12½")	340 mm (13½")
 Effective seat depth	250 mm (9¾")	290 mm (11½")	350 mm (13¾")	400 mm (15¾")
 Back width, top of backrest	200 mm (7¾")	220 mm (8¾")	260 mm (10¼")	280 mm (11")
 Backrest height	350 mm (13¾")	380 mm (15")	420 mm (16½")	470 mm (18½")
 Folded height, seat	240 mm (9½")	260 mm (10¼")	270 mm (10¾")	270 mm (10¾")
 Folded length, seat	370 mm (15")	400 mm (15¾")	430 mm (17")	460 mm (18")
 Backrest angle	30°	30°	30°	30°
 Mass device, seat	2.5 kg (5.5 lb)	3 kg (6.5 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max rated load	39.5 kg (87 lb)	55 kg (121.2 lb)	66 kg (145.5)	76.5 kg (168.6 lb)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1274.2 Pricing and specifications subject to change without notice.

kg (lb)	1	2	3	4
 Mass device, frame	11,5 kg (25.3 lb)	11,5 kg (25.3 lb)	12,5 kg (27.5 lb)	12,5 kg (27.5 lb)
 Mass device, seat	2,5 kg (5.5 lb)	3 kg (6.6 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Mass device	14 kg (30.8 lb)	14,5 kg (31.9 lb)	16,5 kg (36.3 lb)	17 kg (37.4 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max. load, operating force	35 kg (77.1 lb)	35 kg (77.1 lb)	55 kg (121.2 lb)	55 kg (121.2 lb)
 Max rated load	51 kg (112.4 lb)	66,5 kg (146.6 lb)	78,5 kg (173 lb)	89 kg (196 lb)

Base Frame and Mobility Base Options

Seat and Back



Seat and Back

Seat and Back, Blue, Size 1 (880101-61)	\$740
Seat and Back, Blue, Size 2 (880102-61)	\$791
Seat and Back, Blue, Size 3 (880103-61)	\$842
Seat and Back, Blue, Size 4 (880104-61)	\$893

Flamingo Mobility Base

Footrest only included with the High-low frame



Flamingo Frame W/ 3" Casters, Size 1 (880201-NL)	\$867
Flamingo Frame W/ 3" Casters, Size 2 (880202-NL)	\$913
Flamingo Frame W/ 3" Casters, Size 3 (880203-NL)	\$995
Flamingo Frame W/ 3" Casters, Size 4 (880204-NL)	\$1066



Potty Frame, Size 1 (880160-1)	\$459
Potty Frame, Size 2 (880160-2)	\$505



High-Low Bath Frame with 4" Casters, Size 1 (880411)	\$2412
High-Low Bath Frame with 4" Casters, Size 2 (880412)	\$2412
High-Low Bath Frame with 4" Casters, Size 3 (880413)	\$2412
High-Low Bath Frame with 4" Casters, Size 4 (880414)	\$2478



Fitting for Direct Toilet Mounting

Direct Toilet Mount, Size 1 (880141)	\$638
Direct Toilet Mount, Size 2 (880142)	\$638
Direct Toilet Mount, Size 3 (880143)	\$638
Direct Toilet Mount, Size 4 (880144)	\$638



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Footrest for Mobility Base

If wanting a footrest for the flamingo frame or Direct Toilet Mounting, please quote the footrests below



Flamingo Foot Support Pair, Size 1 (880211)	\$306
Flamingo Foot Support Pair, Size 2 (880212)	\$332
Flamingo Foot Support Pair, Size 3 (880213)	\$357
Flamingo Foot Support Pair, Size 4 (880214)	\$383



Calf and Footrest Supports

Calf Support (880199) _____ \$73

Footrests Straps (Pair) (85096) _____ \$40



Pelvic and Trunk Positioning

Pelvic Positioning



Seat with Upholstery

Upholstery, Seat w/wide sides, Blue, Size 2 (880166-61)	\$116
Upholstery, Seat w/wide sides, Blue, Size 3 (880167-61)	\$121
Upholstery, Seat w/wide sides, Blue, Size 4 (880168-61)	\$126



Seat Infill

Seat Infill for Seat, Blue, Size 1 (880185-61)	\$60
Seat Infill for Seat, Blue, Size 2 (880186-61)	\$60
Seat Infill for Seat, Blue, Size 3 (880187-61)	\$67
Seat Infill for Seat, Blue, Size 4 (880188-61)	\$67

Splashguard

Splashguard, Low Profile, Size 1, Blue (880195-61)	\$84
Splashguard, Low Profile, Size 2,3,4, Blue (880196-61)	\$84
Splashguard, High Profile, Size 1, Blue (880191-61)	\$89
Splashguard, High Profile, Size 2,3,4, Blue (880192-61)	\$89

Accessories

Seat Belt w/ Plastic Lock (85100)	\$34
Suction Cups for Seat (880150)	\$147

Neoprene Hip Belts

Neoprene Hip Belts, size 1 (85423-1)	\$107
Neoprene Hip Belts, size 2 (85423-2)	\$107
H-Harness with Plastic Lock (85105)	\$67

Trunk Positioning



Sides with Upholstery

Sides for Back W/ Upholstery, Blue (Each), Size 1 (880111-61)	\$141
Sides for Back W/ Upholstery, Blue (Each), Size 2 (880112-61)	\$141
Sides for Back W/ Upholstery, Blue (Each), Size 3 (880113-61)	\$141
Sides for Back W/ Upholstery, Blue (Each), Size 4 (880114-61)	\$155



Back Extension

Trunk Positioning, Back Extension, Blue Size 1 and 2 (880118-61)	\$225
Back Extension, Size 3 and 4 (880119-61)	\$225

Back Reinforcement Kit

Back Reinforcement Kit 3 and 4 (880155)	\$279
---	-------

Cross vest, Water-proof mesh (net)

Cross vest, Water-proof mesh, Size 0 (89625-040)	\$214
Cross vest, Water-proof mesh, Size 1 (89625-140)	\$214
Cross vest, Water-proof mesh, Size 2 (89625-240)	\$214
Cross vest, Water-proof mesh, Size 3 (89625-340)	\$214
Cross vest, Water-proof mesh, Size 4 (89625-440)	\$214

Mounting Rod

Rod needed for Back Extension

Mounting Rod for Back Extension L:8" (880120)	\$50
---	------

Tilt Support Rod

Tilt Support Rod, Size 3 and 4	\$153
--------------------------------	-------

Vest, Water-Proof Mesh (Each)

More options can be found at etac.com

Vest, Water-Proof Mesh (Each), Size 3/4 (89111-07540)	\$219
Vest, Water-Proof Mesh (Each), Size 1 (89111-140)	\$219
Vest, Water-Proof Mesh (Each), Size 2 (89111-240)	\$219
Vest, Water-Proof Mesh (Each), Size 3 (89111-340)	\$214
Vest, Water-Proof Mesh (Each), Size 4 (89111-440)	\$214



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Fixed Side Supports



Fixed Side Supports, Gray Cover, Set, Size 1, 3.5" x 4" (88401-10)	\$255
Fixed Side Supports, Gray Cover, Set, Size 2, 4.75" x 5.5" (88401-20)	\$255

Swing-Away Brackets

Swing-Away Bracket, Right, Each (83300-R)	\$246
Swing-Away Bracket, Right, Offset (83500-R)	\$246
Swing-Away Bracket, Left, Each (83300-L)	\$246
Swing-Away Bracket, Left, Offset (83500-L)	\$246

Headrests and Hardware

Head Support No. 1, Blue PU (6071000)	\$224
Head Support No. 2, Blue PU (6072000)	\$224
Head Support No. 3, Blue PU (6073000)	\$224
Head Support No. 4, Blue PU (6074000)	\$224

Cushion, Curved

Cushion, Curved, Blue, Each, L: 5.75" x H: 2.5" (6075409)	\$116
Cushion, Curved, Blue, Each, L: 6.5" x H: 2.5" (6075408)	\$116
Cushion, Curved, Blue, Each, L: 7.5" x H: 2.75" (6075407)	\$116

Bar, Straight, Stainless

Bar, Straight, Stainless L: 4" (82569)	N/C
Bar, Straight, Stainless L: 6" (82648)	N/C

Protection Pad Covers

For Swing-Away Brackets

Protection Pad Covers

Protection Pad Covers, Gray (89151)	\$279
-------------------------------------	-------

Armrest Options

Armrests



Swing-Away Armrest/Tray, Blue, Size 1 (880131-61)	\$250
Swing-Away Armrest/Tray, Blue, Size 2 (880132-61)	\$250
Swing-Away Armrest/Tray, Blue, Size 3 (880133-61)	\$250
Swing-Away Armrest/Tray, Blue, Size 4 (880134-61)	\$250



Separate Armrests, Set, 9" (880152)	\$195
-------------------------------------	-------

Accessories

Push Handles



Push Brace (One Size) (880110)	\$105
--------------------------------	-------



For use with High-Low Base

Helper Lift Handle, Size 1 and 2 (880110-HL)	\$93
Helper Lift Handle, Size 3 and 4 (880110-HLL)	\$93

Potty Options



Potty Ring W/ Holders, Size 1 (880121-PR)	\$141
Potty Ring W/ Holders, Size 2 (880122-PR)	\$141
Potty Ring W/ Holders, Size 3 (880123-PR)	\$141
Potty Ring W/ Holders, Size 4 (880124-PR)	\$141



Commode pan w/holder

Commode pan w/holder, Size 1 (880121)	\$131
Commode pan w/holder, Size 2 (880122)	\$131
Commode pan w/holder, Size 3 (880123)	\$131
Commode pan w/holder, Size 4 (880124)	\$131