

R82 Heron

CAD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:


US Configuration

Effective: 01.21.25 - Rev 6




A (Back height) _____
 B (Arm height) _____
 C (Lower leg length) _____
 D (Seat depth) _____
 E (Chest width) _____
 F (Seat width) _____
 G (Height from seat to top of head) _____
 H (Shoulder width) _____
 Age _____
 Weight _____

Measurements



	Heron with hydraulics	Heron with power
Effective seat width (A)	550 mm / 21¾"	550 mm / 21¾"
Seat width between hip supports (accessories)	350-490 mm / 13¾ - 19¼"	350-490 mm / 13¾ - 19¼"
Seat depth with reduction kit (accessories)		
Backrest height (C)	555-625 mm / 21¾-24½"	555-625 mm / 21¾-24½"
Lower leg length (D)	400-560 mm / 15¾-22"	400-560 mm / 15¾-22"
Seat surface height at front edge (E)	510-780 mm / 20-30¾"	510-780 mm / 20-30¾"
Back width, top of backrest (F)	290 mm / 11½"	290 mm / 11½"
Overall width (G)	670 mm / 26½"	670 mm / 26½"
Length (H)	710 mm / 28"	710 mm / 28"
Overall height (I)	1080-1350 mm / 44½-53¾"	1080-1350 mm / 44½-53¾"
Seat plane angle	-5-25°	-5-25°
Back rest angle	0-30°	0-30°
Mass device	28 kg / 61.75 lb	31 kg / 68.25 lb
Max rated load	130 kg / 286.6 lb	163 kg / 359.3 lb


Heron Standard Features



- One size
- 4 grey castors with locks
- Stainless steel frame
- Height adjustable frame
- Integrated push handles
- Built-in tilt in space
- Built-in recline
- Active, height adjustable back
- Potty ring

Description	Size	Part No.	CAD Price
<input type="checkbox"/> Heron	One Size	880505-01	\$7,143
<input type="checkbox"/> Heron Power	One Size	880505-01EL	\$9,094


Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Foot Support, Right-Angled Knee (Pr.) Lower leg range: 300mm-380mm	Short	880506-03	\$675
<input type="checkbox"/> Foot Support, Right-Angled Knee (Pr.) Lower leg range: 390mm-550mm	Long	880506-01	\$675
<input type="checkbox"/> Foot Support, Open Knee Angle (Pr.) Lower leg range: 360mm-530mm	One Size	880506-02	\$675
<input type="checkbox"/> Padding for Foot Supports (Pr.)	One Size	880506-12	\$88
<input type="checkbox"/> Calf Support, Strap (Ea.)	4"	880509-01	\$88
<input type="checkbox"/> Heel Stops, Straps (Pr.)	1½"	880509-02	\$164

Description	Size	Part No.	CAD Price
<input type="checkbox"/> Upholstery, Complete (Set)	One Size	880507-01	\$402
<input type="checkbox"/> Upholstery, Back Only (Ea.)	One Size	101204	\$197
<input type="checkbox"/> Upholstery, Seat Only (Ea.)	One Size	117005	\$221
<input type="checkbox"/> Seat Depth Reduction Reduces by 5cm	One Size	880508-06	\$73

Description	Size	Part No.	CAD Price
<input type="checkbox"/> Armrest, Swing-Away (Pr.)	2 x 14½ x 9¾"	880508-03	\$660
<input type="checkbox"/> Fix-Locks, White, Bag of 4 Pcs ** **Required for use of Vests	One Size	890505-36	\$58



Description	Size	Part No.	CAD Price
<input type="checkbox"/> Vest - Mesh	3	89111-340	\$273
<input type="checkbox"/> Vest - Mesh	4	89111-440	\$273
<input type="checkbox"/> Cross Vest - Mesh	3	89625-340	\$278
<input type="checkbox"/> Cross Vest - Mesh	4	89625-440	\$278
<input type="checkbox"/> Neoprene Hip-Belt	2	85423-2	\$132

Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Sides/Wings (Pr.)	4 x 14"	One Size	880508-01	\$117
--	---------	----------	-----------	-------




<input type="checkbox"/> Upholstery for Sides/Wings (Pr.)	4 x 14"	One Size	880507-04	\$117
---	---------	----------	-----------	-------



<input type="checkbox"/> Hip Supports (Pr.) <i>(seat width range: 13 1/2"-22")</i>	8 1/4"W x 4"H	One Size	880508-02	\$430
---	---------------	----------	-----------	-------

<input type="checkbox"/> Hip Support Spacer Blocks		One Size	880508-02LB	\$111
--	--	----------	-------------	-------

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------



<input type="checkbox"/> Swing-Away Bracket, Stainless, Right (Ea.)		One Size	83000-R	\$308
<input type="checkbox"/> Swing-Away Bracket, Stainless, Left (Ea.)		One Size	83000-L	\$308



<input type="checkbox"/> Bar, Straight, Stainless (Ea.)	6 1/4"	One Size	82648	\$0
<input type="checkbox"/> Bar, Offset, 1" Stainless, (Ea.)	4 3/4"	One Size	150238	\$0

<input type="checkbox"/> Cushion, Curved (Ea.)*	3 1/2 x 7 1/2"	One Size	6175407	\$142
<input type="checkbox"/> Cushion, Curved (Ea.)*	2 1/2 x 6 1/2"	One Size	6175408	\$142
<input type="checkbox"/> Cushion, Curved (Ea.)*	2 1/2 x 5 3/4"	One Size	6175409	\$142

***Must order two for set.**

<input type="checkbox"/> Fixed Side Support Pad	6 x 6"	One Size	880508-09	\$250
---	--------	----------	-----------	-------

<input type="checkbox"/> Protection Pads (Ea.)		One Size	89151	\$126
--	--	----------	-------	-------

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

<input type="checkbox"/> Headrest 2, Grey		One Size	6172000	\$280
<input type="checkbox"/> Headrest 3, Grey		One Size	6173000	\$280
<input type="checkbox"/> Headrest 4, Grey		One Size	6174000	\$280

Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable

<input type="checkbox"/> Stainless Turtle Bar System		One Size	98619-24	\$314
--	--	----------	----------	-------



Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

<input type="checkbox"/> Tray; Mounted in Armrest	17¼ x 7½"	One Size	880508-04	\$430
---	-----------	----------	-----------	-------



<input type="checkbox"/> Commode Pan	7¾ x 13½ x 4¾"	One Size	880130	\$70
--------------------------------------	----------------	----------	--------	------

<input type="checkbox"/> Lid for Commode Pan		One Size	880130-L	\$88
--	--	----------	----------	------



<input type="checkbox"/> Potty Ring		One Size	109072	\$65
-------------------------------------	--	----------	--------	------

<input type="checkbox"/> Potty Ring for Heron Power		One Size	128623	\$65
---	--	----------	--------	------



<input type="checkbox"/> Splash Guard	3 x 5¼"	One Size	880508-05	\$95
---------------------------------------	---------	----------	-----------	------



<input type="checkbox"/> Splash Guard, Low		One Size	880508-07	\$87
--	--	----------	-----------	------

Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Shock Absorber, White			One Size	116924-KIT	\$51
--	--	--	----------	------------	------

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in CAD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com

Etac
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5502-ENG-REV.06

Creating Possibilities