

R82 Heron

USD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | |
| Phone: | City/State/Zip: |
| Email | Phone: |


US Configuration

Effective: 01.21.25 - Rev 6



A (Back height) _____
 B (Arm height) _____
 C (Lower leg length) _____
 D (Seat depth) _____
 E (Chest width) _____
 F (Seat width) _____
 G (Height from seat to top of head) _____
 H (Shoulder width) _____
 Age _____
 Weight _____

Measurements




Heron with hydraulics

Heron with power

| | Heron with hydraulics | Heron with power |
|---|-------------------------|-------------------------|
| Effective seat width (A) | 550 mm / 21¾" | 550 mm / 21¾" |
| Seat width between hip supports (accessories) | 350-490 mm / 13¾ - 19¼" | 350-490 mm / 13¾ - 19¼" |
| Seat depth with reduction kit (accessories) | | |
| Backrest height (C) | 555-625 mm / 21¾-24½" | 555-625 mm / 21¾-24½" |
| Lower leg length (D) | 400-560 mm / 15¾-22" | 400-560 mm / 15¾-22" |
| Seat surface height at front edge (E) | 510-780 mm / 20-30¾" | 510-780 mm / 20-30¾" |
| Back width, top of backrest (F) | 290 mm / 11½" | 290 mm / 11½" |
| Overall width (G) | 670 mm / 26½" | 670 mm / 26½" |
| Length (H) | 710 mm / 28" | 710 mm / 28" |
| Overall height (I) | 1080-1350 mm / 44½-53¾" | 1080-1350 mm / 44½-53¾" |
| Seat plane angle | -5-25° | -5-25° |
| Back rest angle | 0-30° | 0-30° |
| Mass device | 28 kg / 61.75 lb | 31 kg / 68.25 lb |
| Max rated load | 130 kg / 286.6 lb | 163 kg / 359.3 lb |


Heron Standard Features



- One size
- 4 grey castors with locks
- Stainless steel frame
- Height adjustable frame
- Integrated push handles
- Built-in tilt in space
- Built-in recline
- Active, height adjustable back
- Potty ring

| Description | Size | Part No. | USD Price |
|--------------------------------------|----------|-------------|-----------|
| <input type="checkbox"/> Heron | One Size | 880505-01 | \$5,728 |
| <input type="checkbox"/> Heron Power | One Size | 880505-01EL | \$6,995 |

| Accessories | Description | Measurements | Size | Part No. | USD Price |
|-------------|-------------|--------------|------|----------|-----------|
|-------------|-------------|--------------|------|----------|-----------|



| | | | |
|---|----------|-----------|-------|
| <input type="checkbox"/> Foot Support, Right-Angled Knee (Pr.) Lower leg range: 12"-15" | Short | 880506-03 | \$540 |
| <input type="checkbox"/> Foot Support, Right-Angled Knee (Pr.) Lower leg range: 15.35"-21.65" | Long | 880506-01 | \$540 |
| <input type="checkbox"/> Foot Support, Open Knee Angle (Pr.) Lower leg range: 14"-21" | One Size | 880506-02 | \$540 |
| <input type="checkbox"/> Padding for Foot Supports (Pr.) | One Size | 880506-12 | \$72 |
| <input type="checkbox"/> Calf Support, Strap (Ea.) | 4" | 880509-01 | \$72 |
| <input type="checkbox"/> Heel Stops, Straps (Pr.) | 1½" | 880509-02 | \$131 |

| Description | Size | Part No. | USD Price |
|---|----------|-----------|-----------|
| <input type="checkbox"/> Upholstery, Complete (Set) | One Size | 880507-01 | \$322 |
| <input type="checkbox"/> Upholstery, Back Only (Ea.) | One Size | 101204 | \$160 |
| <input type="checkbox"/> Upholstery, Seat Only (Ea.) | One Size | 117005 | \$177 |
| <input type="checkbox"/> Seat Depth Reduction Reduces by 2" | One Size | 880508-06 | \$56 |

| Description | Size | Part No. | USD Price |
|--|---------------|-----------|-----------|
| <input type="checkbox"/> Armrest, Swing-Away (Pr.) | 2 x 14½ x 9¾" | 880508-03 | \$530 |
| <input type="checkbox"/> Fix-Locks, White, Bag of 4 Pcs ** **Required for use of Vests | One Size | 890505-36 | \$46 |



| Description | Size | Part No. | USD Price |
|--|------|-----------|-----------|
| <input type="checkbox"/> Vest - Mesh | 3 | 89111-340 | \$219 |
| <input type="checkbox"/> Vest - Mesh | 4 | 89111-440 | \$219 |
| <input type="checkbox"/> Cross Vest - Mesh | 3 | 89625-340 | \$214 |
| <input type="checkbox"/> Cross Vest - Mesh | 4 | 89625-440 | \$214 |
| <input type="checkbox"/> Neoprene Hip-Belt | 2 | 85423-2 | \$107 |

| Accessories | Description | Measurements | Size | Part No. | USD Price |
|-------------|-------------|--------------|------|----------|-----------|
|-------------|-------------|--------------|------|----------|-----------|



| | | | | |
|--|---------|----------|-----------|------|
| <input type="checkbox"/> Sides/Wings (Pr.) | 4 x 14" | One Size | 880508-01 | \$93 |
|--|---------|----------|-----------|------|




| | | | | |
|---|---------|----------|-----------|------|
| <input type="checkbox"/> Upholstery for Sides/Wings (Pr.) | 4 x 14" | One Size | 880507-04 | \$93 |
|---|---------|----------|-----------|------|



| | | | | |
|---|---------------|----------|-----------|-------|
| <input type="checkbox"/> Hip Supports (Pr.) <i>(seat width range: 13 1/2"-22")</i> | 8 1/4"W x 4"H | One Size | 880508-02 | \$344 |
|---|---------------|----------|-----------|-------|

| | | | | |
|--|--|----------|-------------|------|
| <input type="checkbox"/> Hip Support Spacer Blocks | | One Size | 880508-02LB | \$85 |
|--|--|----------|-------------|------|

| Description | Measurements | Size | Part No. | USD Price | HCPCS |
|-------------|--------------|------|----------|-----------|-------|
|-------------|--------------|------|----------|-----------|-------|



| | | | | | |
|---|--|----------|---------|-------|-------|
| <input type="checkbox"/> Swing-Away Bracket, Stainless, Right (Ea.) | | One Size | 83000-R | \$246 | E1028 |
| <input type="checkbox"/> Swing-Away Bracket, Stainless, Left (Ea.) | | One Size | 83000-L | \$246 | E1028 |



| | | | | | |
|---|--------|----------|--------|-----|-------|
| <input type="checkbox"/> Bar, Straight, Stainless (Ea.) | 6 1/4" | One Size | 82648 | \$0 | E1028 |
| <input type="checkbox"/> Bar, Offset, 1" Stainless, (Ea.) | 4 3/4" | One Size | 150238 | \$0 | E1028 |

| | | | | | |
|---|----------------|----------|---------|-------|-------|
| <input type="checkbox"/> Cushion, Curved (Ea.)* | 3 1/2 x 7 1/2" | One Size | 6175407 | \$116 | E0956 |
| <input type="checkbox"/> Cushion, Curved (Ea.)* | 2 1/2 x 6 1/2" | One Size | 6175408 | \$116 | E0956 |
| <input type="checkbox"/> Cushion, Curved (Ea.)* | 2 1/2 x 5 3/4" | One Size | 6175409 | \$116 | E0956 |

***Must order two for set.**

| | | | | |
|---|--------|----------|-----------|-------|
| <input type="checkbox"/> Fixed Side Support Pad | 6 x 6" | One Size | 880508-09 | \$250 |
|---|--------|----------|-----------|-------|

| | | | | |
|--|--|----------|-------|-------|
| <input type="checkbox"/> Protection Pads (Ea.) | | One Size | 89151 | \$100 |
|--|--|----------|-------|-------|

| Description | Measurements | Size | Part No. | USD Price | HCPCS |
|-------------|--------------|------|----------|-----------|-------|
|-------------|--------------|------|----------|-----------|-------|

| | | | | | |
|---|--|----------|---------|-------|-------|
| <input type="checkbox"/> Headrest 2, Grey | | One Size | 6172000 | \$224 | E0955 |
| <input type="checkbox"/> Headrest 3, Grey | | One Size | 6173000 | \$224 | E0955 |
| <input type="checkbox"/> Headrest 4, Grey | | One Size | 6174000 | \$224 | E0955 |


Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable

| | | | | | |
|--|--|----------|----------|-------|-------|
| <input type="checkbox"/> Stainless Turtle Bar System | | One Size | 98619-24 | \$251 | E0955 |
|--|--|----------|----------|-------|-------|



| Description | Measurements | Size | Part No. | USD Price |
|---|----------------|----------|-----------|-----------|
| <input type="checkbox"/> Tray; Mounted in Armrest | 17¼ x 7½" | One Size | 880508-04 | \$344 |
| <input type="checkbox"/> Commode Pan | 7¾ x 13½ x 4¾" | One Size | 880130 | \$56 |
| <input type="checkbox"/> Lid for Commode Pan | | One Size | 880130-L | \$72 |
| <input type="checkbox"/> Potty Ring | | One Size | 109072 | \$50 |
| <input type="checkbox"/> Potty Ring for Heron Power | | One Size | 128623 | \$50 |
| <input type="checkbox"/> Splash Guard | 3 x 5¼" | One Size | 880508-05 | \$77 |
| <input type="checkbox"/> Splash Guard, Low | | One Size | 880508-07 | \$71 |



| Accessories | Description | Measurements | Size | Part No. | USD Price |
|---|--|--------------|----------|------------|-----------|
|  | <input type="checkbox"/> Shock Absorber, White | | One Size | 116924-KIT | \$42 |

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com

Etac
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5500-ENG-REV.06

Creating Possibilities