



R82 Meerkat - E0638

AADL CAD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | |
| Phone: | City/State/Zip: |
| Email | Phone: |


US Configuration

Effective: 01.01.22 - Rev 8

Meerkat (Baseplate and Column)



Accessories

- Chest Support, Adjustable
- Pelvic Support, Adjustable
- Knee Supports
- Sandals
- Rocker
- Tray



- A (Total Height) _____
- B (Armpit Height) _____
- C (Top of Pelvic Height) _____
- D (Chest Circumference) _____
- E (Pelvic Circumference) _____
- F (Knee Circumference) _____
- G (Foot Length) _____
- H (Foot Width) _____
- Age _____
- Weight _____

Measurements



Product Specs

| | Size 1 | Size 2 | Size 3 |
|--------------------|--|--|--|
| Base Width (A) | 20 ³ / ₄ " | 30 ³ / ₄ " | 30 ³ / ₄ " |
| Base Length (B) | 25 ¹ / ₄ " | 33 ³ / ₄ " | 33 ³ / ₄ " |
| Product Height (C) | 14 ³ / ₄ -29 ¹ / ₂ " | 19 ³ / ₄ -39 ¹ / ₂ " | 25 ¹ / ₂ -51 ¹ / ₄ " |
| Product Weight | 19 ³ / ₄ lbs. | 33 ¹ / ₂ lbs. | 34 ³ / ₄ lbs. |

User Measurement Guide for Stander

| | Size 1 | Size 2 | Size 3 |
|-------------------------------|--|--|--|
| Approx. Age | 1-4 Years | 4-10 Years | 6-18 Years |
| Max. User Weight/Load | 44 lbs. | 143 ¹ / ₄ lbs. | 198 ¹ / ₂ lbs. |
| Approx. Total User Height (A) | 27 ¹ / ₂ -41 ¹ / ₄ " | 37 ¹ / ₂ -53 ¹ / ₄ " | 46-70 ³ / ₄ " |
| Max Armpit Height (B) | 15 ³ / ₄ -29 ¹ / ₂ " | 28 ¹ / ₄ -39 ¹ / ₄ " | 33-51 ¹ / ₄ " |
| Approx. Top of Pelvis (C) | 13 ³ / ₄ -21 ³ / ₄ " | 19 ³ / ₄ -31 ¹ / ₂ " | 25 ¹ / ₂ -41 ¹ / ₄ " |

User Measurement Guide for Accessories

| | Small | Medium | Large | XL |
|--------------------------|--|--|--|--|
| Chest Circumference (D) | 17 ³ / ₄ -23 ¹ / ₂ " | 23 ¹ / ₂ -29 ¹ / ₄ " | 27 ¹ / ₂ -37 ¹ / ₂ " | 37 ¹ / ₂ - 43 1/4" |
| Pelvic Circumference (E) | 17 ³ / ₄ -23 ¹ / ₂ " | 23 ¹ / ₂ -29 ¹ / ₄ " | 27 ¹ / ₂ -37 ¹ / ₂ " | 37 ¹ / ₂ - 43 1/4" |
| Knee Circumference (F) | 7 ¹ / ₂ -10" | 10-11 ¹ / ₂ " | 11 ¹ / ₂ -15 ¹ / ₄ " | |
| Knee Width | 2 ¹ / ₄ " | 3 ¹ / ₄ " | 4" | |
| Heel Width (G) | 2 ¹ / ₄ -3" | 2 ¹ / ₂ -3 ¹ / ₄ " | 3 ¹ / ₄ -3 ³ / ₄ " | |
| Forefoot Width (H) | 3 ³ / ₄ -4 ¹ / ₄ " | 4-5 ³ / ₄ " | 4 ³ / ₄ -6 ¹ / ₂ " | |


The Meerkat is measured in standard set-up and with the center column in vertical position. Different settings may result in another set of measurements. To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without.

NOTES:

Meerkat Standard Features














- Base Plate with Wheels, Column and Tool
- Rocker + Handle




| Description | Size | Part No. | CAD Price | HCPCS |
|---|------|-----------|-----------|-------|
| <input type="checkbox"/> Meerkat, Pastel Blue, size 1 Plate, size 1 Column | 1 | 870011-1R | \$1,261 | E0638 |
| <input type="checkbox"/> Meerkat, Strawberry Pink, size 1 Plate, size 1 Column ●● | 1 | 870011-2R | \$1,261 | E0638 |
| <input type="checkbox"/> Meerkat, Grey, size 2 Plate, size 2 Column | 2 | 870022-R | \$1,448 | E0638 |
| <input type="checkbox"/> Meerkat, Grey, size 2 Plate, Grey, size 3 Column | 3 | 870023-R | \$1,521 | E0638 |
| ●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item | | | | |

| Accessories | Description | Measurements | Size | Part No. | CAD Price | HCPCS |
|-------------|---|--|--------|----------|-----------|-------|
| | <input type="checkbox"/> Pelvic Support | 5½" – 6½" (width) 17¾ - 23½" (circumference) | Small | 869240 | \$527 | |
| | <input type="checkbox"/> Pelvic Support | 6½" – 8½" (width) 23½ - 29" (circumference) | Medium | 869241 | \$532 | |
| | <input type="checkbox"/> Pelvic Support | 8½" – 12"(width) 27½ - 37½" (circumference) | Large | 869242 | \$550 | |
| | <input type="checkbox"/> Pelvic Support | 12" – 14¼ (width) 37½ - 43 ¼" (circumference) | XL | 869243 | \$562 | |
| | <input type="checkbox"/> Chest Support | 5½" – 6½" (width) 17¾ - 23½" (circumference) | Small | 869480 | \$527 | |
| | <input type="checkbox"/> Chest Support | 6½" – 8½" (width) 23½ - 29" (circumference) | Medium | 869481 | \$532 | |
| | <input type="checkbox"/> Chest Support | 8½" – 12"(width) 27½ - 37½" (circumference) | Large | 869482 | \$550 | |
| | <input type="checkbox"/> Chest Support | 12" – 14¼ (width) 37½ - 43 ¼" (circumference) | XL | 869483 | \$562 | |

| Accessories | Description | Measurements | Size | Part No. | CAD Price | HCPCS |
|--|---|---|----------|-------------|-----------|------------|
|  | <input type="checkbox"/> Knee Supports | 7½ - 9¾" | Small | 869840 | \$550 | |
| | <input type="checkbox"/> Knee Supports | 9¾ - 11½" | Medium | 869841 | \$564 | |
| | <input type="checkbox"/> Knee Supports | 11½ - 15¼" | Large | 869842 | \$564 | |
|  | <input type="checkbox"/> Sandals for size 1 Plate | 2¼ - 3" (Heel Width) 3¾ - 4¼" (Forefoot Width) | Small | 869490 | \$277 | |
| | <input type="checkbox"/> Sandals for size 2 Plate | 2½ - 3¼" (Heel Width) 4 - 5¾" (Forefoot Width) | Medium | 869491 | \$290 | |
| | <input type="checkbox"/> Sandals for size 2 Plate | 3¼ - 3¾" (Heel Width) | Large | 869492 | \$299 | |
| <hr/> | | | | | | |
| | Description | Measurements | Size | Part No. | CAD Price | HCPCS |
|  | <input type="checkbox"/> Tray (Only Anterior Use) | | One Size | 869485 | \$416 | E1399 |
|  | <input type="checkbox"/> Displacement kit bracket for mounting on column for pelvic and chest support | | One Size | 124025 | \$94 | <u>Qty</u> |
|  | <input type="checkbox"/> Daily adjustment kit applicable on knee support <i>*For width adjustment on knee support</i> | | One Size | 124026 | \$48 | |
|  | <input type="checkbox"/> Daily adjustment kit, height ** <i>**For height adjustment of chest, pelvic or knee support</i> | | One Size | 124027 | \$37 | |
|  | <input type="checkbox"/> Padding for Meerkat Ratchet Straps (Pr.) | | One Size | 117818-PADS | \$65 | |

| Spare Parts | Description | Measurements | Size | Part No. | CAD Price | HCPCS |
|---|--|--------------|--------------|----------|-----------|------------|
|  | <input type="checkbox"/> Chest Belt | | Small | 869240-1 | \$326 | |
| | <input type="checkbox"/> Chest Belt | | Medium | 869241-1 | \$326 | |
| | <input type="checkbox"/> Chest Belt | | Large | 869242-1 | \$326 | |
| <hr/> | | | | | | |
|  | <input type="checkbox"/> Mounting for Ski Strap Buckle, Chest and Pelvic Support | | One Size | 117818 | \$60 | <u>Qty</u> |
| | <input type="checkbox"/> Mounting for Ski Strap Buckle, Knee Support | | One Size | 117818-1 | \$60 | |
| <hr/> | | | | | | |
|  | <input type="checkbox"/> Knee Support, Cup | | Small | 135310 | \$136 | <u>Qty</u> |
| | <input type="checkbox"/> Knee Support, Cup, Right | | Medium | 135311 | \$136 | |
| | <input type="checkbox"/> Knee Support, Cup, Left | | Medium | 135312 | \$136 | |
| | <input type="checkbox"/> Knee Support, Cup, Right | | Large | 135313 | \$136 | |
| | <input type="checkbox"/> Knee Support, Cup, Left | | Large | 135314 | \$136 | |
| <hr/> | | | | | | |
|  | <input type="checkbox"/> Knee Support, Strap and Pad | | Small | 128254-1 | \$150 | <u>Qty</u> |
| | <input type="checkbox"/> Knee Support, Strap and Pad | | Medium | 128255-1 | \$150 | |
| | <input type="checkbox"/> Knee Support, Strap and Pad | | Large | 128256-1 | \$150 | |
| <hr/> | | | | | | |
| | <input type="checkbox"/> Sandal Strap | | Small,Medium | 127085 | \$56 | <u>Qty</u> |
| | Sandal Strap | | Large | 127086 | \$56 | |

| Spare Parts | Description | Measurements | Size | Part No. | CAD Price | HCPCS |
|---|---|--------------|----------|----------|-----------|-------------------|
| | <input type="checkbox"/> Column Mount Complete | | One size | 117832 | \$127 | |
|  | <input type="checkbox"/> Rods and Screws for Rocker | | 1 | 118598-2 | \$163 | |
| | <input type="checkbox"/> Rods and Screws for Rocker | | 2 | 117047-2 | \$163 | |
| | <input type="checkbox"/> Wheels (Airplane) | | 1 | 3035000 | \$59 | |
| | <input type="checkbox"/> Wheels (Cow) | | 1 | 3035010 | \$59 | <u>Qty</u> |
| | <input type="checkbox"/> Wheels (Grey) | | 2 | 3006825 | \$79 | <u> </u> |

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in CAD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5702-ENG-REV.08

