



# R82 Mustang - E8000/E8002

## AADL CAD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration



Effective: 01.01.22 - Rev 12

- Mustang
- Accessories**
- Chest Support
- Head Support
- Back Support
- Arm Support
- Seat/Sling Seat
- Leg Separator
- Ankle Prompts



- A** (Total Height) \_\_\_\_\_
- B** (Shoulder Height) \_\_\_\_\_
- C** (Chest Height) \_\_\_\_\_
- D** (Buttock Height) \_\_\_\_\_
- E** (Chest Width) \_\_\_\_\_
- F** (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_


### Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches	Size 4 Inches
Width, total (G)	20¼"	22½"	25½"	27½"
Width, inside frame (G1)	14¾"	17½"	19½"	22"
Length with 4" wheels (H)	29¼"	32½"	-	-
Length with 8" wheels (H)	32"	35"	39"	43"
Height with 4" wheels (I)	14¾-25¼"	23½-33½"	-	-
Height with 8" wheels (I)	19-25¼"	23½-33½"	27½-41"	34-50¾"
Seat*) height above floor, 4" wheels	6½-13½"	9½-16"	-	-
Seat*) height above floor, 8" wheels	10-17¼"	12¾-19½"	17½-26½"	24-32¾"
Width between handles *)	5½-13½"	5½-13½"	5½-23½"	15½-23½"
Angle range of center bar (a1)	30°	30°	0-25°	0-25°
Turning radius	18	21½	23½	25¼"
User height, floor to axil	14¾-25¼"	23½-33½"	27½-41½"	34-50¾"
Weight **) with 4" wheels	22 lbs.	24 lbs.	-	-
Weight **) with 8" wheels	28½ lbs.	31 lbs.	40 lbs.	42 lbs.
Max. load / user weight	66 lbs.	88 lbs.	132 lbs.	176 lbs.
Approx. min. user height (4" wheels)	26½"	33"	-	-
Approx. min. user height (8" wheels)	28¾"	34½"	40½"	49½"
Approx. max. user height (4" wheels)	36½"	44½"	-	-
Approx. max. user height (8" wheels)	36½"	44½"	54½"	66"


*The Mustang is measured as standard delivered and with the center bar in vertical position. Different settings may result in another set of measurements. To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without.*

## Mustang Standard Features



- Four Sizes with Four Solid Wheels
- Drag Wheel Function
- Reverse Brake On/Off
- Directional Stabilizer
- Angle Adjustment
- Height Adjustment
- Chest Support, Height/Width Adjustable

Description	Size	Part No.	CAD Price	HPCPS
<input type="checkbox"/> Mustang, Yellow, 4" Casters	1	869011	\$2,516	E8002
<input type="checkbox"/> Mustang, Yellow, 8" Casters	1	869015	\$2,516	E8002
<input type="checkbox"/> Mustang, Yellow, 8" Casters Slim Line Wheels	1	869015-SL	\$2,516	E8002
<input type="checkbox"/> Mustang, White, 4" Casters	2	869025-4	\$2,621	E8002
<input type="checkbox"/> Mustang, White, 8" Casters	2	869024	\$2,621	E8002
<input type="checkbox"/> Mustang, White, 8" Casters Slim Line Wheels	2	869024-SL	\$2,621	E8002
<input type="checkbox"/> Mustang, Red, 8" Casters	3	869033	\$3,707	E8002
<input type="checkbox"/> Mustang, Red, 8" Casters Slim Line Wheels	3	869033-SL	\$3,707	E8002
<input type="checkbox"/> Mustang, Anthracite (Grey), 8" Casters	4	869042	\$4,240	E8002
<input type="checkbox"/> Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels	4	869042-SL	\$4,240	E8002



Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

Required, must select one:



	Inner Circum.	Width Range	Size	Part No.	CAD Price
<input type="checkbox"/> Chest Support	15½-25¼"	6 - 8½"	Small	869810	N/C
<input type="checkbox"/> Chest Support	26¼ -30¾"	9½ - 12"	Medium	869820	N/C
<input type="checkbox"/> Chest Support	34¼ -44¾"	11 - 14¼"	Large	869830	N/C


Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------



<input type="checkbox"/> Cover for Chest Support	H:3¼ x W:9¾"	Small	869810-160	\$150
<input type="checkbox"/> Cover for Chest Support	H:4¾ x W:14¾"	Medium	869820-260	\$182
<input type="checkbox"/> Cover for Chest Support ●●	H:5¾ x W:19"	Large	869830-360	\$195

●●Custom order; Approx. lead time: 6-8 weeks (sz Large)

Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------




	Inner Circum.	Width Range	Size	Part No.	CAD Price
<input type="checkbox"/> Pelvic Support	15½-25¼"	6 - 8½"	Small	869610	\$618
<input type="checkbox"/> Pelvic Support	26¼ -30¾"	9½ - 12"	Medium	869620	\$637
<input type="checkbox"/> Pelvic Support ●●	34¼ -44¾"	11 - 14¼"	Large	869630	\$657

●●Custom order; Approx. lead time: 6-8 weeks (sz Large)




<input type="checkbox"/> Extra Pad for Chest/Pelvic Sup.	H:4¾ x W:2¼"	One Size	869225	\$79
--	--------------	----------	--------	------

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------














<input type="checkbox"/> Sling Seat for Small Chest Support	W:7 x L:9½"	1	869570	\$228
<input type="checkbox"/> Sling Seat for Med Chest Support	W:7 x L:12¾"	2	869571	\$228

Note: Sling Seat n/a on size's 3 and 4 Mustang



<input type="checkbox"/> Mustang Seat	H:6¾ x D:10"	1, 2	869110	\$332
<input type="checkbox"/> Mustang Seat	H:7"x D:12¾"	3, 4	869120	\$351
<input type="checkbox"/> Mustang Comfort Seat	L:8½"	1, 2	869115	\$482
<input type="checkbox"/> Mustang Comfort Seat	L:11¾"	3, 4	869125	\$501

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Hip Support on Curved Bar ♦		1, 2	869222	\$501
	<input type="checkbox"/> Hip Support on Curved Bar ♦		3, 4	869232	\$501
	♦Requires Mustang Seat or Mustang Comfort Seat				
	<input type="checkbox"/> Hip Sup. On Straight Bar ♦	W:7¼"-10"	1, 2	869221	\$501
	<input type="checkbox"/> Hip Sup. On Straight Bar ♦	W:10-12¾"	3, 4	869231	\$501
	♦Requires Mustang Seat or Mustang Comfort Seat				
	<input type="checkbox"/> Sacrum Support ♦♦	H:7", D:5-8¼"	Small	869250	\$209
	<input type="checkbox"/> Sacrum Support ♦♦	H:8", D:6½"-9½"	Large	869255	\$286
	♦Requires Mustang Seat or Mustang Comfort Seat				
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Large)				
	Note: Measurements from Seat				
	<input type="checkbox"/> Abdominal Support	H:5 x W:7½"	One Size	869260	\$120
<hr/>					
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>CAD Price</b>
	<input type="checkbox"/> Hand Bar for Anterior Use	H:41/4 x W:8¼"	Small	869360	\$273
	<input type="checkbox"/> Hand Bar for Anterior Use ♦♦	H:41/4 x W:5½"	Large	869365	\$273
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Large)				
	<input type="checkbox"/> Armrest for Anterior Use w/ Handles		Small	869350	\$533
	<input type="checkbox"/> Armrest for Anterior Use w/ Handles		Large	869355	\$657
	<input type="checkbox"/> Arm Prompts for Anterior Use		Small	869310	\$501
	<input type="checkbox"/> Arm Prompts for Anterior Use		Medium	869320	\$501
	<input type="checkbox"/> Arm Prompts for Anterior Use		Large	869325	\$527
	<input type="checkbox"/> Arm Prompts for Posterior Use		Small	869330	\$527
	<input type="checkbox"/> Arm Prompts for Posterior Use ♦♦		Medium	869340	\$527
	<input type="checkbox"/> Arm Prompts for Posterior Use ♦♦		Large	869345	\$527
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Medium and Large)				
	<input type="checkbox"/> Padded Forearm Positioner (Ea.)		One Size	869370	\$89
<hr/>					
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>CAD Price</b>
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$299
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$299
	♦Required with Anatomic Head Support				
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C
	♦Turtle Bar System Short or Long Required; 98619-25 or 98619-26				
	<input type="checkbox"/> Headband System ♦		One Size	99500	\$299
	♦For Anatomic Head Support Only				

QTY



- Connecting Rod Between Hip and Head Support ♦ One Size 869226 \$73  
♦ **Required to connect Headrest to Hip or Sacrum Supports**

Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



- Leg Separator for sz 1 w/ 4" Wheels 1 869410 \$377
- Leg Separator for sz 1 w/ 8" Wheels 1 869410-8 \$377
- Leg Separator for sz 2 w/ 4" Wheels 2 869430 \$377
- Leg Separator for sz 2 w/ 8" Wheels 2 869420 \$377
- Leg Separator for sz 3 3 869430 \$377
- Leg Separator for sz 4 4 869440 \$455



- Ankle Prompts for sz 1 ♦ 1 869450 \$228
- Ankle Prompts for sz 2 w/ 8" Wheels 2 869455 \$228
- Ankle Prompts for sz 2 w/ 4" Wheels 2 869455-4 \$228
- Ankle Prompts for sz 3 3 869460 \$247
- Ankle Prompts for sz 4 4 869465 \$247

♦ **Can be used on Mustang sz 1 with both 4" and 8" Wheels**



- Leg Pad One Size 869470 \$176

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------



- Attendant Guide Bar L: 27¼" 1, 2 869510 \$260  
**Note: N/A on Mustang sz 3 and 4**



- Storage Pouch •• H: 5¾ x W: 5" 1, 2 869560 \$126  
•• **Custom order; Approx. lead time: 6-8 weeks**



- Storage Bag •• H: 12 x W: 15½" 3, 4 869565 \$253  
•• **Custom order; Approx. lead time: 6-8 weeks**

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------



- 4" Wheel and Holder, Left Side 1 79436-L-KIT \$370
- 4" Wheel and Holder, Right Side 1 79436-R-KIT \$370
- 8" Wheel and Holder, Left Side 2, 3, 4 79438-L-KIT \$533
- 8" Wheel and Holder, Right Side 2, 3, 4 79438-R-KIT \$533
- 8" Slim-Line Wheel and Holder, Left Side 2, 3, 4 79439-L-KIT \$530
- 8" Slim-Line Wheel and Holder, Right Side 2, 3, 4 79439-R-KIT \$530
- 8" Slim-Line Wheels and Bushings, Kit of 4 2, 3, 4 79439-SLKIT \$530

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Prices in CAD & is subject to change without notice**

**For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537**

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
www.etac.us.com

OF5802-ENG-REV.12

