

# R82 Quokka™

## CAD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_


Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	Phone:

### US Configuration


Effective: 01.01.22 Rev 5

Quokka Child Restraint System




A (Back Height) \_\_\_\_\_  
 B (Arm Height) \_\_\_\_\_  
 C (Lower Leg Length) \_\_\_\_\_  
 D (Seat Depth) \_\_\_\_\_  
 E (Chest Width) \_\_\_\_\_  
 F (Seat Width) \_\_\_\_\_  
 G (Height: Seat to Top of Head) \_\_\_\_\_  
 H (Shoulder Width) \_\_\_\_\_  
 Age \_\_\_\_\_  
 Weight \_\_\_\_\_

### Measurements



	<b>Quokka</b>
	<b>Inches</b>
Seat Width Maximum (A)	11"
Seat Depth (B)	11"
Height of Shoulder Strap Slots (C)	7 ½", 9 ½", 11 ½", 13 ½"
Back Height (Top of Head) Maximum (D)	21"
Head Support Height Range Minimum	13 ½"
Head Support Height Range Maximum	19 ½"
Head Support Height	6 ½"
Exterior Shell Width (W), Depth (D), Height (H)	16" x 23" x 23 ½"
Interior Shell Width (W), Depth (D), Height (H)	11" x 11" x 21"
Minimum User Weight	5 lbs.
Maximum User Weight	39.6 lbs. (including cast weight)
Maximum User Height	41"
Product Weight	15.8 lbs.

### Quokka Standard Features



- Quokka Seat, One Size, Black Shell
- Head Support (Red)
- 5-Point Safety Harness
- Harness Padded Covers
- Seat Cover (Red/Black Trim)
- Spica Positioning Wedge
- LATCH System (Lower Anchors and Tether Strap)

Description	Size	Part No.	CAD Price
<input type="checkbox"/> Quokka Seat, (Red/Black Trim)	One Size	905371-2	\$911

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Seat Saver, Incontinence Pad		One Size	907029	\$127
Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Harness Buckle & Strap * <i>*Replacement for standard part</i>	6"	Standard	905386	\$26
	<input type="checkbox"/> Harness Buckle & Strap	10"	XL	905387	\$26
	<input type="checkbox"/> Replacement Seat Cover		One Size	905376-2	\$215
	<input type="checkbox"/> Replacement 5-Point Harness		One Size	905378	\$111
	<input type="checkbox"/> Extra Tether Strap		One Size	905375	\$60
	<input type="checkbox"/> Extra Lower Anchor Strap		One Size	905374	\$60
	<input type="checkbox"/> Spica Wedge Positioning System		One Size	905373	\$39
	<input type="checkbox"/> Replacement Headrest		One Size	905377-2	\$65
	<input type="checkbox"/> Replacement Shoulder Protectors 7"			905383	\$56
	<input type="checkbox"/> Replacement Shoulder Protectors 4"			905384	\$56
	<input type="checkbox"/> Belt Locking Clip Set		Pair	905382	\$19

**IMPORTANT: Due to safety regulations and restrictions, no returns or exchanges are accepted of car seat products.**

**Prices in CAD & is subject to change without notice**

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or Fax: 1-800-760-0537**

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF6002-ENG-REV.05

