

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:



Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Rabbit Up



A (Total Height) _____
B (Armpit Height) _____
C (Top of Pelvic Height) _____
D (Buttock Height) _____
E (Chest Width) _____
F (Pelvic Width) _____
Age _____
Weight _____

Features, incl. in standard price


Stander Frame
Centre Bar
Brake on Castors
Angling Spindle




Mobile Stander Frame
Centre Bar
Driving Wheels
Angling Spindle
Closed Pushrim
0° Camber std. Angle
Spoke protectors
Central brake
Brake on Casters




Measurements




Size 1



Size 2



Size 3



Size 4

Product Specs as a standing frame

Overall width (G)	20¼"	21¾"	23¾"	23¾"
Length (H)	27¾"	31¼"	31¼"	36½"
Total mass	19.8 lbs.	23.1 lbs.	24.2 lbs.	29.3 lbs.

Product Specs as a mobile stander

	w/26" wheels	w/32" wheels	w/32" wheels	w/36" wheels
Overall width (G)	25"	26¾"	28¾"	28¼"
Length (H)	24½"	29½"	29½"	35½"
Total mass	29.75 lbs.	34.1 lbs.	35.25 lbs.	40.75 lbs.

	Size 1	Size 2	Size 3	Size 4
Overall measurements				
Width, inside (G1)	12¾"	14¾"	16¾"	16¾"
Overall height (I) *	26 ¼" – 39 ¼"	33 ¾" – 46 ¾"	39" – 51 ½"	47 ¼" – 60 ¼"
Height from foot plate (I1) *	21 ¾" – 29 ¼"	29 ½" – 42 ½"	32 ½" – 46 ¾"	37 ¾" – 56"
Center bar angle (a1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Foot plate angle (a2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	6°	6°	4°
Center bar length	19¾"	27½"	31½"	39¼"

User Measurement Guide for Stander

Recommended user height	30" - 39"	39" - 54"	43" - 58"	48" - 71"
Max armpit height	29"	42"	46"	56"
User weight max	88.1 lbs.	110.2 lbs.	132.2 lbs.	154.3 lbs.
Recommended user age	1 - 3 yr./old	4 - 10 yr./old	5 - 12 yr./old	7 - 18 yr./old

* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

Rabbit Stander Options

Make a selection of either the Prone Stander or the Mobile Prone Stander

Rabbit-Up Prone Stander



Rabbit-Up, Stander, Yellow, Size 1 (E0638) (8621220-00)	\$1868
Rabbit-Up, Stander, Blue, Size 2 (E0638) (8622220-00)	\$1915
Rabbit-Up, Stander, Grey, Size 3 (E0638) (8623220-00)	\$1962
Rabbit-Up, Stander, Silver, Size 4 (E0638) (8624220-00)	\$2106

Rabbit Mobile Prone Stander



22" Wheels, Yellow, Size 1 (E0642) (8621222-21)	\$3582
24" Wheels, Yellow, Size 1 (E0642) (8621223-21)	\$3582
26" Wheels, Yellow, Size 1 (E0642) (8621224-21)	\$3582
29" Wheels, Blue, Size 2 (E0642) (8622225-21)	\$4486
32" Wheels, Blue, Size 2 (E0642) (8622226-21)	\$4486
29" Wheels, Grey, Size 3 (E0642) (8623225-21)	\$4534
32" Wheels, Grey, Size 3 (E0642) (8623226-21)	\$4534
32" Wheels, Silver, Size 4 (E0642) (8624226-21)	\$4677
36" Wheels, Silver, Size 4 (E0642) (8624227-21)	\$4819

Chest Supports

Supports and Pad



Chest Support and Pad

Chest Support and Pad, W: 7.75-10.75", Size 1 (861704)	\$280
Chest Support and Pad, W: 7.75-10.75", Size 2 and 3 (861702)	\$280
Chest Support and Pad, W: 9.75-12.5", Size 4 (861703)	\$308



Filler Pad

Not required for Surrounding Chest Support

Filler Pad, Small (861421)	\$113
Filler Pad, Large (861422)	\$113



Extension for Center Bar

Available with Chest and Surrounding Chest Supports

Extension for Center Bar (861480)	\$113
-----------------------------------	-------



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Rabbit Up



Surrounding Chest Support and Options



Surrounding Chest Support

Surrounding Chest Support, C:18.5-25.25", Small (861810)	\$667
Surrounding Chest Support, C: 26.5-30.75", Medium (861820)	\$694
Surrounding Chest Support, C: 34.25-44.75", Large (861830)	\$735



Cover for Surrounding Chest Support

Cover for Surrounding Chest Support, C:18.5-25.25", Small (869810-160)	\$164
Cover for Surrounding Chest Support, C: 26.5-30.75", Medium (869820-260)	\$200
Cover for Surrounding Chest Support, C: 34.25-44.75", Large (869830-360)	\$215

Extra Pad for Surrounding Support

Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225)	\$83
---	------



Trays



Trays

Tray, Plastic, Grey (862515)	\$525
Tray With Storage Box (862501)	\$538
Plastic Tray with Bowl and Cover (862525)	\$765

Mounting Hardware required for Tray with Bowl

Tray Accessories

Height Adjustment of Tray (862510)	\$308
------------------------------------	-------

Required with Tray with Storage Box

Mounting Hardware

Required Selection for Tray with Bowl

Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit)	\$378
--	-------

Abdominal Support



Abdominal Supports

Abdominal Support (Each) (861411)	\$219
-----------------------------------	-------

Head Supports

Head Support, Head Support Cover, Connecting System for Head Support, Headrest Mounting Hardware, and Hip Support w/ Pommel are required for a complete Head Support

Head Support



Anatomic Head Support, Small (5099001)	\$348
Anatomic Head Support, Medium (5099002)	\$348
Anatomic Head Support, Large (5099003)	\$348

Connecting System for Head Support




Connecting System for Head Support, Size 1 and 2 (861424-3)	\$352
Connecting System for Head Support, Size 3 and 4 (861424-4)	\$352

Head Support Cover


Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium (9945-83)	N/C
Anatomic Head Support Cover, Large (9967-83)	N/C

Headrest Mounting Hardware



Headrest Mounting Hardware, Short (98619-25)	\$350
Headrest Mounting Hardware, Long (98619-26)	\$358

Pelvic and Hip Control Options



Hip Support w/ Pommel

Front mount to center bar. Required with Head Support.

Hip Support w/ Pommel, Size 1 (861501)	\$688
Hip Support w/ Pommel, Size 2 and 3 (861502)	\$688
Hip Support w/ Pommel, Size 4 (861503)	\$763



Hip Support, With Rotation

Hip Support, With Rotation, Size 2 and 3 (861612)	\$484
Hip Support, With Rotation, Size 4 (861613)	\$525



Hip Support, Front Mount


Hip Support, Front Mount, Size 1 (861601)	\$320
Hip Support, Front Mount, Size 2 and 3 (861602)	\$257
Hip Support, Front Mount, Size 4 (861603)	\$362



Hip Support, Surrounding

Front mount to center bar

Hip Support, Surrounding, Size 1 (861610)	\$754
Hip Support, Surrounding, Size 2 and 3 (861620)	\$754
Hip Support, Surrounding, Size 4 (861630)	\$770



Pommel With Back Support


Pommel With Back Support (861401)	\$618
-----------------------------------	-------

Knee Supports



Knee Support

Knee Supports, Size 1 (Pair) (861301)	\$381
Knee Supports, Size 2 and 3 (Pair) (861302)	\$484
Knee Supports, Size 4 (Pair) (861303)	\$484




Hinged Knee Supports

Hinged Knee Supports, Size 1 (Pair) (861311)	\$593
Hinged Knee Supports, Size 2 and 3 (Pair) (861312)	\$593
Hinged Knee Supports, Size 4 (Pair) (861313)	\$593

Footplate and Sandals Options

Flip-up Footplate required for all Sandals, except the Individual Footrest w/ Hinged Sandals

Footplates and Options




Flip-Up Footplate

Flip-Up Footplate, Size 1 and 2 (861211)	\$299
Flip-Up Footplate, Size 2, 3 and 4 (861212)	\$299



Formed Footplates

Formed Footplates, Size 0 (Pair) (86149)	\$226
Formed Footplates, Size 1 and 2 (Pair) (86147)	\$226
Formed Footplates, Size 3 (Pair) (86148)	\$226



Footplate Straps

Foot Straps, Plastic Lock (Pair) (85096)	\$59
Footplate Straps, Velcro (Pair) (85102)	\$70



Individual Footrest w/ Hinged Sandals

Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201)	\$538
Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202)	\$538

Sandals and Options



Sandals

Sandals for Rabbit Up/ Toucan, Small (869490-RT)	\$326
Sandals for Rabbit Up/ Toucan, Medium (869491-RT)	\$342
Sandals for Rabbit Up/ Toucan, Large (869492-RT)	\$351



Heel Stops

Heel Stops, Small (Pair) (861251)	\$85
Heel Stops, Large (Pair) (861252)	\$99



Sandals

Hinged Sandals w/ Straps, Size 1 (Pair) (861221)	\$185
Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222)	\$185
Hinged Sandals w/ Straps, Size 4 (Pair) (861223)	\$297

Conversion Kits

Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit. Conversion kit is used to convert Prone Stander into Mobile Prone Stander.

Driving Wheel Conversion Kits



Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209)	\$788
Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210)	\$788
Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211)	\$788
Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212)	\$1214
Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213)	\$1214
Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214)	\$1214



Central Brake System

Central Brake System, Size 1 (134215)	\$395
Central Brake System, Size 2 and 3 (134216)	\$395
Central Brake System, Size 4 (134218)	\$395



Camber Kit Rabbit Up

Camber Kit Rabbit Up 0° (134200)	\$71
Camber Kit Rabbit Up +3/-3° (134201)	\$71
Camber Kit Rabbit Up +6/-6° (134202)	\$71