



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Rabbit Up

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:



Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Rabbit Up



A (Total Height) _____
B (Armpit Height) _____
C (Top of Pelvic Height) _____
D (Buttock Height) _____
E (Chest Width) _____
F (Pelvic Width) _____
Age _____
Weight _____

Features, incl. in standard price


Stander Frame
Centre Bar
Brake on Castors
Angling Spindle




Mobile Stander Frame
Centre Bar
Driving Wheels
Angling Spindle
Closed Pushrim
0° Camber std. Angle
Spoke protectors
Central brake
Brake on Casters




Measurements




Size 1



Size 2



Size 3



Size 4

Product Specs as a standing frame

Overall width (G)	20¼"	21¾"	23¾"	23¾"
Length (H)	27¾"	31¼"	31¼"	36½"
Total mass	19.8 lbs.	23.1 lbs.	24.2 lbs.	29.3 lbs.

Product Specs as a mobile stander

	w/26" wheels	w/32" wheels	w/32" wheels	w/36" wheels
Overall width (G)	25"	26¾"	28¾"	28¼"
Length (H)	24½"	29½"	29½"	35½"
Total mass	29.75 lbs.	34.1 lbs.	35.25 lbs.	40.75 lbs.

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1257.3 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Rabbit Up

	Size 1	Size 2	Size 3	Size 4
Overall measurements				
Width, inside (G1)	12¾"	14¾"	16¾"	16¾"
Overall height (I) *	26 ¼" – 39 ¼"	33 ¾" – 46 ¾"	39" – 51 ½"	47 ¼" – 60 ¼"
Height from foot plate (I1) *	21 ¾" – 29 ¼"	29 ½" – 42 ½"	32 ½" – 46 ¾"	37 ¾" – 56"
Center bar angle (a1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Foot plate angle (a2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	6°	6°	4°
Center bar length	19¾"	27½"	31½"	39¼"

User Measurement Guide for Stander

Recommended user height	30" - 39"	39" - 54"	43" - 58"	48" - 71"
Max armpit height	29"	42"	46"	56"
User weight max	88.1 lbs.	110.2 lbs.	132.2 lbs.	154.3 lbs.
Recommended user age	1 - 3 yr./old	4 - 10 yr./old	5 - 12 yr./old	7 - 18 yr./old

* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

Rabbit Stander Options

Make a selection of either the Prone Stander or the Mobile Prone Stander

Rabbit-Up Prone Stander



Rabbit-Up, Stander, Yellow, Size 1 (E0638) (8621220-00)	\$1648
Rabbit-Up, Stander, Blue, Size 2 (E0638) (8622220-00)	\$1690
Rabbit-Up, Stander, Grey, Size 3 (E0638) (8623220-00)	\$1731
Rabbit-Up, Stander, Silver, Size 4 (E0638) (8624220-00)	\$1858

Rabbit Mobile Prone Stander



22" Wheels, Yellow, Size 1 (E0642) (8621222-21)	\$3159
24" Wheels, Yellow, Size 1 (E0642) (8621223-21)	\$3159
26" Wheels, Yellow, Size 1 (E0642) (8621224-21)	\$3159
29" Wheels, Blue, Size 2 (E0642) (8622225-21)	\$3957
32" Wheels, Blue, Size 2 (E0642) (8622226-21)	\$3957
29" Wheels, Grey, Size 3 (E0642) (8623225-21)	\$4000
32" Wheels, Grey, Size 3 (E0642) (8623226-21)	\$4000
32" Wheels, Silver, Size 4 (E0642) (8624226-21)	\$4126
36" Wheels, Silver, Size 4 (E0642) (8624227-21)	\$4253

Chest Supports

Supports and Pad



Chest Support and Pad

Chest Support and Pad, W: 7.75-10.75", Size 1 (861704)	\$246
Chest Support and Pad, W: 7.75-10.75", Size 2 and 3 (861702)	\$246
Chest Support and Pad, W: 9.75-12.5", Size 4 (861703)	\$271



Filler Pad

Not required for Surrounding Chest Support

Filler Pad, Small (861421)	\$98
Filler Pad, Large (861422)	\$98



Extension for Center Bar

Available with Chest and Surrounding Chest Supports

Extension for Center Bar (861480)	\$98
-----------------------------------	------



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Rabbit Up





Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Rabbit Up

Surrounding Chest Support and Options



Surrounding Chest Support

Surrounding Chest Support, C:18.5-25.25", Small (861810)	\$589
Surrounding Chest Support, C: 26.5-30.75", Medium (861820)	\$613
Surrounding Chest Support, C: 34.25-44.75", Large (861830)	\$649



Cover for Surrounding Chest Support

Cover for Surrounding Chest Support, C:18.5-25.25", Small (869810-160)	\$139
Cover for Surrounding Chest Support, C: 26.5-30.75", Medium (869820-260)	\$171
Cover for Surrounding Chest Support, C: 34.25-44.75", Large (869830-360)	\$182



Extra Pad for Surrounding Support

Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225)	\$75
---	------

Trays



Trays

Tray, Plastic, Grey (862515)	\$463
Tray With Storage Box (862501)	\$475
Plastic Tray with Bowl and Cover (862525)	\$672
<i>Mounting Hardware required for Tray with Bowl</i>	

Tray Accessories

Height Adjustment of Tray (862510)	\$271
<i>Required with Tray with Storage Box</i>	

Mounting Hardware

Required Selection for Tray with Bowl

Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit)	\$320
--	-------

Abdominal Support



Abdominal Supports

Abdominal Support (Each) (861411)	\$193
-----------------------------------	-------

Head Supports

Head Support, Head Support Cover, Connecting System for Head Support, Headrest Mounting Hardware, and Hip Support w/ Pommel are required for a complete Head Support

Head Support



Anatomic Head Support, Small (5099001)	\$285
Anatomic Head Support, Medium (5099002)	\$285
Anatomic Head Support, Large (5099003)	\$285

Connecting System for Head Support




Connecting System for Head Support, Size 1 and 2 (861424-3)	\$310
Connecting System for Head Support, Size 3 and 4 (861424-4)	\$310

Head Support Cover


Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium (9945-83)	N/C
Anatomic Head Support Cover, Large (9967-83)	N/C

Headrest Mounting Hardware



Headrest Mounting Hardware, Short (98619-25)	\$250
Headrest Mounting Hardware, Long (98619-26)	\$256

Pelvic and Hip Control Options



Hip Support w/ Pommel


Front mount to center bar. Required with Head Support.

Hip Support w/ Pommel, Size 1 (861501)	\$607
Hip Support w/ Pommel, Size 2 and 3 (861502)	\$607
Hip Support w/ Pommel, Size 4 (861503)	\$672



Hip Support, With Rotation

Hip Support, With Rotation, Size 2 and 3 (861612)	\$426
Hip Support, With Rotation, Size 4 (861613)	\$463



Hip Support, Front Mount


Hip Support, Front Mount, Size 1 (861601)	\$283
Hip Support, Front Mount, Size 2 and 3 (861602)	\$319
Hip Support, Front Mount, Size 4 (861603)	\$352



Hip Support, Surrounding

Front mount to center bar

Hip Support, Surrounding, Size 1 (861610)	\$667
Hip Support, Surrounding, Size 2 and 3 (861620)	\$667
Hip Support, Surrounding, Size 4 (861630)	\$679



Pommel With Back Support

Pommel With Back Support (861401)	\$524
-----------------------------------	-------

Knee Supports



Knee Support

Knee Supports, Size 1 (Pair) (861301)	\$336
Knee Supports, Size 2 and 3 (Pair) (861302)	\$426
Knee Supports, Size 4 (Pair) (861303)	\$426




Hinged Knee Supports

Hinged Knee Supports, Size 1 (Pair) (861311)	\$524
Hinged Knee Supports, Size 2 and 3 (Pair) (861312)	\$524
Hinged Knee Supports, Size 4 (Pair) (861313)	\$524

Footplate and Sandals Options

Flip-up Footplate required for all Sandals, except the Individual Footrest w/ Hinged Sandals

Footplates and Options




Flip-Up Footplate

Flip-Up Footplate, Size 1 and 2 (861211)	\$265
Flip-Up Footplate, Size 2, 3 and 4 (861212)	\$265



Formed Footplates

Formed Footplates, Size 0 (Pair) (86149)	\$198
Formed Footplates, Size 1 and 2 (Pair) (86147)	\$198
Formed Footplates, Size 3 (Pair) (86148)	\$198



Footplate Straps

Foot Straps, Plastic Lock (Pair) (85096) (K0040)	\$42
Footplate Straps, Velcro (Pair) (85102)	\$59



Individual Footrest w/ Hinged Sandals

Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201)	\$475
Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202)	\$475



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Rabbit Up

Sandals and Options



Sandals

Sandals for Rabbit Up/ Toucan, Small (869490-RT)	\$233
Sandals for Rabbit Up/ Toucan, Medium (869491-RT)	\$244
Sandals for Rabbit Up/ Toucan, Large (869492-RT)	\$251



Sandals

Hinged Sandals w/ Straps, Size 1 (Pair) (861221)	\$162
Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222)	\$162
Hinged Sandals w/ Straps, Size 4 (Pair) (861223)	\$283



Heel Stops

Heel Stops, Small (Pair) (861251)	\$61
Heel Stops, Large (Pair) (861252)	\$71

Conversion Kits

Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit. Conversion kit is used to convert Prone Stander into Mobile Prone Stander.

Driving Wheel Conversion Kits



Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209)	\$694
Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210)	\$694
Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211)	\$694
Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212)	\$1071
Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213)	\$1071
Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214)	\$1071



Central Brake System

Central Brake System, Size 1 (134215)	\$349
Central Brake System, Size 2 and 3 (134216)	\$349
Central Brake System, Size 4 (134218)	\$349



Camber Kit Rabbit Up

Camber Kit Rabbit Up 0° (134200)	\$63
Camber Kit Rabbit Up +3/-3° (134201)	\$63
Camber Kit Rabbit Up +6/-6° (134202)	\$63