

R82 Stingray - E1233

USD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

US Configuration


Effective: 01.01.23 - Rev 10

- Frame (Stingray, Trekker, High-Low:v)
- Stingray Seat
- Cushions
- Head Support
- Side Support
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Transport Fittings



- A** (Back Height) _____
- B** (Arm Height) _____
- C** (Lower Leg Length) _____
- D** (Seat Depth) _____
- E** (Chest Width) _____
- F** (Seat Width) _____
- G** (Height: Seat to Top of Head) _____
- H** (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



	Size 1 Inches	Size 2 Inches
Effective Seat Width (A)	7-11 ³ / ₄ "	9-13 ³ / ₄ "
Seat Width Between Knee Supports	11"	12 ³ / ₄ "
Seat Width Between Hip Supports (Accessories)	7", 8 ¹ / ₂ ", 10 ¹ / ₄ "	9", 10 ¹ / ₂ ", 12"
Effective Seat Depth (B)	7 - 11 ³ / ₄ "	9 ³ / ₄ - 15 ¹ / ₂ "
Back Height at Shoulders (C)	10 ¹ / ₂ - 16"	14 - 19 ³ / ₄ "
Backrest Height (C1)	16 - 21 ¹ / ₂ "	19 ¹ / ₂ - 25 ¹ / ₄ "
Footrest to Seat Distance (D)	5 ¹ / ₂ - 17 ¹ / ₂ "	5 ¹ / ₂ - 17 ¹ / ₂ "
Seat Surface Height at Front Edge (E)	20 ¹ / ₄ "	20 ¹ / ₄ "
Back Width, Inside (F)	11 ³ / ₄ "	14"
Back Width, Outside (F1)	17 ¹ / ₂ "	20 ¹ / ₄ "
Overall Width (G)	23"	23"
Length with Turnable Front Wheels (H1)	28 ³ / ₄ "	28 ³ / ₄ "
Length with Large Front Wheels	32 ³ / ₄ "	32 ³ / ₄ "
Overall Height, Push Brace 15° (I)	39", 40", 41 ¹ / ₄ "	39", 40", 41 ¹ / ₄ "
Foot Plate, Width x Depth	12" x 7 ³ / ₄ "	12" x 7 ³ / ₄ "
Folded Length	35"	35"
Folded Height	19"	19"
Total Mass, Frame	29 ¹ / ₂ lbs.	29 ¹ / ₂ lbs.
Total Mass, Seat	12 lbs.	14 ¹ / ₂ lbs.
Max Load	88 lbs.	88 lbs.
Back Recline / Seat Plane Angle	Back: 170° / Seat: 45°	Back: 170° / Seat: 45°
Approx. Minimum User Height	29 ¹ / ₂ "	35 ¹ / ₂ "
Approx. Maximum User Height	36"	47 ¹ / ₂ "

NOTES:

Stingray Standard Features

- Tilt-in-Space - 45°
- Back Recline - 170°
- Turnable Seat - 180°
- One-Size Frame
- Accepts Two Seat Sizes for Maximum Growth Potential
- Adjustable Seat Depth
- Adjustable Back Height
- Suspension



Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Solid Tires, Individual Wheel Locks	8724212	\$2,550	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Hand Brakes ●●	8724211	\$2,942	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks	8724222	\$2,655	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks & Hand Brakes	8724221	\$3,045	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Quick Release Rear Axels	8722212	\$2,810	E1233

●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item

Trekker Frame Standard Features

- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- Single Foot Activated Wheel Lock Bar
- 7½" x 2" Front and 11½" x 2½" Rear Solid Tires
- Quick Release Rear Wheels
- Made in the USA
- Transit Option - *Transit Ready Frame Only*



Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Trekker Frame (TR)	C904661	\$2,560	E1234
<input type="checkbox"/> Trekker Frame - Transit Ready (TR)	C905362	\$2,774	E1234
<input type="checkbox"/> Stingray Seat Module - Interface Adaptor Kit for Trekker Base ♦	C905363	N/C	

♦Required with Trekker Frame

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Transit Option (TR)	C904576	\$214	
<input type="checkbox"/> Using other manufacturer's Head Support	C995681	N/C	

♦Requires Head Support & Pelvic Positioning Belt

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Height Adjustable Push Handle w/Outside Mount (TR)	C904688	N/C	


^Trekker Frame comes standard with push buttons mounted inside Height Adjustable Push Handle


Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> 12½ x 2" Rear Solid Knobby Tire (TR) ●●	C904117	\$213	
<input type="checkbox"/> 7½" x 2" Front & 11½" x 2½" Rear Pneumatic Knobby Tires (TR) ●●	C904118	\$319	

●●Custom order; Approx. lead time: 6-8 weeks

Trekker Frame Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/>	IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required		One Size	C904566	\$209	K0105
<input type="checkbox"/>	O2 Tank Holder (TR) *°		One Size	C904565	\$182	E2208
<input type="checkbox"/>	LTV Ventilator Holder (TR) °		One Size	C904671	\$469	E1029
<input type="checkbox"/>	Ventilator/Suction/Hard Tray (TR)		One Size	C904595	\$469	E1029
<input type="checkbox"/>	Curb Tipper (TR)		One Size	C905310	\$102	
<input type="checkbox"/>	Rear Anti-Tipper (TR)		One Size	C904649	\$120	E0971
<input type="checkbox"/>	Under Seat Storage Basket (TR)		One Size	C904605	\$125	
<input type="checkbox"/>	Medical Necessity Storage Basket (TR)		One Size	C904606	\$307	
	* Must Be Removed During Transit					
	° Manufacturer Installation Required					

High-Low Frame (Trekker)	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:xo* 5" castor, height adj. & tilt on push handle		1,2	89010522	\$1,899	
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	890009-04	\$190	
	<input type="checkbox"/> Push Brace - Angle Adjustable		1,2	8910555	\$402	
	*Requires Push Brace (890009-04 or 8910555)					

High-Low Frame (Stingray)	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:xo* 5" castor, height adj. & tilt on push handle		1	89040522	\$2,034	
	<input type="checkbox"/> High-Low:xo* 5" castor, height adj. & tilt on push handle		1, 2	89050522	\$2,287	
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	890009-04	\$190	
	<input type="checkbox"/> Push Brace, Angle Adjustable		One Size	8910555	\$402	
	*Requires Push Brace (890009-04 or 8910555)					



Stingray Seat	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Seat and Back size 1		1	8711120		
	<input type="checkbox"/> Seat size 1		1	x8711000	\$427	E2292
	<input type="checkbox"/> Back size 1		1	x8712000	\$527	E2291
	<input type="checkbox"/> Seat and Back size 2		2	8712220		
	<input type="checkbox"/> Seat size 2		2	x8721000	\$427	E2292
	<input type="checkbox"/> Back size 2		2	x8722000	\$527	E2291
	<input checked="" type="checkbox"/> Fully Reclining Back (Required)		One Size	x8702000	\$422	E1014

Adjustable seat depth and height adjustable back

Description	Measurements	Size	Part No.	USD Price	HCPCS
Required, must select a color					
X = Color					
<input type="checkbox"/> Blue = 1			<input type="checkbox"/> Red = 2		
<input type="checkbox"/> Pink = 5			<input type="checkbox"/> Grey = 4		
<input type="checkbox"/> Seat and Back Cushion (Pr.)		1	875051x		N/C
<input type="checkbox"/> Seat and Back Cushion (Pr.)		2	875052x		N/C
*Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item					

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Footrest Rod, Elevating ♦	8¼"	Short	x8711100	\$115	K0108
<input type="checkbox"/> Footrest Rod, Elevating ♦	11"	Long	x8712200	\$115	K0108
♦Required with Angle Adjustable Footplate (x8710010)					
<input type="checkbox"/> Angle Adjustable Footplate ♦		One Size	x8710010	\$99	K0040
♦Requires Footrest Rod Short or Long (x8711100 or x8712200)					

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Transit Option*		One Size	8720902	\$209	
	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Hip Supports, Black (Pr.)	6½" x 7"	1	8740230	\$235	E0956
	<input type="checkbox"/> Hip Supports, Black (Pr.)	7" x 8½"	2	8740240	\$235	E0956
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		1	8740250	\$235	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		1	x8740250	\$235	E1028
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		2	8740260	\$235	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		2	x8740260	\$235	E1028
	♦ Hip Support and Multi-Adjustable Hardware must be ordered together					
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	4" x 4"	1	8720301	\$0	E0956
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	5" x 5¾"	2	8720302	\$0	E0956
	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Covers for Knobs on Side Supports (Pr.) ●●		One Size	8720360	\$32	
	●● Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Medial Thigh Support (Ea.) ♦	3½" x 1½" x 2¾"	1	8720501	\$140	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		1	x8720501	\$172	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 3½"	2	8720502	\$140	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		2	x8720502	\$172	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 4¾"	3	8720503	\$140	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		3	x8720503	\$172	E1028
	♦ Medial Thigh Support and Adjustable, Removable Hardware must be ordered together					
	<input type="checkbox"/> Anti-Thrust Wedge, Black (Pr.)	4" x 6"	1, 2, 3	8910400-3	\$48	K0108
	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Contoured Headrest, Black		1	8740810	\$209	E0955
	<input type="checkbox"/> Contoured Headrest, Black		2	8740820	\$209	E0955
	<input type="checkbox"/> Sculpted Head Support ①	9½" x 4½"	Small	8720831	\$214	E0955
	<input type="checkbox"/> Sculpted Head Support ①	11" x 5¾"	Large	8720832	\$214	E0955
	① Preferred option for Stingray seating					
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$266	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C	
	♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26					
	Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable					
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$240	E1028
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$240	E1028
	♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003					

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> 2-Part Headrest		One Size	8720870	\$282	E0955
	<input type="checkbox"/> Headrest Fittings ♦ (Mounts directly onto backrest)		One Size	8720850	\$53	K0108
	<input type="checkbox"/> Universal Headrest Mount ♦Required when using other manufactures Head Support Allow most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation		One Size	87590	\$64	E1028

Other types of headrest available, please see Head Supports Order Form


Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------



■ XX = 60 (Black)



<input type="checkbox"/> Cross Vest (Ea.)	S-P: 7¾–12½"	XS	85720-XX	\$198	E0960
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 10¼–13½"	S	85720-S	\$198	E0960
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 13½ - 17¾"	M	85720-M	\$198	E0960
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 17¾–20¾"	L	85720-L	\$198	E0960
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 20¾ - 24 ½"	XL	85720-XL	\$198	E0960



<input type="checkbox"/> Vest (Ea.)	1¼" x 2½" x 5½" x 3½" x 11"	0	89111-0XX	\$209	E0960
<input type="checkbox"/> Vest (Ea.)	1¼" x 2½" x 5¾" x 4¼" x 12"	¼	89111-075XX	\$209	E0960
<input type="checkbox"/> Vest (Ea.)	2¼" x 2½" x 7" x 5" x 13¼"	1	89111-1XX	\$209	E0960
<input type="checkbox"/> Vest (Ea.)	2½" x 3½" x 9¼" x 5½" x 15¼"	2	89111-2XX	\$209	E0960
<input type="checkbox"/> Vest (Ea.)	2¾" x 4¾" x 11 x 5¾" x 18¼"	3	89111-3XX	\$209	E0960
<input type="checkbox"/> Vest (Ea.)	2¾" x 6¾" x 13¾" x 7½" x 19"	4	89111-4XX	\$209	E0960



<input type="checkbox"/> Torso vest (Ea.)	S-P: 7¾–10¾"	XS	85705-XS	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 10¾–12½"	S	85705-S	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 12½–16"	M	85705-M	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 16–18½"	L	85705-L	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 18½–20½"	XL	85705-XL	\$209	E0960
<input type="checkbox"/> Torso vest (Ea.)	S-P: 20½–23½"	XXL	85705-XXL	\$209	E0960




<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 6¾–12½"	XS	85725-XS	\$140	E0978
<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 8¼–13¾"	S	85725-S	\$140	E0978
<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 9¾–16¼"	M	85725-M	\$140	E0978
<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 12½–17¼"	L	85725-L	\$140	E0978
<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 16¼–20½"	XL	85725-XL	\$140	E0978



<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$101	E0978
---	------------	---	---------	-------	-------



<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦ ♦Required for use of Vests		One Size	862109-35	\$44	
---	--	----------	-----------	------	--






<input type="checkbox"/> Foot Straps, Plastic Lock (Pr.)	Max: 22", A: 13½-20½"	One Size	85096	\$38	K0040
<input type="checkbox"/> Foot Straps, Velcro (Pr.)		One Size	85102	\$53	K0040

Other types of vests and belts available, please see Vest and Belts Order Form

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Folding Top (Canopy), Black		1	8740110	\$344	
	<input type="checkbox"/> Folding Top (Canopy), Black		2	8740120	\$344	
	<input type="checkbox"/> Foam Back Plate ●●	5" x 12"	1	8776325	\$24	
	<input type="checkbox"/> Foam Back Plate ●●	5½" x 13¾"	2	8776326	\$24	
	Future Growth - Fills out the gap between seat and the back					
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Tray with Edge, Grey ◆	17.3" x 13.3"	One Size	8776350	\$105	E0950
	<input type="checkbox"/> Adjustable Removable Tray Hardware ◆		One Size	x8776350	\$172	E1028
	◆ Tray and Adjustable Removable Hardware must be ordered together					
	<input type="checkbox"/> Rain Cover ◆		1	8720141	\$58	
	<input type="checkbox"/> Rain Cover ◆		2	8720142	\$58	
	◆ Requires Folding Top (Canopy)					
	<input type="checkbox"/> Handle Grip	14"	1	8720401	\$209	
	<input type="checkbox"/> Handle Grip	16"	2	8720402	\$209	
	<input type="checkbox"/> Backpack, Black ●●		Small	87043183-1	\$58	
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Medical Essential Carrier		One Size	8720600	\$101	
	<input type="checkbox"/> Vent Tray ◆*	13" x 17"	One Size	8776290	\$480	E1029
	◆ Requires Central Wheel Locks					
	*44 lbs. weight capacity					
	<input type="checkbox"/> O2 Tank Holder, Horizontal Mounted *			8776295	\$365	E2208
	*Cannot be used during transit and manufacturer installation is required					
	<input type="checkbox"/> O2 Tank Holder, Vertical Mounted *			8776296	\$107	E2208
	*Cannot be used during transit and manufacturer installation is required					
	<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required		One Size	C904566	\$209	K0105
	<input type="checkbox"/> Mounting bracket for Stingray backrest*		One Size	8776291	\$147	
	*Required with IV Pole #C904566					

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	Required, must select a color					
	X = Color					
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2				
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4				
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		1	875051x	\$452	
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		2	875052x	\$452	
	●●Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		1	8750410	\$58	
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		2	8750420	\$58	
	<input type="checkbox"/> Hardware Kit for Hip Supports, Fixed		One Size	872023-KIT	\$37	
	<input type="checkbox"/> Directional Arrow, Replacement		One Size	149889	\$6	

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Release Strap for Back Angling			1	87102250	\$48
	<input type="checkbox"/> Release Strap for Back Angling			2	87102255	\$48
	<input type="checkbox"/> Footrest Rod, Elevating	8¼"	Short	8775317	\$115	
	<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	8775318	\$115	
	<input type="checkbox"/> Angle Adjustable Footplate		One Size	8775010	\$272	

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

The R82 Stingray mobility base was tested with a footrest and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

Prices in USD & is subject to change without notice

For Quotes and Orders please send completed Order Form to: sales.us@etac.com or Fax: 1-800-760-0537

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF6300-ENG-REV.09


Creating Possibilities