



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Flamingo Curo

Company:	Ship To:		
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	State:
ATP:	Zip:	Phone:	
Marked For:	Email/Fax:		
Associated Quote/Order#:		Effective: 4/13/2026	

## Base Frame and Mobility Base Options

### Flamingo Curo Seat and Back



Seat and Back, Size 1 (882001)	\$675
Seat and Back, Size 2 (882002)	\$725
Seat and Back, Size 3 (882003)	\$795

### Seat and Back Upholstery

*Upholstery must match seat size selection*



#### Seat Base Upholstery Inlay

Seat Base Upholstery Inlay, Blue, Size 1 (882071)	\$149
Seat Base Upholstery Inlay, Blue, Size 2 (882072)	\$162
Seat Base Upholstery Inlay, Blue, Size 3 (882073)	\$175



#### Back Upholstery Inlay

Back Upholstery Inlay, Blue, Size 1 (882075)	\$153
Back Upholstery Inlay, Blue, Size 2 (882076)	\$160
Back Upholstery Inlay, Blue, Size 3 (882077)	\$175

### Flamingo Curo Base Frame

*Frames must match seat size selection*

#### Flamingo Curo Compact Frame



Flamingo Curo Compact Frame with 3" Casters, Size 1 (882011)	\$995
Flamingo Curo Compact Frame with 3" Casters, Size 2 (882012)	\$1030
Flamingo Curo Compact Frame with 3" Casters, Size 3 (882013)	\$1075



#### Potty Frame

Potty Frame, Size 1 (882050)	\$325
------------------------------	-------

*Available only with Size 1 Seat*

#### Flamingo Curo High-Low Frame



Flamingo Curo High-Low Frame with 3" Casters, Size 1 (882031)	\$3150
Flamingo Curo High-Low Frame with 3" Casters, Size 2 (882032)	\$3250
Flamingo Curo High-Low Frame with 3" Casters, Size 3 (882033)	\$3350



#### Hardware for Toilet Mounting

Direct Toilet Mounting, Size 1 (882061)	\$270
Direct Toilet Mounting, Size 2 (882062)	\$285
Direct Toilet Mounting, Size 3 (882063)	\$295



#### Suction Cups

Suction Cups (880150)	\$175
-----------------------	-------

*Requires Commode Pan / Potty Ring Holders*



#### Casters

4" Casters (880100-S)	\$356
-----------------------	-------

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1425.0 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Flamingo Curo

## Footrest

*Footrest options must match seat size selection*



### Foot Support with Quick Release

Removeable Foot Support, Size 1 (882151)	\$362
Removeable Foot Support, Size 2 (882152)	\$384
Removeable Foot Support, Size 3 (882153)	\$405



### Calf Support

Calf Support, Size 1 (882135)	\$152
Calf Support, Size 2 (882136)	\$155
Calf Support, Size 3 (882137)	\$159

## Pelvic and Trunk Positioning

### Pelvic Positioning



### Seat Infill

*Must match seat size selection*

Seat Infill for Seat, Blue PU, Size 1 (882121)	\$150
Seat Infill for Seat, Blue PU, Size 2 (882122)	\$150
Seat Infill for Seat, Blue PU, Size 3 (882123)	\$150



### Splashguard

Splashguard, Small (882145)	\$130
Splashguard, Large (882147)	\$130



### Hip Guide

*Requires Armrests or Hip Guide Attachment Rod*

Hip Guide, Small (882111)	\$175
Hip Guide, Large (882112)	\$180



### Hip Belt

Neoprene Hip Belt, Size 1 (85423-1)	\$107
Neoprene Hip Belt, Size 2 (85423-2)	\$107
Neoprene Hip Belt, Size 3 (85423-3)	\$114



Hip Guide Attachment Rod (882116)	\$30
-----------------------------------	------

*Not compatible with Armrests*

### Trunk Positioning



### Shoulder Support with PU

*Not compatible with Side Supports. Must match seat size selection.*

Shoulder Support Inlay, Blue PU, Size 1 (882081)	\$225
Shoulder Support Inlay, Blue PU, Size 2 (882082)	\$243
Shoulder Support Inlay, Blue PU, Size 3 (882083)	\$265



### Lower Back Support

*Must match seat size selection*

Lower Back Support, Size 1 (882131)	\$165
Lower Back Support, Size 2 (882132)	\$165
Lower Back Support, Size 3 (882133)	\$165



### Fixed Side Supports

Fixed Side Supports, Small (882090)	\$180
Fixed Side Supports, Medium (882091)	\$195
Fixed Side Supports, Large (882092)	\$215



### Swing-Away Side Supports

Swing-Away Side Supports, Small (882095)	\$325
Swing-Away Side Supports, Medium (882096)	\$345
Swing-Away Side Supports, Large (882097)	\$365



### Cross Vest

Cross Vest, Net, Size 0 (89625-040)	\$214
Cross Vest, Net, Size 1 (89625-140)	\$214
Cross Vest, Net, Size 2 (89625-240)	\$214
Cross Vest, Net, Size 3 (89625-340)	\$214
Cross Vest, Net, Size 4 (89625-440)	\$214



### Vest

Vest, Net, Size 0 (89111-040)	\$219
Vest, Net, Size 3/4 (89111-07540)	\$219
Vest, Net, Size 1 (89111-140)	\$219
Vest, Net, Size 2 (89111-240)	\$219
Vest, Net, Size 3 (89111-340)	\$219
Vest, Net, Size 4 (89111-440)	\$219



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Flamingo Curo

## Headrests and Hardware



Head Support, Small, Curved, Standard Hardware (9931056-61)	\$465
Head Support, Medium, Curved, Standard Hardware (9932056-61)	\$485
Head Support, Small, Curved, Extended Hardware (9931066-61)	\$475
Head Support, Medium, Curved, Extended Hardware (9932066-61)	\$495

## Armrest Options



### Armrests

Armrests, Small (882113)	\$235
Armrests, Large (882114)	\$255



### Tray

*Requires Armrests to mount. Must match seat size selection.*

Tray, Blue, Size 1 (882141)	\$250
Tray, Blue, Size 2 (882142)	\$275
Tray, Blue, Size 3 (882143)	\$295

## Accessories

### Handles and Push Braces



### Push Brace

Push Brace (880110)	\$198
---------------------	-------



### Push Handle

*Requires High-Low Frame*

Helper Handle, Size 1 (880110-HL)	\$78
Helper Handle, Sizes 2 and 3 (880110-HLL)	\$78

### Potty Options



### Commode Pan

*Requires Commode Pan / Potty Ring Holders*

Commode Pan, Size 1 (870800)	\$20
Commode Pan, Sizes 2 and 3 (870801)	\$45



### Potty Ring

*Requires Commode Pan / Potty Ring Holders*

Potty Ring, Size 1 (870805)	\$45
Potty Ring, Sizes 2 and 3 (870807)	\$55



### Commode Pan Lid

Commode Pan Lid, Size 1 (870800-L)	\$84
Commode Pan Lid, Sizes 2 and 3 (870802-L)	\$84



### Extended Potty Ring

*Requires Commode Pan / Potty Ring Holders*














Potty Ring, Extended, Sizes 2 and 3 (870807-E)	\$55
--	------



### Commode Pan / Potty Ring Holders















Commode Pan / Potty Ring Holders, Size 1 (880125-1)	\$30
Commode Pan / Potty Ring Holders, Size 2 (880125-2)	\$35
Commode Pan / Potty Ring Holders, Size 3 (880125-3)	\$40

## Measurement Guide









Seat size 1, 2 & 3			
mm (inch) kg (lb)	 1	 2	 3
 Backrest recline angle	(-20°)-(30°)	(-20°)-(30°)	(-20°)-(30°)
 Overall backrest height**	425-495 mm (16¾-12½")	470-550 mm (18½-21¾")	515-615 mm (20¼-24¼")
 Effective backrest height**	375/410/445 mm (14¾/16¼/17½")	405/445/485 mm (16/17½/19")	450/500/550 mm (17¾/19¾/21¾")
 Effective back width (with side supports)	140-240 mm (5½-9½")	170-270 mm (6¾-10¾")	220-320 mm (8¾-12½")
 Effective back width (without side supports)	300 mm (11¾")	335 mm (13¼")	390 mm (15¼")
 Effective seat depth**	240-300 mm (9½-11¾")	285-345 mm (11¼-13½")	335-425 mm (13¼-16¾")
 Overall length	360-375 mm (14¼-14¾")	395-430 mm (15½-17")	435-495 mm (17¼-19½")
 Effective seat width (without upper leg supports)	265 mm (10½")	300 mm (11¾")	355 mm (14")
 Effective seat width (with upper leg supports)	185/225/265 mm (7¼/8¾/10½")	220/260/300 mm (8¾/10¼/11¾")	275/315/355 mm (10¾/12½/14")
 Effective foot support to seat distance	220-285 mm (8¾-11¼")	260-345 mm (10¼-13½")	300-400 mm (11¾-15¾")
 Mass device*	5.2 kg (11.5 lb)	6.1 kg (13.4 lb)	7.5 kg (16.5 lb)
 Maximum user mass	<35 kg (<77 lb)	<55 kg (<121 lb)	<70 kg (<154 lb)
 Max rated load	42.2 kg (93 lb)	63.1 kg (139.1 lb)	79.5 kg (175.3 lb)

\*seat only (without accessories) \*\*including PUR

# Flamingo Curo

Seat & Compact frame				Seat & High-low frame			
mm (inch) kg (lb)	 1	 2	 3	 1	 2	 3	
 Overall height*	905-1000 mm (35¾-39¾")	950-1055 mm (37½-41½")	995-1120 mm (39¾-44")	895-1075 mm (35¾-42¾")	930-1205 mm (36½-47½")	975-1270 mm (38½-50")	
 Overall height, to pushbrace*	850-1020 mm (33½-40¾")	890-1070 mm (35-42¾")	930-1130 mm (36½-44½")	900-1180 mm (35½-46½")	945-1235 mm (37¼-48½")	990-1350 mm (39-53¼")	
 Overall length with leg rest*	935 mm (36¾")	955 mm (37½")	1110 mm (43¾")	930 mm (36½")	1065 mm (42")	1200 mm (47¼")	
 Overall width*	470 mm (18½")	510 mm (20")	560 mm (22")	485 mm (19")	580 mm (22¾")	580 mm (22¾")	
 Mass device*	12.4 kg (27.3 lb)	13.4 kg (29.5 lb)	14.9 kg (32.8 lb)	14.8 kg (32.6 lb)	17 kg (37.5 lb)	18.1 kg (39.9 lb)	
 Overall minimum turning diameter*	850 mm (33½")	875 mm (34½")	905 mm (35¾")	785 mm (31")	955 mm (37½")	955 mm (37½")	
 Seat surface height at front edge***	480-580 mm (19-22¾")	480-580 mm (19-22¾")	480-580 mm (19-22¾")	530-740 mm (20¾-29¾")	530-740 mm (20¾-29¾")	530-740 mm (20¾-29¾")	
 Seat plane angle (tilt)**	(-10°)-(30°)	(-10°)-(30°)	(-10°)-(30°)	(-5°)-(20°)	(-5°)-(20°)	(-5°)-(20°)	

Seat & Potty frame		mm (inch) kg (lb)					
	 Overall height	 Overall length	 Overall length with leg rest***	 Overall width	 Mass device	 Effective foot support to seat distance***	 Seat surface height at front edge***
 1	630-750 mm (24¾-29½")	555 mm (21¾")	870 mm (34¾")	430 mm (17")	7.1 kg (16 lb)	220-285 mm (8¾-11¼")	290-330 mm (11½-13")

\*seat and High-low frame with 75 mm wheels    \*\*(+/-5°)    \*\*\*including PUR