



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 3/12/2025

x:panda Shape - Advanced

	Infant	1	2
 Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
 Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
 Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
 Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
 Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
 Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
 Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Shape Base Options

Please see measurement guide for growth information. A selection must be made in each section below.

Seat Base

Includes base and plate

Seat Base Size

Size 1 Advanced Seat (E2292)	\$2851
Size 2 Advanced Seat (E2292)	\$2851

Recline Options

Selection required

Reclining Dynamic Back (K0108)	N/C
Reclining Adjustable Dynamic Back (K0108)	N/C



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Back and Cushions

Upper Back Size

"Select Upper Back Base then select Upper and Lower Back Cushion and Covers. Includes quick release detachment (#9050374)"



Size 1 Base

Size 1 Base Upper - H: 7" (905073) (E2291) _____ N/C
Total back height range 12.25" - 16.25"
Size 1 Base Upper - H: 8.75" (905074) (E2291) _____ N/C
Total back height range 13.75" - 18.50"

Size 2 Base

Size 2 Base Upper - H: 7" (905073) (E2291) _____ N/C
Total back height range 13.5" - 17.25"
Size 2 Base Upper - H: 8.75" (905074) (E2291) _____ N/C
Total back height range 15" - 19.75"
Size 2 Base Upper - H: 10.75" (905075) (E2291) _____ N/C
Total back height range 17" - 21.75"

Upper Back Cushion and Cover

Upper back cushion height must match Upper Back Base height

Upper Back Cushion

Upper Back Cushion - W: 7" x H: 7" (902103) _____ N/C
Upper Back Cushion - W: 8.75" x H: 7" (902104) _____ N/C
Upper Back Cushion - W: 10.25" x H: 7" (902105) _____ N/C
Upper Back Cushion - W: 8.75" x H: 8.75" (902106) _____ N/C
Upper Back Cushion - W: 10.25" x H: 8.75" (902107) _____ N/C
Upper Back Cushion - W: 11.75" x H: 8.75" (902108) _____ N/C
Upper Back Cushion - W: 13.5" x H: 8.75" (902109) _____ N/C
Upper Back Cushion - W: 10.25" x H: 10.75" (902110) _____ N/C
Upper Back Cushion - W: 11.75" x H: 10.75" (902111) _____ N/C
Upper Back Cushion - W: 13.5" x H: 10.75" (902112) _____ N/C
Upper Back Cushion - W: 15" x H: 10.75" (902113) _____ N/C

Upper Back Cover

Selection required

Red (-50) _____ N/C
Grey (-60) _____ N/C
Black Wipeable (-70) _____ N/C

Lower Back Cushion and Cover

The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.

Lower Back Cushion

Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120) _____ N/C
Size 1 Lower Back Cushion - W: 10.25" x H: 4.75" (902121) _____ N/C
Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902122) _____ N/C
Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123) _____ N/C
Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124) _____ N/C
Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125) _____ N/C
Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126) _____ N/C

Lower Back Cover

Selection Required

Red (-50) _____ N/C
Grey (-60) _____ N/C
Black Wipeable (-70) _____ N/C

Daily Adjustment Kit for Upper Back

Allows for tool-less adjustment of back height. Ideal for frequent adjustment.



Daily Adjustment Kit - Upper Back (905081) _____

\$26

Spacer Blocks

Spacer Block Holder without Cushion

Cushion width must match spacer block holder width

Spacer Block Holder, Size 1 - Back W: 7", 8.75", 10.25" _____ \$224
Spacer Block Holder, Size 2 - Back W: 8.75", 10.25", 11.75", 13.5" _____ \$224

Spacer Block Cushion Cover

Black (-80) _____ N/C
Black Wipeable (-70) _____ N/C

Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011) _____ N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012) _____ N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013) _____ N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014) _____ N/C
Spacer Block Cushion Size 5 - W: 13.5" x H: 2.25" (903015) _____ N/C





Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Cushions and Covers

Seat Cushion

Seat Cushion - Split

Required when using Quick Release Pommel

Split Cushion, Size 1- W: 8.75" x L: 19.75" (901113)	N/C
Split Cushion, Size 1- W: 10.25" x L: 19.75" (901114)	N/C
Split Cushion, Size 1- W: 11.75" x L: 19.75" (901115)	N/C
Split Cushion, Size 2- W: 12.5" x L: 25.5" (901117)	N/C
Split Cushion, Size 2- W: 14.25" x L: 25.5" (901118)	N/C
Split Cushion, Size 2- W: 15.75" x L: 25.5" (901119)	N/C

Color

Required selection for cover color

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$26

Headrest and Hardware Options

Headrest Style

Includes Cushion

Type of Head Support

Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)	\$275
Curved, Small - Black - W: 7", H: 5" (E0955)	\$275
Curved, Medium - Black Wipeable (E0955)	\$300
Curved, Medium - Black (E0955)	\$300
Anatomic, Small - Black Wipeable (E0955)	\$339
Anatomic, Small - Black (E0955)	\$339
Anatomic, Medium - Black Wipeable (E0955)	\$351
Anatomic, Medium - Black (E0955)	\$351
Anatomic, Large - Black Wipeable (E0955)	\$365
Anatomic, Large - Black (E0955)	\$365
Adjustable with Large Laterals - Black (E0955)	\$406
Adjustable with Large Laterals and Speakers - Black (E0955)	\$539
Adjustable with Four Laterals - Black (E0955)	\$504

Mounting Hardware for Curved and Anatomic Headrest

Required with Curved and Anatomic Head support

Headrest Mounting Hardware

Mounting Brackets, Standard (E1028)	\$306
5.5" of forward/backward adjustment	
Mounting Brackets, Extended (E1028)	\$320
8.75" of forward/backward adjustment	
Mounting Brackets, 3D Adjustment (E1028)	\$332
8.75" of forward/backward adjustment	

Mounting Hardware for Adjustable Headrest Style

Required with Adjustable Head Support

Short Mounting Hardware (E1028)	\$306
Suitable for upper back heights of 7", 8.75", 10.75" and 13.5"	
Long Mounting Hardware (E1028)	\$306
Suitable for upper back heights of 10.75" and 13.5"	

Headrest Options



Daily Adjustment Kits

*Creates toolless adjustability, For angle adjustment of head support.
For use with Anatomic and Curved Headrests only.*

Daily Adjustment Kit Standard - Set of Two	\$39
Daily Adjustment Kit Extended - Set of Three	\$59

Additional Headrest Options

Headband System	\$306
For use with Anatomic Headrest	
Universal Headrest Mount	\$82
Required if using a third party headrest	

Arm Supports

Arm Supports

Arm Supports Size

Cover not included, select separately

Size 2 - W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$416
Size 3 - W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$416

Arm Support Cover

Black	\$96
Black Wipeable	\$96

Daily Adjustment Kit

For toolless height adjustment of arm supports

Height Adjustment Kit of Arm Supports (8910615-11)	\$26
--	------



Upper Arm Support - Fixed

Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.

Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$320
For size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$320
For size 2 Base	

Upper Arm Supports - Adjustable

Not compatible with Swing-Away Lateral Supports.

Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$428
For Size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$428
For size 2 Base	



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Knee and Postural Supports

Abduction and Adduction

Inner Knee Support



For Size 1 and 2 Black - L: 7.5", H: 7.5" (905001-80) (E0957)	\$452
For Size 1 and 2 Black Wipeable - L: 7.5", H: 7.5" (905001-70) (E0957)	\$452
For Size 2 Black - L: 9.5", H: 9.5" (905002-80) (E0957)	\$452
For Size 2 Black Wipeable - L: 9.5", H: 9.5" (905002-70) (E0957)	\$452

Outer Knee Support



Standard Set for Size 1, Black - L: 7.5", H: 7.5" (905008-80) (E0956)	\$300
Standard Set for Size 1, Black Wipeable - L: 7", H: 7.5" (905008-70) (E0956)	\$300
Standard Set for Size 2, Black - L: 9.5", H: 9.5" (905009-80) (E0956)	\$300
Standard Set for Size 2, Black Wipeable - L: 9.5", H: 9.5" (905009-70) (E0956)	\$300
Flip Down Set for Size 2, Black - L: 7.75", H: 7.75" (905010-80) (E0956)	\$357
Flip Down Set for Size 2, Black Wipeable - L: 7.75", H: 7.75" (905010-70) (E0956)	\$357

Quick Release Pommel

For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$447
For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$447
For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905078-80) (E0957)	\$447
For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905078-70) (E0957)	\$447
For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905082-80) (E0957)	\$447
For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905082-70) (E0957)	\$447

Postural Supports

Width Adjustable Pelvis Supports

Mounted from the back



Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" (905041-80) (E0956)	\$300
Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" (905041-70) (E0956)	\$300
Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" (905042-80) (E0956)	\$300
Width Adjustable Pelvis Supports, Large Black Wipeable - L: 5.5", H: 5" (905042-70) (E0956)	\$300

Width Adjustable Hip and Pelvis Supports

Mounted from the sides



Width Adjustable Hip Supports, Small Black - L: 4" H: 4.25" (905044-80) (E0956)	\$428
Width Adjustable Hip Supports, Small Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$428
Width Adjustable Hip Supports, Large Black - L: 5.5" H: 5" (905045-80) (E0956)	\$428
Width Adjustable Hip Supports, Large Black Wipeable - L: 5.5", H: 5" (905045-70) (E0956)	\$428

Adjustable Anti-Thrust Wedges

Short - W: 4.5", L: 2.75-4.25", H: .75" (905052) (K0108)	\$141
<i>For Size 1</i>	
Medium - W: 5.25", L: 4.25-6.75", H: 1.25" (905053) (K0108)	\$141
<i>For Size 1</i>	
Long - W: 5.25", L: 4.75-7.75", H: 1.75" (905054) (K0108)	\$141
<i>For Size 1</i>	
Medium - W: 6", L: 5-8.75", H: 1.25" (905055) (K0108)	\$141
<i>For Size 2</i>	
Long - W: 6", L: 5.50-9.5", H: 2" (905056) (K0108)	\$141
<i>For Size 2</i>	

Trays and Accessories



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950)	\$357
<i>For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.</i>	



Shape 6 Tray with Edge Size 2 (8910725-2) (E0950)	\$357
<i>For Size 2. Double rod, Grey ABS.</i>	



Shape 8 Tray with Edge Size 1 (8976350) (E0950)	\$357
---	-------

Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028)	\$459
Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028)	\$327
Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028)	\$459
Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028)	\$327
Arm Pads for Tray (8910730)	\$386
<i>For shape 5 and 6</i>	
Toy Bar (8910740)	\$306



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Accessory Holders and Brackets

Remote Control Holder for Size 1 and 2 - Arm Support Mounted (905030-1)	\$84
Remote Control Holder for Size 3 - Arm Support Mounted (905030-2)	\$84
Remote Control Holder - Back Mounted (905030)	\$161

Mouting Hardware

Bracket for Communication Aid, Size 1 and 2 - Arm Support Mounted (905033-1)	\$84
Bracket for Communication Aid, Size 3 - Arm Support Mounted (905033-2)	\$84

Canopy

Canopy does not come with the brackets. Both the Canopy and brackets must be quoted when mounting to shape.



Canopy Without Brackets

Black Canopy - Size 1 (8740110)	\$357
Black Canopy - Size 2 (8740120)	\$357

Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1)	\$102
Canopy Brackets with Push Handle for Size 2 (905039-2)	\$102

Vest and Belts

More options can be found at etac.com

Vests

Cross Vests



Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960)	\$253
Cross Vest Small - S-P: 10.25" - 13.5" (85720-S) (E0960)	\$253
Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960)	\$253
Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960)	\$253
Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960)	\$253

Torso Vests



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$268
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$268
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$268
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$268

Fixlocks Bracket



Higher Brackets (8910260)	\$96
<i>Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.</i>	

Wider Brackets (905068)	\$71
-------------------------	------



Pelvic Belts and Harnesses

4-Point Belts



4-Point Padded Hip Belts 1" X-Small - H: 6.75" - 12.5" (85725-XS) (E0960)	\$179
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$179
4-Point Padded Hip Belts 1" Medium - H: 9.75" - 16.25" (85725-M) (E0960)	\$179
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$179

2-Point Pelvic Belt



2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$130
---	-------

Pelvic Harness



Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0)	\$326
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1)	\$326
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2)	\$326

Pelvis Belt



Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$326
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-XS) (E0978)	\$326
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$326



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Lateral Support Cushion and Fittings

Lateral Supports

If choosing Laterals, both the Left and Right cushions must be selected.

Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	N/C
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	N/C
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	N/C

Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	N/C
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	N/C
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	N/C

Type of Fitting

Selection required, fittings are sold as a pair.

Fixed Fitting (E1028)	\$212
Swing-Away and Angle Adjustable Fitting (E1028)	\$332

Footrests for Mobility Bases

Footrest attached the mobility bases. If mobility base are desired with a footrest make a selection here and select desired mobility base below.

Footrests

Select footrest type and desired mounting bar

Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5" (95145-1) (K0040)	\$351
Angle and Height Adjustable Footplate, Medium, 14" x 9" (95145-2)(K0040)	\$351
Angle and Height Adjustable Footplate, Large, 16" x 9.5" (95145-3) (K0040)	\$351

Footrest - Multi-Adjustable

Mounting hardware included. Do not quote mounting bar. Not available with Foot Box,

Multi-Adjust Foot Support, Small - 6" - 10.5" (8910852-1)	\$539
Multi-Adjust Foot Support, Medium - 6" - 14.5" (8910852-2)	\$539
Multi-Adjust Foot Support, Large - 6" - 17.5" (8910852-3)	\$539

Mounting Bar

Required for Angle and Height Adjustable Footplates

Short Mounting Bar - 6-10.5" (8910821-1)	N/C
Medium Mounting Bar - 6-14.5" (8910821-2)	N/C
Long Mounting Bar - 6-17.5" (8910821-3)	N/C

Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjust	\$102
--------------------------------	-------

Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$173
Neoprene Sandals - Medium (85696-2)	\$173
Formed Shoeholder - Small (86147)	\$220
Formed Shoeholder - Medium (86148)	\$220
Footstraps w/ fix lock (85096) (E0951)	\$49
Ankle Straps, X-Small - C: 6.25"-9" (85745-XS)	\$247
Ankle Straps, Medium - C: 9"-12.25" (85745-M)	\$247
Ankle Straps, X-Large - C: 11.5"-15" (85745-XL)	\$247

Foot Box Kits

Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi Adjustable Footrests.

Foot Box Kit

Foot Box Kit, Size 1 Plate (945145-1 BOX)	\$336
Foot Box Kit, Size 2 Plate (945145-2 BOX)	\$357
Foot Box Kit, Size 3 Plate (945145-3 BOX)	\$385

Strap for foot box

Required selection

Strap for Foot Box (95145-20)	\$137
-------------------------------	-------

Foot Plate Padding

Required selection

Foot Plate Padding, Size 1 (95145-1PAD)	\$168
Foot Plate Padding, Size 2 (95145-2PAD)	\$231
Foot Plate Padding, Size 3 (95145-3PAD)	\$287

Padding for foot plate rod

Required selection

Padding for Foot Plate Rod (95145-150)	\$90
--	------



Attaches to mounting bar

Calf Support for Centrally Mounted Foot Plate (85131)	\$343
---	-------



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com



x:panda shape
advanced

Chilla Mobility Base

Chilla Frame

Only Available with Size 1 x:panda Shape

Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) \$3682

Shape Adapter

x:panda Shape adapter Plate (881091) \$1225

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) \$314

O2 Holder (881096) \$245

Under Storage Basket (881092) \$429

IV Pole (881097) \$665

Vent Tray (881095) \$462

High-low:xo Mobility Base and Options

High-low:xo

Manual Lifts

Manual Lift - indoor Size 1 (K0108) \$2422

Manual Lift - indoor, Size 2 (K0108) \$2599

Outdoor High-low:xo

Outdoor- manual lift, Size 2 (89120152) \$4080

Outdoor- power lift, Size 2 \$5216

Power Lifts

Indoor- power lift, Size 1 (K0108) \$3882

Indoor- power lift, Indoor, Size 2 (K0108) \$4177

Anti-Tips on the High-low:xo

Anti-tips - Indoor (890007) \$273

Anti-tips - Outdoor (890008) \$293

Transit - Outdoor High-low:xo

Transit (890005) \$319

Push Brace for High-low:xo - Indoor

Required selection

Fixed Angle Push Brace (890009-04) \$242

Angle Adjustable Push Brace (8910555) \$512

Multi Frame Mobility Base and Options

Multi Frame Base Options

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel

Locks \$3888

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake \$4257

Anti- Tip Option

Anti-tips \$240

Transit Option

Transit for Multi Frame \$331

Multi-Frame Accessories

Adaptor for x:panda on Multi Frame - for x:panda sz 1 & 2 \$290

Vent Tray \$611

O2 Holder \$465

Ki Mobility Bases and Options

Ki Mobility Bases

Separate Order form for base must be filled out. Order form is located at Kimobility.com. Footrests are not provided with the base. If footrest are desired, selections must be made above in the "Mobility Base Footrest" section.

Tyke x:panda

Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Tyke for x:panda N/C

Little Wave Flip for x:panda

Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Flip for x:panda N/C