



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote: PO Number: \_\_\_\_\_ City: \_\_\_\_\_ Province: \_\_\_\_\_  
 Clinician / Sales Rep: \_\_\_\_\_ Postal Code: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_  
 Associated Quote/Order#: \_\_\_\_\_ Effective: 11/14/2023

	Infant	1	2
Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

## Base

Please see measurement guide for growth information. A selection must be made in each section below.

### Advanced Seat Size

Includes base and plate

Size 1 Advanced Seat <i>See chart above for size details (SE0001947)</i>	\$2795
Size 2 Advanced Seat <i>See chart above for size details (SE0001947)</i>	\$2795

### Back Recline

Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless	N/C
Adjustable Stepless Recline with Adjustable Dynamic, Lock of Dynamic, Toolless	N/C

## Upper Back

*Choose cushions under Cushions. Includes quick release detachment (#9050374)*

### Size 1 Base

Size 1 Base Upper - W: 8.75" x H: 7" <i>Total back height range 12.25" - 16.25"</i>	N/C
Size 1 Base Upper - W: 10.25" x H: 8.75" <i>Total back height range 13.75" - 18.50"</i>	N/C



### Size 2 Base

Size 2 Base Upper - W: 8.75" x H: 7" <i>Total back height range 13.5" - 17.25"</i>	N/C
Size 2 Base Upper - W: 10.25" x H: 8.75" <i>Total back height range 15" - 19.75"</i>	N/C
Size 2 Base Upper - W: 11.75" x H: 9.5" <i>Total back height range 17" - 21.75"</i>	N/C

## Daily Adjustment Kit - Upper Back

Daily Adjustment Kit - Upper Back <i>Allows for tool-less adjustment of back height. Ideal for frequent adjustment.</i>	\$25
--	------



## Spacer Blocks

*Pad between back and upper back. Used when back height is increased. Not a required selection. Field retrofitable.*

### Spacer Block Holder without Cushion

Spacer Block Holder, Size 1 - Back widths 7", 8.75", 10.25"	\$219
Spacer Block Holder, Size 2 - Back widths 8.75", 10.25", 11.75", 13.5"	\$219
Spacer Block Holder, Size 3 - Back width 11.75", 13.5", 15", 16.5"	\$219

### Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25"	N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25"	N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25"	N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25"	N/C
Spacer Block Cushion Size 5 - W: 13.5" x H: 2.25"	N/C

### Spacer Block Cushion Cover

Black	N/C
Black Wipeable	N/C

## Cushions

### Seat Cushion

*Required selections*

#### Color

Red	N/C
Grey	N/C
Black Wipeable	\$25

### Seat Cushion - Split

*Advanced seat requires split cushion*

Split Cushion, Size 1 - W: 10.25" x L: 19.75"	N/C
Split Cushion, Size 1 - W: 11.75" x L: 19.75"	N/C
Split Cushion, Size 2 - W: 12.5" x L: 25.5"	N/C
Split Cushion, Size 2 - W: 15.75" x L: 25.5"	N/C

### Lower Back Cushion

*The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.*

#### Lower Back Color

Red	N/C
Grey	N/C
Black Wipeable	N/C

#### Size 1 - Lower Back Cushion Size

*Select if using Size 1 Seat*

Size 1 Lower Back Cushion - W: 8.75" x H: 4.75"	N/C
Size 1 Lower Back Cushion - W: 10.25" x H: 4.75"	N/C
Size 1 Lower Back Cushion - W: 11.75" x H: 4.75"	N/C

#### Size 2 - Lower Back Cushion Size

*Select if using Size 2 Seat*

Size 2 Lower Back Cushion - W: 11" x H: 5.5"	N/C
Size 2 Lower Back Cushion - W: 12.5" x H: 5.5"	N/C
Size 2 Lower Back Cushion - W: 14.25" x H: 5.5"	N/C
Size 2 Lower Back Cushion - W: 15.75" x H: 5.5"	N/C

### Upper Back Cushion

*Required selections*

#### Upper Back Cushion Cover Color

Red	N/C
Grey	N/C
Black Wipeable	N/C

#### Upper Back Cushion

Upper Back Cushion - W: 7" x H: 5.5"	N/C
Upper Back Cushion - W: 7" x H: 7"	N/C
Upper Back Cushion - W: 8.75" x H: 7"	N/C
Upper Back Cushion - W: 10.25" x H: 7"	N/C
Upper Back Cushion - W: 8.75" x H: 8.75"	N/C
Upper Back Cushion - W: 10.25" x H: 8.75"	N/C
Upper Back Cushion - W: 11.75" x H: 8.75"	N/C
Upper Back Cushion - W: 13.5" x H: 8.75"	N/C
Upper Back Cushion - W: 10.25" x H: 10.75"	N/C
Upper Back Cushion - W: 13.5" x H: 10.75"	N/C
Upper Back Cushion - W: 15" x H: 10.75"	N/C

## Headrest

### Headrest Style

*Includes Cushion*

Curved, Small - Black Wipeable - W: 7", H: 5"	\$269
Curved, Small - Black - W: 7", H: 5"	\$269
Curved, Medium - Black Wipeable	\$294
Curved, Medium - Black	\$294
Anatomic, Small - Black Wipeable	\$332
Anatomic, Small - Black	\$332
Anatomic, Medium - Black Wipeable	\$344
Anatomic, Medium - Black	\$344
Anatomic, Large - Black Wipeable	\$357
Anatomic, Large - Black	\$357
Adjustable with Large Laterals - Black	\$398
Adjustable with Large Laterals and Speakers - Black	\$528
Adjustable with Four Laterals - Black	\$494

### Mounting Hardware for Curved and Anatomic Headrest Style

*Required with Curved and Anatomic Head Support*

Mounting Brackets, Standard <i>5.5" of forward/backward adjustment</i>	\$300
Mounting Brackets, Extended <i>8.75" of forward/backward adjustment</i>	\$313
Mounting Brackets, 3D Adjustment <i>8.75" of forward/backward adjustment</i>	\$325

### Mounting Hardware for Adjustable Headrest Style

*Required with Adjustable Head Support*

Short Mounting Hardware <i>Suitable for upper back heights of 7", 8.75", 10.75" and 13.5"</i>	\$300
Long Mounting Hardware <i>Suitable for upper back heights of 10.75" and 13.5"</i>	\$300

### Headrest Options

#### Daily Adjustment Kits

*Creates toolless adjustability. For angle adjustment of head support.*

Daily Adjustment Kit Standard - Set of Two	\$38
Daily Adjustment Kit Extended - Set of Three	\$57



#### Additional Headrest Options

Headband System <i>For use with Anatomic Headrest</i>	\$300
Universal Headrest Mount <i>Required if using a third party headrest</i>	\$80

Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

## Arm Supports

### Arm Supports Size

*Cover not included, select separately*

Size 2 - W: 3" L: 10.25" H: 4" - 6.5" (SEAND0005)	\$407
Size 3 - W: 3" L: 12.25" H: 4" - 8.75" (SEAND0005)	\$407



### Arm Support Cover

Black	\$94
Black Wipeable	\$94

### Daily Adjustment Kit

*For toolless height adjustment of arm supports*

Height Adjustment Kit of Arm Supports	\$25
---------------------------------------	------



### Upper Arm Support - Fixed

*Includes Black Wipeable Cover*

Size 1 - L: 5.5" H: 7" For size 1 Base. (SEAND0005)	\$313
Size 2 - L: 7.5" H: 10.75" For size 2 Base. (SEAND0005)	\$313



### Upper Arm Supports - Adjustable

*Includes Black Wipeable Cover*

Size 1 - L: 5.5" H: 7" For size 1 Base. (SEAND0005)	\$419
Size 2 - L: 7.5" H: 10.75" For size 2 Base. (SEAND0005)	\$419

## Abduction and Adduction

### Inner Knee Support

*Includes cover. Used for controlling adduction.*

For Size 1 and 2 Black - L: 7.5", H: 7.5" _____	\$443
For Size 1 and 2 Black Wipeable - L: 7.5", H: 7.5" _____	\$443
For Size 2 Black - L: 9.5", H: 9.5" _____	\$443
For Size 2 Black Wipeable - L: 9.5", H: 9.5" _____	\$443



### Quick Release Pommel

*Includes cover*

For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" _____ (SESND1043)	\$438
For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" _____ (SESND1043)	\$438
For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" _____ (SESND1043)	\$438
For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" _____ (SESND1043)	\$438
For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" _____ (SESND1043)	\$438
For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" _____ (SESND1043)	\$438



### Outer Knee Support

*Includes cover. Used for controlling abduction.*

Standard Set for Size 1, Black - L: 6", H: 3.25" _____	\$130
Standard Set for Size 1, Black Wipeable - L: 6", H: 3.25" _____	\$130
Standard Set for Size 2, Black - L: 7.75", H: 3.5" _____	\$130
Standard Set for Size 2, Black Wipeable - L: 7.75", H: 3.5" _____	\$130
Flip Down Set for Size 2, Black - L: 7.75", H: 7.75" _____	\$130
Flip Down Set for Size 2, Black Wipeable - L: 7.75", H: 7.75" _____	\$130



## Accessories

### Accessory Holders and Brackets

#### Holder for Remote Control - Arm Support Mounted

*For use with powered base*

Remote Control Holder for Size 1 and 2 - Arm Support Mounted _____	\$82
Remote Control Holder for Size 3 - Arm Support Mounted _____	\$82



#### Holder For Remote - Back Mounted

Remote Control Holder - Back Mounted _____	\$157
--	-------



#### Bracket for Communication Aid - Arm Support Mounted

Bracket for Communication Aid, Size 1 and 2 - Arm Support Mounted _____	\$82
Bracket for Communication Aid, Size 3 - Arm Support Mounted _____	\$82



### Canopy

#### Canopy Without Brackets

*Push handle and brackets ordered separately. Not compatible with Size 1 Tyke Base.*

Black Canopy - Size 1 _____ <i>Requires upper back H: 5.5", 3.25"</i>	\$350
Black Canopy - Size 2 _____ <i>Requires upper back H: 8.75", 10.50", 13.50"</i>	\$350



#### Canopy Brackets

Canopy Brackets with Push Handle for Size 1 _____	\$100
Canopy Brackets with Push Handle for Size 2 _____	\$100



**Postural Supports**

**Width Adjustable Pelvis Supports**

*Includes cover, mounted from the back*

- Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" \_\_\_\_\_ \$294  
*Generally used on Size 1, compatible with Size 2*
- Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" \_\_\_\_\_ \$294
- Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" \_\_\_\_\_ \$294  
*Generally used on size 2, compatible on Size 1 with 6.25" hip width*
- Width Adjustable Pelvis Supports, Large Black Wipeable - L: 5.5", H: 5" \_\_\_\_\_ \$294



**Width Adjustable Hip and Pelvis Supports**

*Includes cover, mounted from the sides*

- Width Adjustable Hip Supports, Small Black - L: 4" H: 4.25" \_\_\_\_\_ \$419
- Width Adjustable Hip Supports, Small Black Wipeable - L: 4" H: 4.25" \_\_\_\_\_ \$419
- Width Adjustable Hip Supports, Large Black - L: 5.5" H: 5" \_\_\_\_\_ \$419
- Width Adjustable Hip Supports, Large Black Wipeable - L: 5.5", H: 5" \_\_\_\_\_ \$419



**Adjustable Anti-Thrust Wedges - For Size 1**

- Short - W: 4.5", L: 2.75-4.25", H: .75" \_\_\_\_\_ \$138
- Medium - W: 5.25", L: 4.25-6.75", H: 1.25" \_\_\_\_\_ \$138
- Long - W: 5.25", L: 4.75-7.75", H: 1.75" \_\_\_\_\_ \$138

**Adjustable Anti-Thrust Wedges - For Size 2**

*Can be used with Size 1 seat widths of 11.75"*

- Medium - W: 6", L: 5 -8.75", H: 1.25" \_\_\_\_\_ \$138
- Long - W: 6", L: 5.50-9.5", H: 2" \_\_\_\_\_ \$138

**Trays and Activity Centers**

**Tray**

- Shape 5 Tray with Edge Size 1 \_\_\_\_\_ \$350  
*Double rod, grey ABS*
- Shape 6 Tray with Edge Size 2 \_\_\_\_\_ \$350  
*Double rod, grey ABS. Fits size 2 width 11", 12.5", and 14.25".*
- Shape 8 Tray with Edge Size 1 \_\_\_\_\_ \$350



**Padding for Tray**

- Padding for Size 1 Tray, Shape 5 Plate, Over and Under \_\_\_\_\_ \$450
- Padding for Size 1 Tray, Top, Shape 5 Plate, Over \_\_\_\_\_ \$320
- Padding for Size 2 Tray, Shape 6 Plate, Over and Under \_\_\_\_\_ \$450
- Padding for Size 2 Tray, Top, Shape 6 Plate, Over \_\_\_\_\_ \$320



**Toy Bar**

- Toy Bar \_\_\_\_\_ \$300



**Arm Pads for Tray**

- Arm Pads for Tray \_\_\_\_\_ \$378  
*Mounts on top of tray*

**Vests and Belts**

**Fixlocks**

Higher Brackets <i>Brackets for higher strap guidance. Includes set of fixlocks. Requires Wider Brackets below.</i>	\$94
Wider Brackets	\$69



**Torso Vests**

Torso Vest Small - S-P: 10.75" - 12.5" <i>(SERND0020)</i>	\$262
Torso Vest Medium - S-P: 12.5" - 16" <i>(SERND0020)</i>	\$262
Torso Vest Large - S-P: 16" - 18.5" <i>(SERND0020)</i>	\$262
Torso Vest X-Large - S-P: 18.5" - 20.5" <i>(SERND0020)</i>	\$262



**Cross Vests**

Cross Vest X-Small - S-P: 7.75" - 12.5" <i>(SERND0020)</i>	\$248
Cross Vest Small - S-P 10.25" - 13.5" <i>(SERND0020)</i>	\$248
Cross Vest Medium - S-P: 13.5" - 17.75" <i>(SERND0020)</i>	\$248
Cross Vest Large - S-P: 17.75" - 20.75" <i>(SERND0020)</i>	\$248
Cross Vest X-Large - S-P: 20.75" - 24.5" <i>(SERND0020)</i>	\$248



**4-Point Belts**

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" <i>(SERND0020)</i>	\$175
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" <i>(SERND0020)</i>	\$175
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" <i>(SERND0020)</i>	\$175
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" <i>(SERND0020)</i>	\$175

**2-Point Pelvic Belt, Black**

2-Point Pelvic Belt Black - L: 13-14.5" <i>(SERND0030)</i>	\$127
---	-------

**Pelvic Harness**

Pelvic Harness, 0 - A:6.5" x B:7.5" <i>(SERND0020)</i>	\$319
Pelvic Harness, 1 - A:8.5" x B:8.75" <i>(SERND0020)</i>	\$319
Pelvic Harness, 2 - A:10.5" x B:9.75" <i>(SERND0020)</i>	\$319

**Pelvis Belt**

Pelvis Belt, X-Small - H-H: 61/4"- 73/4" (E0978)	\$319
Pelvis Belt, Small - H-H: 73/4"- 91/2" (E0978)	\$319
Pelvis Belt, Medium - H-H: 91/2"- 11" (E0978)	\$319



## Lateral Supports

### Type of Fitting

Fixed Fitting (SERND0020)	\$207
Swing-Away and Angle Adjustable Fitting (SELND0010)	\$325

### Left Cushion

*Provided with Black Wipeable cover*

Left Rectangular Small - L: 4.25" x H: 3.25"	N/C
Left Rectangular Medium - L: 6" x H: 3.5"	N/C
Left Rectangular Large - L: 7.5" x H: 3.5"	N/C

### Right Cushion

*Provided with Black Wipeable cover*

Right Rectangular Small - L: 4.25" x H: 3.25"	N/C
Right Rectangular Medium - L: 6" x H: 3.5"	N/C
Right Rectangular Large - L: 7.5" x H: 3.5"	N/C

## Mobility Bases and Footrests

### High-low:xo

#### Manual Lifts

High-low:xo Frame Manual Lift, Indoor, Size 1	\$2374
High-low:xo Frame Manual Lift, Indoor, Size 2	\$2548

#### Outdoor High-low:xo

High-low:xo Frame Outdoor, Size 1, Manual Lift - Size 1 and 2	\$4000
High-low:xo Frame Outdoor, Size 1, Power Lift - Size 1 and 2	\$5113

#### Power Lifts

High-low:xo Frame Power Lift, Indoor, Size 1	\$3805
High-low:xo Frame Power Lift, Indoor, Size 2	\$4095

#### Transit - Outdoor High-low:xo

Transit	\$312
---------	-------

#### Push Brace for High-low:xo - Indoor

Fixed Angle Push Brace	\$237
Angle Adjustable Push Brace	\$501

#### Anti-Tips on the High-low:xo

Anti-tips - Indoor	\$267
Anti-tips - Outdoor	\$287

### Footrests for High-low:xo

#### Footrest - Angle and Height Adjustable

Angle and Height Adjustable Footplate, Small, 12" x 7.5"	\$344
Angle and Height Adjustable Footplate, Medium, 14" x 9"	\$344
Angle and Height Adjustable Footplate, Large, 16" x 9.5"	\$344

#### Neoprene Sandals

Neoprene Sandals - Small	\$169
Neoprene Sandals - Medium	\$169

#### Mounting Bar

*Required selection for Angle and Height Adjustable Footrests*

Short Mounting Bar - 6-10.5"	N/C
Medium Mounting Bar - 6-14.5"	N/C
Long Mounting Bar - 6-17.5"	N/C

#### Formed Shoeholder

Formed Shoeholder - Small	\$215
Formed Shoeholder - Medium	\$215

#### Footstraps w/ fix lock

Footstraps w/ fix lock	\$48
------------------------	------

#### Footrest - Multi-Adjustable

Multi-Adjust Foot Support, Small - 6" - 10.5"	\$528
Multi-Adjust Foot Support, Medium - 6" - 14.5"	\$528
Multi-Adjust Foot Support, Large - 6" - 17.5"	\$528

#### Ankle Straps

Ankle Straps, X-Small - C: 6.25"-9"	\$242
Ankle Straps, Medium - C: 9"-12.25"	\$242
Ankle Straps, X-Large - C: 11.5"-15"	\$242

#### Width Extensions

Width Extension - Multi Adjust	\$100
--------------------------------	-------

### Stingray Frame

#### Stingray Tilt in Space

Stingray Tilt-In-Space Frame for x:panda, Size 1	\$3853
Stingray Tilt-In-Space Frame for x:panda, Central Wheel Locks, Size 1	\$4433

#### Stingray Frame Accessories

Transit	\$261
Vertical O2 Tank Holder	\$455
Medical Essential Carrier	\$127
Vent Tray	\$600
IV Pole, Height Adjustable and Collapsible	\$262
IV Pole Holder	\$183

### Footrests for Stingray

#### Footrest - Angle and Height Adjustable

Angle and Height Adjustable Footplate, Small, 12" x 7.5"	\$344
Angle and Height Adjustable Footplate, Medium, 14" x 9"	\$344
Angle and Height Adjustable Footplate, Large, 16" x 9.5"	\$344

#### Neoprene Sandals

Neoprene Sandals - Small	\$169
Neoprene Sandals - Medium	\$169

#### Mounting Bar

*Required selection for Angle and Height Adjustable Footrests*

Short Mounting Bar - 6-10.5"	N/C
Medium Mounting Bar - 6-14.5"	N/C
Long Mounting Bar - 6-17.5"	N/C

#### Formed Shoeholder

Formed Shoeholder - Small	\$215
Formed Shoeholder - Medium	\$215

#### Footstraps w/ fix lock

Footstraps w/ fix lock	\$48
------------------------	------

#### Footrest - Multi-Adjustable

Multi-Adjust Foot Support, Small - 6" - 10.5"	\$528
Multi-Adjust Foot Support, Medium - 6" - 14.5"	\$528
Multi-Adjust Foot Support, Large - 6" - 17.5"	\$528

#### Ankle Straps

Ankle Straps, X-Small - C: 6.25"-9"	\$242
Ankle Straps, Medium - C: 9"-12.25"	\$242
Ankle Straps, X-Large - C: 11.5"-15"	\$242

#### Width Extensions

Width Extension - Multi Adjust	\$100
--------------------------------	-------

### Multi Frame

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel Locks	\$3811
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake	\$4173

#### Anti-tips

Anti-tips	\$235
-----------	-------

#### Transit Option

Transit for Multi Frame	\$324
-------------------------	-------

#### Adaptor for x:panda

Adaptor for x:panda on Multi Frame - for x:panda sz 1 & 2	\$284
---	-------

#### Vent Tray for Multi Frame

Vent Tray	\$599
-----------	-------

#### O2 Holder

O2 Holder	\$455
-----------	-------

### Footrests for Multi Frame

#### Footrest - Angle and Height Adjustable

Angle and Height Adjustable Footplate, Small, 12" x 7.5"	\$344
Angle and Height Adjustable Footplate, Medium, 14" x 9"	\$344
Angle and Height Adjustable Footplate, Large, 16" x 9.5"	\$344

#### Width Extensions

Width Extension - Multi Adjust	\$100
--------------------------------	-------

Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

**Mounting Bar**

*Required selection for Angle and Height Adjustable Footrests*

Short Mounting Bar - 6-10.5" _____	N/C
Medium Mounting Bar - 6-14.5" _____	N/C
Long Mounting Bar - 6-17.5" _____	N/C

**Footrest - Multi-Adjustable**

Multi-Adjust Foot Support, Small - 6" - 10.5" _____	\$528
Multi-Adjust Foot Support, Medium - 6" - 14.5" _____	\$528
Multi-Adjust Foot Support, Large - 6" - 17.5" _____	\$528

**Neoprene Sandals**

Neoprene Sandals - Small _____	\$169
Neoprene Sandals - Medium _____	\$169

**Formed Shoeholder**

Formed Shoeholder - Small _____	\$215
Formed Shoeholder - Medium _____	\$215

**Footstraps w/ fix lock**

Footstraps w/ fix lock _____	\$48
------------------------------	------

**Ki Mobility Bases**

*Separate order form for base must be filled out. Order form is located at [www.kimobility.com](http://www.kimobility.com). Footrests are not provided with base. If footrests are desired selection must be made below.*

Tyke for x:panda _____	
Flip for x:panda _____	

**Footrests for Ki Mobility Bases**

**Footrest - Angle and Height Adjustable**

Angle and Height Adjustable Footplate, Small, 12" x 7.5" _____	\$344
Angle and Height Adjustable Footplate, Medium, 14" x 9" _____	\$344
Angle and Height Adjustable Footplate, Large, 16" x 9.5" _____	\$344

**Mounting Bar**

*Required selection for Angle and Height Adjustable Footrests*

Short Mounting Bar - 6-10.5" _____	N/C
Medium Mounting Bar - 6-14.5" _____	N/C
Long Mounting Bar - 6-17.5" _____	N/C

**Footrest - Multi-Adjustable**

Multi-Adjust Foot Support, Small - 6" - 10.5" _____	\$528
Multi-Adjust Foot Support, Medium - 6" - 14.5" _____	\$528
Multi-Adjust Foot Support, Large - 6" - 17.5" _____	\$528

**Width Extensions**

Width Extension - Multi Adjust _____	\$100
--------------------------------------	-------

**Neoprene Sandals**

Neoprene Sandals - Small _____	\$169
Neoprene Sandals - Medium _____	\$169

**Formed Shoeholder**

Formed Shoeholder - Small _____	\$215
Formed Shoeholder - Medium _____	\$215

**Footstraps w/ fix lock**

Footstraps w/ fix lock _____	\$48
------------------------------	------

**Ankle Straps**

Ankle Straps, X-Small - C: 6.25"-9" _____	\$242
Ankle Straps, Medium - C: 9"-12.25" _____	\$242
Ankle Straps, X-Large - C: 11.5"-15" _____	\$242