



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com



x:panda shape - Infant Seat

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 3/12/2025

	Infant	1	2
Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

## Seat Bases

Please see measurement guide for growth information. A selection must be made in each section below.

### Set Base

Includes base and plate

Size 1 Infant Base

See chart above for size details

\$2545

### Recline Options

Selection Required

Adjustable Stepless Recline, Toolless (K0108)

N/C

Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (K0108)

N/C

### Back

Upper Back with Quick Release Mount

N/C

## Cushions and Covers

### Seat Cushion

#### Cushion Color

*Cushion including seat, back, seat pit, upper arm cushions. Order covers for side support belt separately. Required Selection.*

Pink	N/C
Light Blue	N/C
Black Wipeable	\$21

#### Extra Cushion and Cover

*An extra cushion and cover set, in addition to the installed cushion and cover selected above. Comes uninstalled.*

Pink	\$416
Light Blue	\$416
Black Wipeable	\$479

#### Extra Cover only

*Select for an extra cover (no cushion). Comes uninstalled. Optional selection.*

Pink	\$224
Light Blue	\$224

## Headrests and Hardware Options

### Headrest Style

#### Type of Head Support

*Includes Short Rod length. Includes Cushion.*

Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)	\$275
Curved, Small - Black - W: 7", H: 5" (E0955)	\$275
Anatomic, Small - Black, Wipeable (E0955)	\$339
Anatomic, Small - Black (E0955)	\$339

### Headrest Mounting Hardware

*2 pivot points, standard 80mm rod, height adjustment rod 905037*

Standard Mounting Hardware (E1028)	\$192
------------------------------------	-------



### Headrest Options

*Knob replaces bolt for tool-less adjustment of headrest angle*

Daily Adjustment Kit Standard - Set of Two	\$39
Daily Adjustment Kit Standard - Set of Three	\$59

## Vest, Straps and Belts



### Side Supports Belts

Infant Side Support Belt, Black Wipeable (E0978)	\$269
Infant Side Support Belt, Black (E0978)	\$269

### Torso Vest

Torso Vest X-Small - S-P: 7.75" - 10.75" (E0960)	\$268
Torso Vest Small - S-P: 10.75" - 12.5" (E0960)	\$268

### Side Support Extra Cover

Extra Side Support Cover, Black, Wipeable	\$122
Extra Side Support Cover, Black	\$122

### Cross Vest

Cross Vest X-Small - S-P: 7.75" - 12.5" (E0960)	\$253
Cross Vest Small - S-P: 10.75" - 13.5" (E0960)	\$253

### Padded Hip Belt

4-Point Padded Hip Belt, 1" X-Small - H-H: 6.75" - 12.5" (E0978)	\$179
4-Point Padded Hip Belt, 1" Small - H-H: 8.25" - 13.75" (E0978)	\$179
2-Point Pelvic Belt Black - L: 13-14.5" (E0978)	\$130

## Trays and Accessories

### Tray and Activities Centres



### Tray

Soft tray for Infant Seating, Black - W: 12.5" x L: 7.5"	\$192
--	-------

### Canopy

#### Canopy

Black Canopy - Size 1	\$357
-----------------------	-------

#### Canopy Brackets

Canopy Brackets with Push Handle for Size 1	\$102
---	-------

## Footrests for Mobility Bases

Footrest attached the mobility bases. If mobility base are desired with a footrest make a selection here and select desired mobility base below.

### Footrests

#### Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5"  
(95145-1)(K0040) \$351

#### Mounting Bar

Short Mounting Bar - 6-10.5" (8910821-1) N/C

### Footrest Accessories

Neoprene Sandals - Small \$173  
 Neoprene Sandals - Medium \$173  
 Footstraps w/ fix lock (E0951) \$49  
 Ankle Straps, X-Small - C: 6.25" - 9" \$247

## Chilla Mobility Base

### Chilla Frame

Only available with size 1 Shape

#### Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234 \$3682

#### Shape Adapter

x:panda Shape adapter Plate (881091) K0208 \$1225

### Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) \$314 Vent Tray (881095) \$462  
 IV Pole (881097) \$665 O2 Holder (881096) \$245

## Mobility Bases

### High-low:xo

#### Manual Lifts

Manual Lift - indoor Size 1 (K0108) \$2422

#### Anti-Tips on the High-low:xo

Anti-tips - Indoor \$273  
 Anti-tips - Outdoor \$293

#### Power Lifts

Indoor - power lift, Size 1 (K0108) \$3882

### Push Brace for High-low:xo - Indoor

Required for High-low:xo

Fixed Angle Push Brace \$242  
 Angle Adjustable Push Brace \$512

### Ki Mobility Bases

Separate order form for base must be filled out. Order form is located at [www.kimobility.com](http://www.kimobility.com). Footrests are not provided with base. If footrests are desired selection must be made below.

Tyke for x:panda N/C  
 Flip for x:panda N/C