



by Etac  
Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com



x:panda shape - Simple Seat

Company:		Ship To:	
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	Province:
Clinician / Sales Rep:		Postal Code:	Phone:
Marked For:		Email/Fax:	
Associated Quote/Order#:		Effective: 3/12/2025	

## x:panda Shape - Simple Seat

	Infant	1	2
Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

## Shape Base Options

Please see measurement guide for growth information. A selection must be made in each section below.

### Shape Base Options

Selection required

#### Seat Base

Size 1 Simple Seat (E2292)	\$2545
Size 2 Simple Seat (E2292)	\$2545

### Seat Plate

Selection required

#### Seat Plate Size 1

Size 1, Seat W: 10", Lower Back W: 8.75"	N/C
Size 1, Seat W: 11.75", Lower Back W: 11.75"	N/C

#### Seat Plate Size 2

Size 2, Seat W: 12.5", Lower Back W: 12.5"	N/C
Size 2, Seat W: 15.75", Lower Back W: 15.75"	N/C

### Recline Options

Selection required

#### Back Recline

Adjustable Stepless Recline, Toolless (K0108)	N/C
Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (K0108)	N/C

## Back and Cushions

### Upper Back Size

Select upper back base then select Upper and Lower back cushion and cover. Includes quick release detachment (#9050374)



#### Size 1 Base

Size 1 Base Upper - H: 7" (905073) (E2291) \_\_\_\_\_ N/C  
Height range: 12.25" - 16.25"  
Size 1 Base Upper - H: 8.75" (905074) (E2291) \_\_\_\_\_ N/C  
Height range: 13.75" - 18.50"

#### Size 2 Base

Size 2 Base Upper - H: 7" (905073) (E2291) \_\_\_\_\_ N/C  
Height range: 13.5" - 17.25"  
Size 2 Base Upper - H: 8.75" (905074) (E2291) \_\_\_\_\_ N/C  
Height range: 15" - 19.25"  
Size 2 Base Upper - H: 10.75" (905075) (E2291) \_\_\_\_\_ N/C  
Height range: 17" - 21.75"

### Upper Back Cushion and Cover

Upper back cushion height must match Upper Back Base height

#### Upper Back Cushion Size

Upper Back Cushion-W: 7"xH: 5.5" (902101) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 7"xH: 7" (902103) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 8.75"xH: 7" (902104) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 10.25"xH: 7" (902105) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 8.75"xH: 8.75" (902106) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 10.25"xH: 8.75" (902107) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 11.75"xH: 8.75" (902108) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 13.25"xH: 8.75" (902109) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 10.25"xH: 10.75" (902110) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 11.75"xH: 10.75" (902111) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 13.75"xH: 10.75" (902112) \_\_\_\_\_ N/C  
Upper Back Cushion-W: 15"xH: 10.75" (902113) \_\_\_\_\_ N/C

#### Upper Back Cushion Cover

Selected required

Red (-50) \_\_\_\_\_ N/C  
Grey (-60) \_\_\_\_\_ N/C  
Black Wipeable (-70) \_\_\_\_\_ N/C

### Lower Back Cushion

The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.

#### Lower Back Cushion

Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120) \_\_\_\_\_ N/C  
Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902121) \_\_\_\_\_ N/C  
Size 1 Lower Back Cushion - W: 11" x H: 5.5" (902122) \_\_\_\_\_ N/C  
Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123) \_\_\_\_\_ N/C  
Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124) \_\_\_\_\_ N/C  
Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125) \_\_\_\_\_ N/C  
Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126) \_\_\_\_\_ N/C

#### Lower Back Cover

Selection required

Red (-50) \_\_\_\_\_ N/C  
Grey (-60) \_\_\_\_\_ N/C  
Black Wipeable (-70) \_\_\_\_\_ N/C

### Daily Adjustment Kit - Upper Back

Allows for tool-less adjustable of back height. Ideal for frequent adjustment.



#### Daily Adjustment Kit

Daily Adjustment Kit - Upper Back (905081) \_\_\_\_\_ \$26

## Spacer Blocks

### Standard Spacer Block Holder without Cushion

*Cushion widths must match pacer block holder width*

Spacer Block Holder, Size 1 - Back widths 7", 8.75", 10.25" (903100-1)	\$224
Spacer Block Holder, Size 2 - Back widths 8.75", 10.25", 11.75", 13.5" (903100-2)	\$224

### Standard Spacer Block Cushion Cover

Black (-80)	N/C
Black Wipeable (-70)	N/C

### Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011)	N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012)	N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013)	N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014)	N/C
Spacer Block Cushion Size 5 - W: 13.25" x H: 2.25" (903015)	N/C



## Cushions and Covers

### Seat Cushions

#### Seat Cushion Size - One Piece

One Piece Cushion, Size 1 - W: 9.75" x L: 19.75" (901101)	N/C
One Piece Cushion, Size 1 - W: 11.75" x L: 19.75" (901103)	N/C
One Piece Cushion, Size 2 - W: 12.50" x L: 25.5" (901105)	N/C
One Piece Cushion, Size 2 - W: 15.75" x L: 25.5" (901107)	N/C

#### Seat Cushion Size - Split

*Split cushion require or use with pommel*

Split Cushion, Size 1 - W: 10.25" x L: 19.75" (901114)	N/C
Split Cushion, Size 1 - W: 11.75" x L: 19.75" (901115)	N/C
Split Cushion, Size 2 - W: 12.5" x L: 25.5" (901117)	N/C
Split Cushion, Size 2 - W: 15.75" x L: 25.5" (901119)	N/C

### Seat Cover

*Required selection for cover color*

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$26

## Headrest and Hardware Options

### Headrest and Hardware Options

*Includes Cushion*

#### Type of Head Support

Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)	\$275
Curved, Small - Black - W: 7", H: 5" (E0955)	\$275
Curved, Medium - Black Wipeable (E0955)	\$300
Curved, Medium - Black (E0955)	\$300
Anatomic, Small - Black Wipeable (E0955)	\$339
Anatomic, Small - Black (E0955)	\$339
Anatomic, Medium - Black Wipeable (E0955)	\$351
Anatomic, Medium - Black (E0955)	\$351
Anatomic, Large - Black Wipeable (E0955)	\$365
Anatomic, Large - Black (E0955)	\$365
Adjustable w/ Large Laterals - Black (E0955)	\$406
Adjustable with Large Laterals and Speakers - Black (E0955)	\$539
Adjustable with Four Laterals - Black (E0955)	\$504

#### Additional Headrest Options

Headband System	\$306
Universal Headrest Mount	\$82

### Mounting Hardware for Curved and Anatomic Headrest

#### Headrest Mounting Hardware

Mounting Brackets, Standard (E1028)	\$306
Mounting Brackets, Extended (E1028)	\$320
Mounting Brackets, 3D Adjustment (E1028)	\$332

### Mounting Hardware for Adjustable Headrest Style

Short Mounting Hardware (E1028)	\$306
Long Mounting Hardware (E1028)	\$306

### Headrest Options



#### Daily Adjustment Kits

*Required with Curved and Anatomic Head support*

Daily Adjustment Kit Standard - Set of Two	\$39
Daily Adjustment Kit Extended - Set of Three	\$59

## Arm Supports

### Arm Support Size

#### Arm Supports Size

*Cover not included, select separately*

Size 2 - W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$416
Size 3 - W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$416

#### Arm Supports Cover

Black (-80)	\$96
Black Wipeable (-70)	\$96

#### Daily Adjustment Kit

*For toolless height adjustment of arm supports*

Height Adjustment Kit of Arm Supports (8910615-11)	\$26
--	------



#### Upper Arm Support - Fixed

*Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.*

Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$320
For size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$320
For Size 2 Base	

#### Upper Arm Supports - Adjustable

*Not compatible with Swing-Away Lateral Supports*

Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$428
For Size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$428
For size 2 Base	

## Knee and Postural Supports

### Abduction and Adduction



#### Outer Knee Support

Standard Set for Size 1, Black - L: 6", H: 3.25" (905005-70) (E0956)	\$300
Standard Set for Size 1, Black Wipeable - L: 6", H: 3.25" (905005-80) (E0956)	\$300
Standard Set for Size 2, Black - L: 7.75", H: 3.5" (905006-80) (E0956)	\$300
Standard Set for Size 2, Black Wipeable - L: 7.75", H: 3.5" (905006-70)(E0956)	\$133

#### Quick Release Pommel

For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$447
For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$447
For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905082-80) (E0957)	\$447
For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905082-70) (E0957)	\$447
For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905078-80) (E0957)	\$447
For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905078-70)(E0957)	\$447

### Postural Supports

#### Width Adjustable Pelvis Supports

*Mounted from the back*


Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" (905041-80) (E0956)	\$300
Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" (905041-70)(E0956)	\$300
Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" (905042-80) (E0956)	\$300
Width Adjustable Pelvis Supports, Large Black Wipeable - L: 5.5", H: 5" (905042-70) (E0956)	\$300

#### Width Adjustable Hip and Pelvis Supports

*Mounted from the sides*


Width Adjustable Hip Supports, Small Black - L: 4" H: 4.25" (905044-80) (E0956)	\$428
Width Adjustable Hip Supports, Small Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$428
Width Adjustable Hip Supports, Large Black - L: 5.5" H: 5" (905045-80) (E0956)	\$428
Width Adjustable Hip Supports, Large Black Wipeable - L: 5.5", H: 5" (905045-70) (E0956)	\$428

## Trays and Accessories

### Trays and Activity Centers

#### Tray



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950)	\$357
For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.	
Shape 6 Tray with Edge Size 2 (8910725-2) (E0950)	\$357
For Size 2. Double rod, Grey ABS.	
Shape 8 Tray with Edge Size 1 (8976350) (E0950)	\$357



Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028)	\$459
Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028)	\$327
Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028)	\$459
Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028)	\$327

## Canopy and Mounting Hardware

*Canopy does not come with brackets. Both the canopy and bracket must be quoted when mounting to the shape.*



### Canopy Without Brackets

Black Canopy - Size 1 (8740110)	\$357
Black Canopy - Size 2 (8740120)	\$357



## Additional Tray Accessories

Toy Bar (8910740) <i>For shape 5 and 6</i>	\$306
---	-------

*For Shape 5 and 6*

Arm Pads for Tray (8910730)	\$386
-----------------------------	-------



### Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1)	\$102
Canopy Brackets with Push Handle for Size 2 (905039-2)	\$102

## Accessory Holders and Brackets

Remote Control Holder Size 1 - Arm Support Mounted (905030-1)	\$84
Remote Control Holder Size 2 - Arm Support Mounted (905030-2)	\$84
Remote Control Holder - Back Mounted (905030)	\$161
Bracket for Communication Aid - Arm Support Mounted (905033-1)	\$84

## Vests and Belts

*More options can be found at etac.com*

### Vests

#### Torso Vests



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$268
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$268
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$268
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$268

#### Cross Vests



Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960)	\$253
Cross Vest Small - S-P 10.25" - 13.5" (85720-S) (E0960)	\$253
Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960)	\$253
Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960)	\$253
Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960)	\$253

### Fixlocks Bracket



Higher Brackets (8910260) <i>Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.</i>	\$96
Wider Brackets (905068)	\$71

## Pelvic Belts and Harnesses

### 4-Point Belts

4-Point Padded Hip Belts 1" X-Small - H: 6.75" - 12.5" (85725-XS) (E0960)	\$179
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$179
4-Point Padded Hip Belts 1" Medium - H: 9.75" - 16.25" (85725-M) (E0960)	\$179
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$179

### Pelvic Harness

Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0)	\$326
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1)	\$326
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2)	\$326

### Pelvis Belt

Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$326
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-XS) (E0978)	\$326
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$326

### 2-Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$130
---	-------

## Lateral Support Cushion and Fittings

*if choosing Laterals, both the left and right cushions must be selected*

### Left Lateral Supports

#### Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	N/C
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	N/C
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	N/C

### Right Lateral Supports

#### Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	N/C
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	N/C
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	N/C

### Fitting for Lateral Supports

*Selection required, fittings are sold as a pair*

#### Type of Fitting

Fixed Fitting (E1028)	\$212
Swing-Away and Angle Adjustable Fitting (E1028)	\$332

## Footrests for Mobility Bases

*Footrests attaches to the mobility bases. If mobility base are desired with a footrest, make a selection here and then choose mobility base below.*

### Footrest

*Select footrest type and desired mounting bar*

#### Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5" (95145-1) (K0040)	\$351
Angle and Height Adjustable Footplate, Medium, 14" x 9" (95145-2) (K0040)	\$351
Angle and Height Adjustable Footplate, Large, 16" x 9.5" (95145-3) (K0040)	\$351

#### Mounting Bar

Short Mounting Bar - 6-10.5"	N/C
Medium Mounting Bar - 6-14.5"	N/C
Long Mounting Bar - 6-17.5"	N/C

#### Width Extensions

*Available only with Multi-Adjustable Foot Supports*

Width Extension - Multi Adjust	\$102
--------------------------------	-------

#### Footrest - Multi-Adjustable

*Mounting hardware comes standard. Not compatible with the Footbox.*

Multi-Adjust Foot Support, Small - 6" - 10.5" (8910852-1)	\$539
Multi-Adjust Foot Support, Medium - 6" - 14.5" (8910852-2)	\$539
Multi-Adjust Foot Support, Large - 6" - 17.5" (8910852-3)	\$539

### Foot Box Kits

*Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi-Adjustable Footrest.*

#### Foot Box Kit

Foot Box Kit, Size 1 Plate (945145-1 BOX)	\$336
Foot Box Kit, Size 2 Plate (945145-2 BOX)	\$357
Foot Box Kit, Size 3 Plate (945145-3 BOX)	\$385

#### Strap for foot box

*Required with all sizes*

Strap for Foot Box (95145-20)	\$137
-------------------------------	-------

*Required selection*

Foot Plate Padding, Size 1 (95145-1PAD)	\$168
Foot Plate Padding, Size 2 (95145-2PAD)	\$231
Foot Plate Padding, Size 3 (95145-3PAD)	\$287

#### Padding for foot plate rod

*Selection required*

Padding for Foot Plate Rod (95145-150)	\$90
--	------

### Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$173
Neoprene Sandals - Medium (85696-2)	\$173
Formed Shoeholder - Small (86147)	\$220
Formed Shoeholder - Medium (86148)	\$220
Footstraps w/ fix lock (85096) (E0951)	\$49
Ankle Straps, X-Small - C: 6.25"-9" (85745-XS)	\$247
Ankle Straps, Medium - C: 9"-12.25" (85745-M)	\$247
Ankle Straps, X-Large - C: 11.5"-15" (85745-XL)	\$247



*Attaches to Mounting Bar*

Calf Support for Centrally Mounted Foot Plate (85131)	\$343
---	-------

## Chilla Mobility Base

### Chilla Frame

Only available with size 1 shape

### Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) \$3682

### Shape Adapter

x:panda Shape adapter Plate (881091) \$1225

### Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) \$314 Vent Tray (881095) \$462

IV Pole (881097) \$665 O2 Holder (881096) \$245

## High-low:xo Mobility Base and Options

### High-low:xo

#### Manual Lifts

High-low:xo Frame Manual Lift, Indoor, Size 1 \$2422  
 High-low:xo Frame Manual Lift, Indoor, Size 2 \$2599

#### Power Lifts

High-low:xo Frame Power Lift, Indoor, Size 1 \$3882  
 High-low:xo Frame Power Lift, Indoor, Size 2 \$4177

#### Outdoor High-low:xo

High-low:xo Frame Outdoor, Size 2, Manual Lift \$4080  
 High-low:xo Frame Outdoor, Size 2, Power Lift \$5216

#### Anti-Tips on the High-low:xo

Anti-tips - Indoor \$273  
 Anti-tips - Outdoor \$293

#### Transit - Outdoor High-low:xo

Transit \$319

### Push Brace for High-low:xo - Indoor

Required selection for High-Low base

Fixed Angle Push Brace \$242  
 Angle Adjustable Push Brace \$512

## Multi Frame Mobility Base and Options

### Multi Frame

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel  
 Locks \$3888  
 Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake \$4257  
 Adaptor for x:panda on Multi Frame - for x:panda sz 1 & 2 \$290

#### Base options

Anti-tips \$240  
 Transit for Multi Frame \$331  
 Vent Tray \$611  
 O2 Holder \$465

## Ki Mobility Bases

Separate order form for base must be filled out. Order form located at Kimobility.com. Footrest are not provide with base. if footrests are desired, make a selection below.

### Ki Mobility Bases

Price for base calculated with Ki Mobility's order form

Flip for x:panda N/C  
 Tyke for x:panda N/C