



**etac®**

# Etac Rex

Wall-mounted toilet arm support

# Rex toilet arm support

Rex is a stable and safe arm support that gives you safety close to the toilet. It offers extra support for individuals who need it, promoting independence and confidence in daily routines.

## Support designed for you

Rex is a foldable, space-saving and reliable wall-mounted toilet arm support. It is available in three different lengths and in variants with or without supporting legs. Choose a model that suits you, your bathroom and your requirements. It can easily be fitted onto the wall and blends well into your bathroom environment.

## Suitable for all bathrooms

The clever design ensures that the arm supports take up very little room when folded up. They are close to the wall and do not get in the way. Furthermore, thanks to their angled shape, it is always easy to reach and grip onto them. They are also designed to guarantee a gentle and smooth folding down process.



Rex comes with a 5-year warranty and a 10-year service life.

### Rex toilet arm support

Rex 60 cm	81705010
Rex 70 cm	81705020
Rex 70 cm with supporting leg	81705030
Rex 85 cm	81705040
Rex 85 cm with supporting leg	81705050

### Materials

Aluminium, polyurethane, ABS plastic, polypropylene.  
Supporting legs: Aluminium, polyamide, TPE, stainless steel.

### Weight

60 cm: 2.7 kg (59 lbs)  
70 cm: 2.9 kg (64 lbs)  
85 cm: 3.3 kg (73 lbs)  
Supporting leg: 0.7 kg (15.4 lbs)

### Maximum user weight

Rex 60/70 cm: 135 kg (298 lbs)  
Rex 85 cm: 100 kg (220 lbs)  
With supporting legs: 150 kg (330 lbs)

### Colour

White with volcano grey details

### Maintenance



Clean the product with a solvent-free cleaning agent (pH level 5-9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.






## Features



### Comfortable grip

The triangular shape of the arm support provides a stable and natural grip.




### Easy to clean

With its clean, smooth surfaces the arm support is both stylish and easy to clean.



### Stable and safe

Top quality materials and construction provides security when it is most needed. Etac Rex is also available with support legs if extra stability is required.



### Smart design

The arm supports take up very little room when folded up, while angled shape makes them easy to reach, grip and fold down.



## Accessories



### Toilet paper holder

81706001

Toilet paper holder in stainless steel which is attached to the arm support



### Mounting set

81706004

Can be used for mounting in concrete, plywood or wood



### Storage basket

81706003

Storage basket in polyester/mesh which is attached to the arm support.



Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit [www.etac.com](http://www.etac.com)

**Etac**

+46 371 587 330  
info@etac.se  
[www.etac.com](http://www.etac.com)

**etac**<sup>®</sup>  
Creating Possibilities