



**etac®**

# Etac Swift

Shower stool/chair

# Swift

## shower stool & shower chair

For over 20 years, users around the world have relied on Etac Swift shower stools and chairs to help them feel safe and confident in the bathroom. Developed in cooperation with rehab professionals, industrial designers and users, Swift shower stools and chairs are designed to ensure wellbeing, support and comfort.

After 50 years of creating possibilities in the home healthcare and rehabilitation markets, we pride ourselves on a unique range of high-quality products. With its award-winning design and carefully crafted features, **the Swift Shower chair is a true Etac original!**



### Designed for safe and sustainable use

A truly great product will not only perform well for the primary user in its intended environment. It should offer practical functionality throughout the product lifetime. The Swift shower stool/chair comes in a convenient flatpack for easy storage, shipment and carrying. The product itself is very easy to assemble without the need for any tools.

At Etac, products are designed with safety, hygiene and sustainability in mind. The Swift shower stool/chair is thoroughly tested in accordance with European and international standards and withstands rigorous use and reconditioning.



*Swift comes with a 5-year warranty and 10-year service life.*

# Discover the benefits

With sturdy height-adjustable legs and detachable back and arm supports that gently click into place, the Swift shower stool/chair has been designed for safe use and easy adaptation.



Swift  
shower stool

Swift  
shower stool with arm supports

## Easily adapted

Detachable, modular back and arm supports that gently click into place without the need for tools.

## Sturdy grip

Without arm supports, the side openings enable a stable grip and offer excellent support.

## Smart drainage

Specially-designed perforations that allow excess water to quickly drain away.

## Non-slip surface

The surface structure ensures a safe, non-slip seat even when wet.

## Flexible seat

Reversible seat with a generous recess.

## Height adjustable

Sturdy legs with a convenient built-in height adjustment.

## Dependably stable

Responsive, non-slip ferrules that provide stability on any floor.



Swift  
shower chair

# Swift shower stool/chair

A shower chair with many options.

## One Swift - three options

- The basic version of Swift comes without arm supports. The side openings enable a stable grip and provide excellent support.
- By adding arm supports, you get extra support when sitting down and standing up.
- By adding both arm and back support, you get a Swift shower chair, for a comfortable and safe shower experience.

## Two height options

Swift shower stool/chair are available in two different sizes, both of which are fully adjustable. The lower model can be adjusted into four different height positions and the higher model can be adjusted into seven positions.

The height-adjustable function can also be used to tilt the Swift shower chair slightly forward. This is particularly useful for anyone who has difficulty getting up and down from a seated position.



Etac Swift	White	Grey	Lagoon green	Blue
- shower stool	81701610	81701410	81701710	81701510
- shower stool with arm supports	81701620	81701420	81701720	
- shower chair with back and arm supports	81701630	81701430	81701730	81701530
- shower stool Low		81701450		
- shower chair Low with back and arm supports		81701460		

### Colours



### Dimensions

Total height with back support: 80–95 cm (31½–37½")  
 Total height Swift Low with back support: 72.5–80 cm (28½–31½")  
 Total width, stool: 54 cm (21¼")  
 Total width, with arm support: 54.5 cm (21½")  
 Total depth: 50 cm (19¾")  
 Total depth Swift Low: 46.5 cm (18")  
 Width between arm support: 45 cm (17¾")  
 Seat height: 42–57 cm (16½–22½")  
 Seat height Swift Low: 34.5–42 cm (13½–16½")  
 Seat area: 54 x 41 cm (21¼ x 16¼")

### Weight

Swift shower stool: 3.1 kg (6.8 lbs)  
 Swift Low shower stool: 2.9 kg (6.4 lbs)  
 Swift shower chair: 4.7 kg (10.4 lbs)  
 Swift Low shower chair: 4.5 kg (8.8 lbs)

### Max. user weight

Max user weight: 160 kg (352 lbs)

### Materials

Seat: Polypropylene  
 Legs: Aluminium, polyamide  
 Ferrules: TPE  
 Seat pads: Polyethylene

### Maintenance



Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85 °C for 3 minutes.



## Features



### Height adjustable

It is easy to adjust the height. Just click and press. The oval shape prevents the telescopic legs from turning when making the adjustment. There are no loose parts to misplace.



### Dependably stable

Responsive, non-slip ferrules adapt to the floor and provide stability on the most uneven of surfaces.



### Easily adapted

Modular design with detachable back and arm supports that gently click into place without the need for tools. Easy to customise for different user needs.



### Ease of access

A reversible seat with a generous recess and a smooth surface provides good access and comfort.

## Colour variants

Swift shower stools and chairs are available in four different colours. Grey, lagoon green, blue and white. All to give the user great options and possibility to match the shower stool/chair with the rest of the bathroom.







## Accessories



### Back support

Detachable, modular back support that gently clicks into place without the need for tools.

Grey ..... 84006181  
Lagoon green..... 84006198  
Blue..... 84006183  
White ..... 84006196



### Pair of arm supports

Detachable, modular arm supports that gently click into place without the need for tools.

Grey ..... 84006184  
Lagoon green..... 84006199  
Blue..... 84006186  
White ..... 84006197



### Soap dish

Attaches to the side of the seat. Perfect place to store shower gel or shampoo. Also has a built-in shower holder.

Grey ..... 80209266



### Seat pad

The non-slip pads provide an instant warm feeling. Sitting is more comfortable and not as slippery. Easy to wipe down and simple to click into place.

Grey ..... 84005073



### Back support pad

The non-slip pads provide an instant warm feeling. Easy to wipe down and simple to click into place.

Grey ..... 84005072





Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit [www.etac.com](http://www.etac.com)

Etac



+46 371 587 330



[info@etac.se](mailto:info@etac.se)



[www.etac.com](http://www.etac.com)

 **etac**<sup>®</sup>  
Creating Possibilities