



Etac presents:

S.P.A. Collection

Safe. Personal. Assistance.

 **etac**[®]

Creating Possibilities



Lean on us. Stay yourself.

For most of us, there comes a time – due to age or some other reason – when we need a little extra support. You are not as supple as you used to be, and cannot move like you used to. You feel as if you have to act a certain way. This is not the case, however.

Etac S.P.A. Collection comprises *Rex toilet arm support*, *Flex grab rail* and *Relax shower seat*. All three have a unique function but work together to ensure your bathroom is as comfortable, attractive and safe as possible. They give you the chance to continue being yourself.



Etac Rex toilet arm support

Secure support – when and where you need it most.

Retaining your independence in the bathroom is an important part of continuing to be yourself. Etac Rex is a stable and safe arm support designed to create a secure environment around the toilet.

The modern design suits any bathroom, even those with limited space.

Support designed for you

Etac Rex is available in three different lengths and can be ordered with or without supporting legs. Choose a model that suits you, your bathroom and your requirements. Can easily be fitted onto the wall and blends well into your bathroom environment.

Suitable for all bathrooms

The clever design ensures that the arm supports take up very little room when folded up. They are close to the wall and do not get in the way. Furthermore, thanks to their angled shape, it is always easy to reach and grip onto them. They are also designed to guarantee a gentle and smooth folding down process.

Rex features

Comfortable grip

The triangular shape of the arm support provides a stable and natural grip.

Stable and safe


Top quality materials and construction provides security when it is most needed. Etac Rex is also available with support legs if extra stability is required.

Easy to clean

With its clean, smooth surfaces, the arm support is both stylish and easy to clean.

Lightweight

The arm supports take up very little room when folded up, while the angled shape makes them easy to reach, grip and fold down.





Etac Relax shower seat

**No unnecessary space.
No unnecessary worry.**

Taking a shower is a great experience and we think that everyone should be able to enjoy that daily luxury and continue to be themselves. Etac Relax is a wall-mounted shower seat that fits most bathrooms without affecting your safety or comfort.

It is always close at hand, but as it can be folded upwards it never takes up any unnecessary space.

Scandinavian design for the modern bathroom

Etac Relax has been designed to ensure a secure and enjoyable shower experience. Design, function and safety are extremely important to us, and this is presented here in the form of a discreet, wall-mounted shower seat.

Adapted to your needs

How do you move about in the shower, how much support do you need, and which model best meets your needs? Etac Relax comes in two colours and six different models, which have all been designed to create options specifically for you.

Relax features

Large, stable seat

Large and stable seat when folded down and extremely discreet when folded up.

Easy-grip arm support

The shape and material of the folding arm supports provide a safe and secure grip.

Functional design

The shape of the seat offers the option of varied sitting positions, resulting in improved access.

Easy to fit

As there are only three fixing points, installation is both quick and easy.



Etac **Flex** grab rail



Flexibility and freedom

Etac Flex grab rail is the perfect choice for those who need a little extra support in the bathroom. Thanks to its flexible modular function, it is easy to adjust the handle to suit both you and your bathroom environment.

Easy to attach and easy to grip

Attach with screws or adhesive. Choose the method that best suits you and your bathroom. The pure Scandinavian design blends in with other products and you can choose between two colours - white and grey.

The support you need with the design you want

The material of the grab rail and its well-designed shape and surface provide a stable grip even if your hands or the handles are wet.

Flex features

Easy to adjust

The modular system can easily be assembled into a length and shape that suits you and your bathroom.

Secure grip


The well-designed shape and surface of the grab rail provide a stable grip even if your hands are wet.

Creating Possibilities

You can connect as many grab rails as you want and choose the angles that suit you. This creates possibilities.

Attach in two ways

Attach with screws or adhesive. The latter is perfect for those who cannot or do not want to drill in the bathroom.



Etac Rex toilet arm support

Etac Rex is available in three different lengths and in variants with or without supporting legs. The accessories toilet paper holder and storage basket can be attached to the arm support.



Maximum user weight

Etac Rex 60/70 cm: 135 kg Etac

Rex 85 cm: 100 kg

Med stödben: 150 kg

Etac Relax shower seat

Relax shower seat comes in six different models. If extra support is needed, choose a variant with supporting legs.



Maximum user weight

Med/utan stödben: 150/125 kg

Etac Flex grab rail

Flex is easy to customize. Choose the colour, angle and length that fits your needs.



Maximum load

100kg



Read more

Visit www.etac.com for more information about the products..





Safe. Personal. Assistance.

Etac S.P.A. Collection comprises Rex toilet arm support, Flex grab rail and Relax shower seat. All three have a unique function but work together to ensure your bathroom is as comfortable, attractive and safe as possible. And to ensure you can continue being yourself.

Lean on us. Stay yourself.

Visit Etac.com for more information or get in touch and we will tell you more.

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 58 73 30
Fax +46 371 58 73 90
info@etac.se www.etac.com


Creating Possibilities