



Molift RgoSling MediumBack Plus

A comfortable sling with full-body support for users weighing from 160 kg.

There are several considerations when caring for a bariatric patient, such as weight, body shape and medical condition. It is also important to enable care with compassion and dignity.

RgoSling MediumBack Plus is designed for users weighing over 160 kg. The sling allows the lifting of heavy users in a safe and comfortable way.


RgoSling MediumBack Plus provides full body support. It has a split leg opening and six straps with fixed loops at four heights. If head support is needed, a separate headrest can be used in combination with the sling.

The sling is intended to be used in combination with a 4-point sling bar.



Full-body support

The sling is supporting the whole body.



Height options

Straps with fixed loops and multiple height options



Color-coded loops

Colour coded hoisting loops helps to find an individual fit.



Comfortable

Padded comfortable leg-support with dual straps



The lifting straps of the leg support and the back have multiple height options to be able to make an individual assessment depending on the patients needs. The leg support can be positioned either crossed as this photo shows or on the same side as shown on the front page.

Molift RgoSling MediumBack Plus

Item no.

L
XL

1721940
1721950



Material: Polyester, PE cellfoam
Maximum user weight: 500 kg

Washing at higher temperatures wears out the material faster. The sling must be inspected regularly, preferably prior to each lifting operation but especially after it has been washed. Make sure not to damage or remove labels when cleaning.

For the latest news and continuously updated product information – visit: www.etac.com