



HoverMatt® PROS Air Sling

Versatile patient handling – with wipeable surface for easy cleaning.

The HoverMatt PROS Air Sling combines air-assisted repositioning with sling functionality for safe and efficient patient handling. It enables smooth turning, boosting, transferring and lifting with less effort, helping reduce caregiver strain.

Breathable, wipeable materials minimise shear and friction while protecting sensitive skin. Used in conjunction with the HoverCover, it facilitates micro-climate management and enhances patient comfort.

Designed for use in ICUs and acute care, its single-patient-use design supports infection control and eliminates the need for laundering.



Versatile functionality

Supports turning, repositioning, transferring and lifting.




Efficient workflow

Reduces caregiver effort, promoting staff safety and efficiency.



Gentle movement

Helps protect skin and promote patient comfort.



Single-patient use

Minimises cross-contamination risk; no laundering required.



Product

HoverMatt PROS Air Sling incl. HoverCover

Pcs

5

Item no.

PROS-HMSL-CS



Maximum user weight:

Lift 273 kg

Air-assist 544 kg

Length:

1980 mm

Width:

1230 mm

Height:

180 mm

Material (Latex-free):

Nylon twill

Construction:

Sewn

Compatible with the HTAir Supply, item.no HTAIR-2300

If soiled, the HoverMatt PROS Air Sling may be wiped down with disinfectant wipes or a cleaning solution used by your hospital for medical equipment disinfection. A 10:1 bleach solution (10 parts water: one part bleach) may also be used. NOTE: Cleaning with bleach solution may discolor fabric. To help keep the PROS Air clean, HoverTech recommends the use of the HoverCover™ Disposable Absorbent Cover. Whatever the patient is lying on to keep the hospital bed clean may also be placed on top of the PROS Air. The single-patient use PROS Air Sling is not intended to be laundered.

For the intended purpose of the products and latest information, please visit www.etac.com

Etac

 +46 371 587 330

 info@etac.se

 www.etac.com


Creating Possibilities