



Etac Swift

Shower stool/chair

 **etac**[®]

Etac Swift shower stool/chair

Designed for safe, sustainable use

For the past 20 years, users around the world have relied on Etac Swift shower stools and chairs to help them feel safe and confident in the bathroom.

At Etac, products are designed with safety, hygiene and sustainability in mind. The Swift shower stool/chair is thoroughly tested in accordance with European and international standards and withstands rigorous use and reconditioning. Its estimated lifetime is 10 years and it comes with a 5-year warranty.

Adaptable support

The basic version of Swift comes without armrests. The side openings enable a stable grip and provide excellent support. By adding armrests, you get extra support when sitting down and standing up. For a complete Swift shower chair, add both armrest and backrest. Should extra comfort be needed, seat and back pads are available.

With its award-winning design and carefully crafted features, the Swift Shower chair is a true Etac original!



Etac Swift shower chair



Etac Swift shower stool



Etac Swift shower stool with armrests

Etac Swift	White	Grey	Lagoon green	Blue
- shower stool	81701610	81701410	81701710	81701510
- shower stool with armrests	81701620	81701420	81701720	
- shower chair with back and armrests	81701630	81701430	81701730	81701530
- shower stool Low		81701450		
- shower chair Low with back and armrests		81701460		

Colours



Dimensions

Total height with backrest: 80–95 cm (31½–37½")
 Total height Swift Low with backrest: 72.5–80 cm (28½–31½")
 Total width, stool: 54 cm (21¼")
 Total width, with armrests: 54.5 cm (21½")
 Total depth: 50 cm (19¾")
 Total depth Swift Low: 46.5 cm (18")
 Width between armrests: 45 cm (17¾")
 Seat height: 42–57 cm (16½–22½")
 Seat height Swift Low: 34.5–42 cm (13½–16½")
 Seat area: 54 x 41 cm (21¼ x 16¼")

Weight

Swift shower stool: 3.1 kg (6.8 lbs)
 Swift Low shower stool: 2.9 kg (6.4 lbs)
 Swift shower chair: 4.7 kg (10.4 lbs)
 Swift Low shower chair: 4.5 kg (8.8 lbs)

Materials

Seat: Polypropylene
 Legs: Aluminium, polyamide
 Ferrules: TPE
 Seat pads: Polyethylene

Max. user weight

Max user weight: 160 kg (352 lbs)

Maintenance



Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85 °C for 3 minutes.



Features



Height adjustable

It is easy to adjust the height. Just click and press. The oval shape prevents the telescopic legs from turning when making the adjustment. There are no loose parts to misplace.



Dependably stable

Responsive, non-slip ferrules adapt to the floor and provide stability on the most uneven of surfaces.



Easily adapted

Modular design with detachable back and armrests that gently click into place without the need for tools. Easy to customise for different user needs.



Ease of access

A reversible seat with a generous recess and a smooth surface provides good access and comfort.

Two heights and four colour options

Swift shower stools/chairs have a max user weight of 160 kg and are available in two different sizes, both of which are fully adjustable. The lower model can be adjusted into four different height positions and the higher model can be adjusted into seven positions.

The height-adjustable function can also be used to tilt the shower chair slightly forward. This is particularly useful for anyone who has difficulty getting up and down from a seated position. Swift shower stool/chairs are available in four colours, grey, blue, lagoon green and white.



Accessories



Backrest

Detachable, modular backrest that gently clicks into place without the need for tools.

Grey	84006181
Lagoon green.....	84006198
Blue.....	84006183
White.....	84006196



Pair of armrests

Detachable, modular armrests that gently click into place without the need for tools.

Grey	84006184
Lagoon green.....	84006199
Blue.....	84006186
White.....	84006197



Soap dish

Attaches to the side of the seat. Perfect place to store shower gel or shampoo. Also has a built-in shower holder.

Grey	80209266
------------	----------



Seat pad

The non-slip pads provide an instant warm feeling. Sitting is more comfortable and not as slippery. Easy to wipe down and simple to click into place.

Grey	84005073
------------	----------



Backrest pad

The non-slip pads provide an instant warm feeling. Easy to wipe down and simple to click into place.

Grey.....	84005072
-----------	----------



Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 - 58 73 30
info@etac.se www.etac.com


Creating Possibilities