

Molift Smart 150

Compact smartness - light, foldable and transportable without tools

Molift Smart 150 is designed to be easily wheeled and transported to where it is required. It is designed for home care, for travelling and easy storage. The Molift Smart is foldable in one piece without using any tools.

Due to its low weight, small footprint and unique construction, Molift Smart 150 is easy to manoeuvre even in narrow spaces. The hoisting range provides an excellent maximum hoisting height and and can also lift from the floor. The 4-point sling bar ensures a comfortable and spacious hoisting position for the user.

Molift Smart 150 has been awarded for Design Excellence by the Norwegian Design Council.





Foldable without tools Easy handling simplifies everyday life and transport when traveling.



4-point sling bar Ensures the user a comfortable and spacious hoisting position.



Easy to manoeuvre Low weight and small footprint facilitates use in narrow spaces.



Lift from floor Great hoisting interval that allows lifting from the floor.



Accessories



Molift Soft Travel Bag

The Molift Soft Travel Bag protects and enables easy transporting of the hoist, battery, charger and sling during travel.



Molift Travel Suitcase

This hard shell suitcase provides additional protection when transporting by air. It is TSA approved for air travel and ideal for the client with an active life.



Product name

Molift Smart 150	1020
Molift Soft Travel Bag	3049177
Molift Travel Suitcase	0990200

Material

Aluminium, steel

Weight

Total Weight Weight when separated

Dimensions

Length Leg height Dimensions folded Hoisting height Hoisting interval

1200 mm 110 mm L1160 x W475 x H360 mm 270-1665 mm 1410 mm

Part 1: 13 kg, Part 2: 13 kg

26 kg

Item no. Technical

Service Software Battery Battery charge time Hoisting speed Rated performance

Included NiMH 14.4 V - 2.6 Ah 50 hoists/75 kg/500 mm



approx. 3 hours 60 mm/sec

CE

Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85°C for 3 minutes.

For the latest news and continuously updated product information - visit: www.etac.com

Etac AB Kista Science Tower SE-164 51 Kista, Sweden Tel +46 371 - 58 73 30 info@etac.se www.etac.com

