

CHILD RESTRAINT SYSTEM

CARROT 3

BOOSTER SEAT



- **Adaptive Positioning**
High density foam inserts for optimal positioning and maximum comfort
- **Grows With Your Child**
Able to accommodate growth from 79-165 lbs / 55"-69"
- **"Free Angle" Recline up to 180°**
Conforms to the vehicle's seat
- **LATCH Equipped**
Lower anchors and tethers for children for easy and safe installation
- **Offers a Large Configuration**
To meet the needs of children and adults who require back height support and a wide seat

Shoulder Protector Wings



Gives lateral shoulder support and protection
(Width Adjustable)

Wide Seat



The Carrot 3 Booster Seat is for users who require a wide seat
(16.5")

4-Point Positioning H-Harness



To maintain midline positioning of the trunk for increased postural control

Booster Kit



Option to transform an existing Carrot 3 into a booster in 3 easy steps


Lightweight • Compact Folding • Convenient

Etac - United States
2830 California Street
Torrance, CA 90503
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com




Convaid
by Etac

Carrot 3 Booster Kit (Option to transform Carrot 3 Seat to Booster Seat)



1 Booster Seat Base


2 4-Point Positioning H-Harness



3 H-Harness Attaching Bracket

4 Backrest Attaching Metal Rods


Carrot 3 Booster Kit Required Extensions



1 Shoulder Protector Wings

2 2" Backrest Extension

Accessories



Booster Support Tray
(Crash Tested)



Seat Angle Adjustment Wedge
(Height 4.25")



3D Body Balance Pads



Backrest Extension
(4")

Upholstery Colors



Panther Black



Cherry Red




Royal Blue

Carrot 3 Specifications with Extensions

Carrot 3 Car Seat

Carrot 3 Booster Seat



*Height is measured from the top of the headrest to the seat.

Booster Seat Standard Specifications

Measurements:	Seat Width	Seat Depth	Backrest Height	Head Support Height	Recline Adjustment	Product Weight	Product Overall Width
INCHES	16.5	17	28.8 - 32.88	8.5	90° - 180°	21 lbs	20 - 22
	User Weight Range	User Height Range	Shoulder Protector Wings Height	Backrest Height	Backrest Extensions Height	Shoulder Wing Width Adjustment	
INCHES	79 lbs - 165 lbs	55 - 69	6	11.8	2 or 4	15.25 - 19.25	

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately. Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.