



UPRIGHT EZ RIDER

- **10° Fixed Tilt**
Provides upright posture to enhance functional positioning
- **Multiple Upholstery Styles & Colors**
Removable, washable Cordura or Textilene available in various colors
- **5 Years of Growth Capabilities**
Adjustable seat depth extends chair life for a growing child
- **Self-Tensio®**
Built in self-tensioning seating system increases sitting tolerance



Positioning options & accessories sold separately

Transit




Meets or exceeds WC19 safety standards for safe transportation

Planar




Solid seat & back provides in-depth seating & therapeutic positioning

Convertible



Removable Mag Rear Wheels allows independent mobility

Push Handle



Two Piece Push Handle (Option)
or
One Piece Height Adjustable Push Handle (Standard)

Lightweight • Compact Folding • Convenient

Convaidd Products, LLC
2830 California Street
Torrance, California 90503
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com



Convaidd
by Etac

Accessories

- Extended Headrest Cover with windows **A** (Canopy)
- Extended Headrest Cover without windows (Canopy)
- Headrest Cover (Canopy)
- Mosquito Net **B** (Canopy Required)
- Rain Cover (Canopy Required)
- Upper Extremity Support Surface Tray **C** (Clear Tray) (Armrests Required)
- Utility Bag **D**
- Wheelchair Travel Bag
- Rear Anti-Tip Tubes **E**
- Caster Locks (For 2" Wide Front Wheels)
- Mag Wheel Bag **F** (EZ CV Required)



Frame Options

- One Piece Height Adjustable Push Handle (Standard)
- Two Piece Push Handle
- Height Adjustable Flip-up Armrests
- 4" Extended Two Piece Push Handle

Transit Options

- 5-Point Harness (EZ12 Only)
- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt
- Upper Extremity Support Surface Tray (Foam Tray) (Armrests Required)

Upholstery Options



(Zippered Pockets in Cordura Seat & Back Upholstery)

EZ Rider 18



Crash tested model shown

Standard Features

- Two Piece Push Handle
- 8" Front, 11" Rear Solid Tires
- Cordura or Textile Upholstery (Cordura Shown)
- Toggle Wheel Locks

Head Supports



- Headrest Extension (Required for Transit) - Curved Cordura (Shown) - Rectangular Textile

- Occi Headwings

- Padded Headwings

Trunk Positioning



- H-Harness with Padded Covers (Required for Transit)

- Soft Adjustable Lateral Support w/ Scoliosis Strap (Single or Double Flap)

- Reducer Seat Insert (Shown)
- Full Torso Support Vest

Pelvic Positioning



- Support Cushion

- Align Cushion

- Position Cushion

- Solid Seat Insert (Required for all Cushions)

- Incontinence Cushion Liner (Solid Seat Insert & Cushion Required)

Upholstery Colors

Textile

Navy Blue

Teal Green

Violet Purple

Wild Raspberry

Cordura

Panther Black

Cherry Red

Sassy Purple

Princess Pink

Royal Blue


Forest Green

Chocolate Brown

Aqua Blue

Apple Green

Foot & Leg Positioning



- Angle Adjustable Footplates (Shown)
- Foot Positioners (Required for Transit)
- Footplate Securement Strap

- Lateral Thigh Support Adductor (Shown)
- Medial Thigh Support Abductor
- Calf Panel

Compact Fold



Closure Strap

Tire Options

EZ -EZ Rider EZ CV - EZ Rider Convertible Option

Front Tires

- 7.5" x 2" Front Pneumatic (Rolko) ----- EZ 12/14/16
- 8" x 2" Front Solid Ribbed ----- EZ 18
- 8" x 2" Front Pneumatic Ribbed ----- EZ 18

Rear Tires

- 12.5" x 2.5" Rear Solid Knobby (Rolko) ----- EZ 12/14/16
- 12.5" x 2.5" Rear Pneumatic Knobby (Rolko) ----- EZ 12/14/16
- 12.5" x 2.5" Rear Solid Knobby ----- EZ 18
- 12.5" x 2.5" Rear Pneumatic Knobby ----- EZ 18
- 20" Mag Rear Solid Tires ----- EZ CV 12
- 22" Mag Rear Solid Tires ----- EZ CV 14
- 24" Mag Rear Solid Tires ----- EZ CV 16/18