

FIXED-TILT CRUISER



- 30° Fixed Tilt**
 Increased tilt angle improves postural control through upper body positioning..
- Multiple Upholstery Styles & Colors**
 Removable, washable Cordura or Textilene available in various colors
- 5 Years of Growth Capabilities**
 Adjustable seat depth extends chair life for a growing child
- Self-Tensio®**
 Built in self-tensioning seating system increases sitting tolerance



Positioning options & accessories sold separately

Transit



Meets or exceeds WC19 safety standards for safe transportation

Planar



Solid seat & back provides in-depth seating & therapeutic positioning

Scout



All-terrain option built for maximum comfort and durability

Push Handle



Two Piece Push Handle (Option) or One Piece Height Adjustable Push Handle (Standard)

Lightweight • Compact Folding • Convenient

Convoid Products, LLC
 2830 California Street
 Torrance, California 90503
 +1 844-US-MOBILITY (844-876-6245)
 www.etac.us.com



Convoid
 by Etac

Accessoires

- Extended Headrest Cover with windows (Canopy) **A**
- Extended Headrest Cover without windows (Canopy)
- Headrest Cover (Canopy)
- Mosquito Net (Canopy Required) **B**
- Rain Cover (Canopy Required)
- Attendant Hand Brake **C**
- Upper Extremity Support Surface Tray (Clear Tray) **D**
- Utility Bag
- Saddle Bags (Pair)
- Wheelchair Travel Bag
- Under Seat Storage Basket **E**
- Rear Anti-Tip Tubes **F**
- Curb Tipper **G**
- Caster Locks (For 2" Wide Front Wheels) **H**



Head Supports

- Headrest Extension (Required for Transit)
 - Curved Cordura (Shown)
 - Rectangular Textilene
- Occi Headwings
- Padded Headwings

Trunk Positioning

- H-Harness with Padded Covers (Required for Transit)
- Soft Adjustable Lateral Support w/ Scoliosis Strap (Single or Double Flap)
- Reducer Seat Insert (Shown)
- Full Torso Support Vest

Pelvic Positioning

- Support Cushion
- Align Cushion
- Position Cushion
- Solid Seat Insert (Required for all Cushions)
- Incontinence Cushion Liner (Solid Seat Insert & Cushion Required)

Frame Options

- One Piece Height Adjustable Push Handle (Standard)
- Two Piece Push Handle
- Anatomic Back (For Additional Lumbar Support)

Transit Options

- 5-Point Harness (CX10/12 Only)
- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt

Upholstery Options



(Zippered Pockets in Cordura Seat & Back Upholstery)

Cruiser 18



Crash tested model shown

Standard Features

- Two Piece Push Handle
- 8" Front, 11" Rear Solid Tires
- Cordura or Textilene Upholstery (Cordura Shown)
- Toggle Wheel Locks

Upholstery Colors

- Navy Blue
- Teal Green
- Violet Purple
- Wild Raspberry

Cordura

- Panther Black
- Cherry Red
- Sassy Purple
- Princess Pink
- Royal Blue
- Forest Green
- Chocolate Brown
- Aqua Blue
- Apple Green

Foot & Leg Positioning

- Medial Thigh Support Abductor
- Lateral Thigh Support Adductor (Shown)
- Calf Panel
- Foot Positioners (Required for Transit)
- Footplate Securement Strap
- Angle Adjustable Footplates

Compact Fold



Closure Strap

Tire Options

cx - Cruiser **cx sc** - Cruiser Scout Option

Front Tire

- 7.5" x 2" Front Pneumatic (Rolko) ----- **cx** 10/12/14/16
- 8" x 2" Front Solid Ribbed----- **cx sc** 12/14/16/18 **cx** 18
- 8" x 2" Front Pneumatic Ribbed----- **cx sc** 12/14/16/18 **cx** 18

Rear Tire

- 12.5" x 2.5" Rear Solid Knobby (Rolko)----- **cx** 10/12/14/16
- 12.5" x 2.5" Rear Pneumatic Knobby (Rolko)----- **cx** 10/12/14/16
- 12.5" x 2.5" Rear Solid Knobby----- **cx** 18
- 12.5" x 2.5" Rear Pneumatic Knobby----- **cx** 18