

*Fixed - Tilt*

# METRO



- **Contoured Upholstery**

Provides a closer fit for greater support and comfort

- **Adjustable Trunk Support**

Promotes superior lumbar positioning

- **Self-Tensio®**

Built in self-tensioning seating system increases sitting tolerance



## Transit



Adheres to WC19 safety standards for safe transportation  
(optional)

## Mesh Upholstery Center



Increases ventilation to keep the occupant cooler

## "Quick Release" Rear Wheels



Easy to remove for storage or transport

## Two-Piece Push Handle



Quick and easy to fold, store, and transport

**Lightweight • Compact Folding • Convenient**

Convaid Products, LLC  
2830 California Street  
Torrance, California 90503  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)



 **Convaid**  
by Etac

Accessories



3-Pt Positioning Belt with  
Depth Adjustable Crotch Strap  
(Transit Required)



Travel Bag



Under Seat Storage  
Basket

Tire Options

- 12.5" Rear Solid Knobby Tire
- 7.5" x 2" Front 12.5" Rear Pneumatic Knobby Tire
- 8" x 2" Front Ribbed, 12.5" Rear Solid or Pneumatic Knobby Tire

Transit Options

- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt


Measurement	Chair Model	Seat Width*	Overall Width	Overall Height	Overall Length	Seat Depth*	Seat Back Height	Seat to Footrest*	Seat to floor	Seat Angle	Headrest Extension	Weight Capacity Regular/Transit	Weight of Chair Δ
in.	MT 12	12	23.5	36.8	36.5	11	21	8 - 15	24	30°	6	75 / 66 lbs	23 lbs
cm	MT 12	30.5	60	93.3	92.7	28	53.3	20.3 - 38	61		15	34 / 30 kg	10.4 kg
in.	MT 14	14	25.5	37.5	34.5	13	23	8 - 14.5	24	30°	8	100 / 100 lbs	24 lbs
cm	MT 14	35.5	64.8	95.3	87.6	33	58.4	20.3 - 36.8	61		20	45.5 / 45.5kg	11 kg
in.	MT 16	16	27.5	38	35	14	25	13 - 18	23	30°	8	170 / 170 lbs	27 lbs
cm	MT 16	40.5	70	96.5	89	35.6	63.5	33 - 45.7	58.4		20	77 / 77 kg	12.3 kg
in.	MT 18	18	28.5	38	42	16	25	3 - 24	23.5	30°	8	250 / 170 lbs	29 lbs
cm	MT 18	45.5	72.4	96.5	106.7	40.6	63.5	7.3 - 61	59.7		20	114 / 77 kg	13.2 kg

in. - Inch, cm - Centimeter \* Plus or Minus .5 inches (1.27 cm) due to fabric. Δ- Weight of chair is without front rigging

Upholstery Colors



Cherry Red



Royal Blue



Metro Transit