

Convaid Flyer[®]

A Truly Complex Rehab Chair. The lightest and most portable pediatric chair on the market is now offered with a complete line-up of complex rehab wheelchair positioning accessories and options.

Hardware Kit for Vest & 4-Point Belts

Quick adjustability and removal of vest and belts for accessibility and transfers

Removable Adjustable Armrest by R82

Quick removal for convenience and transfer. Angle, height and width adjustable

Swing-Away Trunk Laterals by R82

Easy swing-away mechanism for transfers and accessibility

Medial Thigh Support

Quick removal for transfers and accessibility

Elevating Leg Rest Cover

NEW

Quick Release Axles


Remove with the wheel to make the Flyer even more compact and lighter for easy transport

NEW


65° Leg Supports

New 65° removable leg supports accommodate longer leg lengths and avoid caster interference

Mobility Base Color Options



Royal Blue



Cosmic Silver



Cherry Red



Plum Purple

NEW