

# Upright VIVO



- **10° Fixed Tilt**  
Promotes an upright position for individuals with good trunk strength and upper body stability
- **95° - 110° of Adjustable Recline (Size 16 Only)**  
Opens the seat-to-back angle for better digestion, improved respiration, and pressure relief
- **5 Years of Growth Capabilities**  
Extend the life of the chair for a growing child
- **Self-Tensio®**  
Built in self-tensioning seating system increases sitting tolerance



## Frame Options



Height  
Adjustable Armrest

## Adjustable Back Upholstery



Provides Contoured  
Positioning

## Zippered Pocket in Upholstery



Allow additional foam  
padding to be added for  
firmer support

**Lightweight • Compact Folding • Convenient**

Convaid Products, LLC  
2830 California Street  
Torrance, California 90503  
+1 844-US-MOBILITY (844-876-6245)  
www.etac.us.com



**Convaid**  
by Etac

# Accessories

- Headrest Cover (Canopy)
- Grab Bar
- Upper Extremity Support Surface Tray
- Travel Bag

## Pelvic Positioning



3-Pt Positioning Belt  
with Depth Adjustable  
Crotch Strap

## Trunk Positioning



H-Harness with  
Padded Covers



Full Torso Support Vest

## Frame Options

- Height Adjustable Armrests (Pair)

## Upholstery Colors

Cordura



Panther  
Black

Cherry  
Red



Royal  
Blue

Princess  
Pink




Sassy  
Purple

Forest  
Green



Apple  
Green

Aqua  
Blue



Chocolate  
Brown

## Foot & Leg Positioning



Medial Thigh Support  
(Abductor)



Lateral Thigh Support  
(Adductor)

- Footplate  
Securement Strap

- Foot Positioners



Calf Panel

## Head Supports

- Occi Headwings
- Headrest Extension

## Standard Specifications

Measurement	Chair Model	Seat Width*	Seat Depth*	Long Seat Depth	Seat Back Height	Seat to Footplate	Seat to Floor	Long Back Height	Seat to Back Angle	Seat Angle	Overall Height	Push Handle Adjustments	Overall Length	Overall Width	Folded Length	Folded Width	Folded Height	Weight Capacity	Weight of Chair $\Delta$
in	VV 12	12	8 - 10	11 - 13	21.5	7.9 - 11	17.3	N/A	95° - 110°	10°	36.6	N/A	31.5	19.7	38.2	14.6	17	44 lbs	14.3 lbs
cm	VV 12	30.5	20.3 - 25.4	27.9 - 33	54.6	20 - 28	43.9	N/A	---	---	92.96	N/A	80	50	97	37	43.2	20 kg	6.5 kg
in	VV 14	14	10 - 12	12 - 15	23	9.8 - 14	18.5	26	95° - 110°	10°	41.2	N/A	39	21.5	46.6	11.4	17.4	88.2 lbs	17.2 lbs
cm	VV 14	35.6	25.4 - 30.5	30.5 - 38.1	58.4	25 - 35.6	47	66.04	---	---	104.6	N/A	99.1	54.6	118.4	29	44.2	40 kg	7.8 kg
in	VV 16 <sup>o</sup>	16	13 - 15	16 - 18	27	7 - 21	22	N/A	95° - 120°	12°	N/A	32.5 - 47.3	44.5	24.75	50	13	15	121 lbs	27 lbs
cm	VV 16 <sup>o</sup>	40.6	33 - 38	40.6 - 45.8	68.6	17.8 - 53.3	55.9	N/A	---	---	N/A	82.6 - 120	113	62.9	127	33	38.1	55 kg	12.2 kg

<sup>o</sup> One Piece Push Handle comes standard on VV16 only

\* Plus or Minus .5" due to fabric

$\Delta$  Weight of chair is without front rigging