


Slings

Recycling Instructions (EN)

molift[®]
by Etac

RI30001 Rev C 2022-11-30



All Nova and Molift sling models are made of different types of textiles and can be disposed of as general waste.



Details like buckles and ladder locks are made of recyclable plastics and can be removed from the sling and recycled separately. Scissor may be required on some models.



Some details are made of metals. These can be removed from the sling and recycled separately. Scissor may be required.

Some countries and regions may have collection of textiles for recycling. Always follow local rules and regulations.



Etac A/S
Parallevej 3
DK-8751 Gedved
www.etac.com