


Etac Slings

Recycling Instructions English

 **etac**[®]

RI30001 Rev B 2015-06-01



All Nova and Molift sling models are made of different types of textiles and can be disposed of as general waste.



Details like buckles and ladderlocks are made of recyclable plastics (e.g. nylon) and can be removed from the sling and recycled separately. Scissor may be required on some models.



Some details are made of metals. These can be removed from the sling and recycled separately. Scissor may be required.

Some countries and regions may have collection of textiles for recycling. Always follow local rules and regulations.

Etac Sverige AB
Box 203, 334 24
Anderstorp, Sweden
Tel 46 371-58 73 30
Fax 46 371-58 73 90
info@etac.se www.etac.se

Etac AS,
Etac Supply Gjøvik
Hadelandsveien 2,
2816 Gjøvik, Norway
Tel +47 4000 1004
molift@etac.com www.molift.com