


Rufus

Duschstol hopfällbar med urtag

Art. nr: 8120 8012



Pos	Benämning	Art nr
1	Teleskopbenhylsa	84001150
2	Skruvsats Rufus 4 st	84001168
3	Nivåfäste Rufus 2 st	84001177
4	Doppsko 20.5 mm 4 st	84001193
5	Bricka Rufus 10 st	84001197
6	Rörändpropp 10 st	84001202
7	Clips 12 st	84001210
8	Sits med urtag Rufus 68 cm	84001219
9	Sits standard Rufus med urtag 130 st	87001017
10	Ryggbricka Rufus 200 st	87001020

78347A

05-08-23

Etac AB, Box 203, SE-334 24 Anderstorp

Tel. +46/371-58 73 00, Fax. +46/371-58 73 90

www.etac.com