


# Rufus

Showerchair Rufus foldb.recess

Art. No: 8120 8012



Pos	Description	Art no
1	Telescopic legcasing 12pcs.....	84001150
2	Screw kit Rufus 4 pcs .....	84001168
3	Seat level socket Rufus 2 pcs .....	84001177
4	Ferrule 20,5mm/4pcs .....	84001193
5	Washer Rufus 10 pcs .....	84001197
6	End cover 10 pcs.....	84001202
7	Clip 12 pcs .....	84001210
8	Seat Rufus 68 cm recess .....	84001219
9	Seat standard with recess .....	87001017
10	Backrest washer .....	87001020

78357A

06-02-06

Etac, Box 203, SE-334 24 Anderstorp

Tel. +46/371-58 73 00, Fax. +46/371-58 73 90

www.etac.com