




Spare part list

Molift Transfer Pro

molift®
by Etac

SPL70199 - 2024-09-30 - Rev. 3.0



www.etac.com



Fig.1




Fig.2




Fig.3




Fig.4




Fig.5




Fig.	Pos.	Title	Item no.
1	1	Front caster	1630301
	2	Rear caster	1630302
	3	Brake	1630304
	4	Foot plate	1630289
	5	Leg adjustment covers – RH/LH	1630309
2	6	Seat pad RH	1630306
	6	Seat pad LH	1630307
3	7	Leg support	1630308
	8	Short user guide (no illustration)	1630305
	9	Mounting kit for Transfer Pro (no illustration)	1630010
4	10	Handle	1630310
5	11	Leg adjustment mechanism	1630311