




Spare part list – EN

Molift Partner 255

molift[®]
by Etac

SPL03401 - 2026-09-11 - Rev. 5.0



www.etac.com



Fig.1




Fig.3




Fig.2




Fig.	Pos.	Title	Item no.
1	1	Front wheel kit ø75 Mover 205/300/Partner 255	0950210
	2	Rear wheel ø100 w. brake Mover 205/300/Partner 255	0950211
	3	Locking handle kit chassis Mover 205/300/Partner 255	0950249
	4	Chassis contact kit Mover 205/300/Partner 255	0950246
	5	Special tool Molift Smart 150/Mover/Partner 255 (no illustration)	0920501
2	6	Leg adjust mechanism repair kit with motor	1420540
	7	Leg adjust mechanism repair kit without motor	1420530
	8	Motor, leg adjust mechanism	1320155
3	9	Spare handle kit Mover 205/300/Partner 255	0950239
	10	Control box mounting kit	0950247
	11	Spare control box, BLE	0950252
	11	Spare ctrl box w/o chrg. BLE (only available for Japan)	0950252-J
	12	Channel1 cable f/BLE contr.box	1860116
	13	Channel2 cable f/BLE contr.box	1860117
	14	Cable RJ50 f/handctrl. 2&4 btn (no illustration)	0950248

Fig.4




Fig.5




Fig.7




Fig.8




Fig.9



Fig.10



Fig.11



Fig.12

U In : 120-240 V AC, 50-60 Hz
I In : Max. 350 mA, P In: 35-55 VA

Fig.	Pos.	Title	Item no.
4	15	Cable for Molift Service Tool (USB/RJ50)	2510250
	16	Prog.box kit f/BLE control box (only for service technician)	0950254
5	17	Spare cable f/prog. box, RJ50 (only for service technician)	0950152
	18	USB A to USB B cable 1M (only for service technician)	0950153
6	19	Spiral cable USA plug (no illustration)	2920092
	20	Spiral cable EU plug (no illustration)	2920093
	21	Spiral cable GB plug (no illustration)	2920094
	22	Spiral cable AU plug (no illustration)	2920095
7		Product label	1860111
8		Label Safety Check	1100303
9		Label Safety Check	1110306
10		Label Molift by Etac	0950128
11		Label charging instructions (not available for Japan)	0950123
12		Label elec.data f/BLE ctrl.box (not available for Japan)	0950124

**Etac Ltd.**

Unit 60, Hartlebury Trading Estate,
Hartlebury, Kidderminster,
Worcestershire, DY10 4JB
+44 121 561 2222

**Promefa SA**

Avenue de Morges 33
1027 Lonay, CH
+41 44 872 97 79

**Etac A/S**

Parallevej 3
DK-8751 Gedved



+45 796 858 33



patienthandling@etac.com



www.etac.com

molift[®]
by Etac