




Spare part list – EN

Molift Mover 180

molift[®]
by Etac

SPL08201 - 2025-03-12 - Rev. 1.0



www.etac.com



Fig.1




Fig.2




Fig.3




Fig.	Pos.	Title	Item no.
1	1	Front wheel ø75 kit Smart 150/Mover 180	0950201
	2	Front wheel ø100 kit Mover 180	0950207
	3	Front wheel bracket kit Mover 180	0950206
	4	Rear wheel ø75 kit w. brake	0950205
	5	Rear wheel ø100 kit w. brake	0950208
	6	Chassis contact kit Mover 180	0950245
	7	Special tool Molift Smart/Mover (no illustration)	0920501
2	8	Leg adjust mechanism repair kit with motor	1420540
	9	Leg adjust mechanism repair kit without motor	1420530
	10	Motor for leg adjust	1320155
3	11	Actuator kit Mover 180	0950229
	12	Spare handle kit Mover 180	0950238
	13	Control box COL50 mounting kit	0950247
	14	CH1 COL50 cable	1860116
	15	CH2 COL50 cable	1860117
	16	COL50 control box	0950252
	17	Spare cable hand control HB32/HB31 (no illustration)	0950248

Fig.6



Fig.8




Fig.9



Fig.10



Fig.11




Fig.4




Fig.	Pos.	Title	Item no.
4	18	Programming PB31 box for COL50 (only for service technician)	0950151
	19	Spare cable for PB31 - COL50 (only for service technician)	0950152
	20	USB A to USB B cable 1M (only for service technician)	0950153
5	21	Spiral cable USA plug (no illustration)	2920092
	22	Spiral cable EU plug (no illustration)	2920093
	23	Spiral cable GB plug (no illustration)	2920094
	24	Spiral cable AU plug (no illustration)	2920095
6		Product label	1860111
7		Label Safety Check	1100303
8		Label Safety Check	1110306
9		Label Molift by Etac	0950128
10		Label charging instructions	0950123
11		Label COL50 electrical data	0950124



Etac Ltd.

Unit 60, Hartlebury Trading Estate,
Hartlebury, Kidderminster,
Worcestershire, DY10 4JB
+44 121 561 2222



Sodimed SA

Ch. Praz-Devant 12
CH-1032 Romanel-sur-Lausanne
+41 58 911 06 06



Etac A/S

Parallevej 3
DK-8751 Gedved



+45 796 858 33



patienthandling@etac.com



www.etac.com

molift[®]
by Etac