

Wombat Solo

Exploded views / Part lists

Rev. date (2025.11)

On the following pages you will find the parts used to assemble this product.
Please note that some spare parts cannot be purchased individually but are part of a complete kit. Please also note that for some spare parts extra delivery time may occur.




Contents

Wombat Solo 1-2	4
Item numbers for Wombat Solo 1-2.....	5
Wombat Solo 1-2 lift on frame	6
Item numbers for Wombat Solo 1-2 lift on frame	7
Wombat Solo 1-2 lift/tilt on frame.....	8
Item numbers for Wombat Solo 1-2 lift/tilt on frame	9
Wombat Solo 1-2 lift frame tilt on brace	10
Item numbers for Wombat Solo 1-2 lift frame tilt on brace	11
Wombat Solo 1-2 lift and tilt on brace.....	12
Item numbers for Wombat Solo 1-2 lift and tilt on brace.....	13
Wombat Solo 1-2, no tilt	14
Item numbers for Wombat Solo 1-2, no tilt	15
Wombat Solo back size 1-2	16
Item numbers for Wombat Solo back size 1-2	16
Wombat Solo seat size 1-2.....	17
Item numbers for Wombat Solo seat size 1-2.....	18
Clutch house complete Wombat Solo size 1+2	19
Item numbers for Clutch house complete Wombat Solo size 1+2	20
Wombat solo size 3 lift/tilt on frame	21
Item numbers for Wombat solo size 3 lift/tilt on frame	22
Wombat Solo size 3 lift on frame.....	23
Item numbers for Wombat Solo size 3 lift on frame.....	24
Wombat Solo size 3 lift and tilt on brace	25
Item numbers for Wombat Solo size 3 lift and tilt on brace.....	26
Wombat Solo size 3	27
Item numbers for Wombat Solo size 3	28
Wombat Solo size 3, no tilt.....	29
Item numbers for Wombat Solo size 3, no tilt.....	30
Wombat Living back size 3.....	31
Item numbers for Wombat Living back size 3.....	32
Wombat seat 3 for solo	33
Item numbers for Wombat seat 3 for solo.....	33
Clutch house complete Wombat size 3	34
Item numbers for Clutch house complete Wombat size 3	35

Frame assembly	36
Item numbers for Frame assembly	36
Footrest frame, size 1-2-3	37
Item numbers for Footrest frame, size 1-2-3	38
Footrest seat, size 1-2-3	39
Item numbers for Footrest seat, size 1-2-3	40
Fixed side supports, with covers (95121, 95121-S & 95121-23)	41
Item numbers for Fixed side supports, with covers (95121)	42
Item numbers for Fixed side supports, with covers (95121-S)	42
Item numbers for Fixed side supports, with covers (95121-23)	42


Wombat Solo 1-2



Item numbers for Wombat Solo 1-2

Position	No.	Description	Quantity	Unit of Measure	Ean code
1	134307-100CBSF	Solo Frame std m.100mm wheel+CB	1	PCS	5707292875129
	134307-100CBS	Solo Frame std m.100mm wheel+CB	1	PCS	5707292813060
	134307-075SF	Solo Frame m. 75 mm wheel	1	PCS	5707292875068
	134307-075S	Solo Frame std	1	PCS	5707292811134
	134307-100SF	Solo Frame std m. 100 mm wheel	1	PCS	5707292875099
2*	134302-TLSW	Cloth house complete size 1+2	1	PCS	5707292837325
	134302-TLS	Cloth house complete size 1+2	1	PCS	5707292811370
3*	139245-2	Wombat Solo seat size 2	1	PCS	5707292868909
	139245-1	Wombat Solo seat size 1	1	PCS	5707292868893
4*	121092-6	Wombat size 2 eith high back	1	PCS	5707292784834
	121092-5	Wombat size1 with high back	1	PCS	5707292784827
	121092-21	Wombat back size 2	1	PCS	5707292622822
	121092-11	Wombat back size 1	1	PCS	5707292807533
	121092-1	Wombat back size 1	1	PCS	5707292622808
5*	139214-2	Footrest Simple frame size 2	1	PCS	5707292869623
	135020-2	Footrest size 2 Wombat Solo	1	PCS	5707292813671
	139214-1	Footrest Simple frame size 1	1	PCS	5707292869616
6*	135020-2	Footrest size 2 Wombat Solo	1	PCS	5707292813671
	0609826	Stainless insex M6x22	4	PCS	5707292394477
	135020-1	Footrest size 1 Wombat Solo	1	PCS	5707292813619
7	0609826	Stainless insex M6x22	4	PCS	5707292394477
*Choose size, color, type					


Wombat Solo 1-2 lift on frame



Item numbers for Wombat Solo 1-2 lift on frame

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	121092-11	Wombat back size 1	PCS	1	5707292807533
	121092-6	Wombat size 2 eith high back	PCS	1	5707292784834
	121092-5	Wombat size1 with high back	PCS	1	5707292784827
	121092-21	Wombat back size 2	PCS	1	5707292622822
2*	139245-1	Wombat Solo seat size 1	PCS	1	5707292868893
	139245-2	Wombat Solo seat size 2	PCS	1	5707292868909
3*	134302-LFTS	Cloth house complete size1+2	PCS	1	5707292811479
	134302-LFTSW	Cloth house complete str 1+2	PCS	1	5707292837332
4	0387000	End nipple for outer casing	PCS	1	5707292136091
5	0388006	Screwnipple 5x7,5mm	PCS	1	5707292136114
6	0380060	Flexible v-brake pipe.	PCS	1	5707292425195
7	0371853	Strained wires,1,5 mm	PCS	1	5707292135889
8	1760005	Compensating nut 5mm	PCS	1	5707292141491
9	0380000	Cable tighter M5x23	PCS	1	5707292135933
10	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
11	127592	Cable fix	PCS	1	5707292810540
12	0609737	Hexalobular (6 Lobe) flat head	PCS	2	5707292828743
13*	134307-075LFS	Solo stel m 75 mm wheel	PCS	1	5707292813008
	134307-100CBLFF	Solo stel m 100 mm wheel+CB	PCS	1	5707292875136
	134307-075LFSF	Solo frame 75 mm wheel	PCS	1	5707292875075
14	136393	Kabelbøjle	PCS	1	5707292871169
15	139214-1	Footrest Simple frame size 1	PCS	1	5707292869616
16*	135020-2	Footrest size 2 Wombat Solo	PCS	1	5707292813671
	135020-1	Footrest size 1 Wombat Solo	PCS	1	5707292813619
17	0609826	Stainless insex M6x22	PCS	4	5707292394477
*Choose size, color, type					


Wombat Solo 1-2 lift/tilt on frame



Item numbers for Wombat Solo 1-2 lift/tilt on frame

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	121092-5	Wombat size1 with high back	PCS	1	5707292784827
	121092-21	Wombat back size 2	PCS	1	5707292622822
	121092-11	Wombat back size 1	PCS	1	5707292807533
2*	139245-1	Wombat Solo seat size 1	PCS	1	5707292868893
	139245-2	Wombat Solo seat size 2	PCS	1	5707292868909
3*	134302-LFBTFBW	Cloth house compl str1+2, SWIV	PCS	1	5707292837349
	134302-LFBTFB	lift frame tilt frame	PCS	1	5707292811486
4	0387000	End nipple for outer casing	PCS	2	5707292136091
5	0388006	Screwnipple 5x7,5mm	PCS	2	5707292136114
6	0380060	Flexible v-brake pipe.	PCS	2	5707292425195
7	0371853	Strained wire, 1,5 mm	PCS	1	5707292135889
8	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
9	1760005	Compensating nut 5mm	PCS	2	5707292141491
10	0380000	Cable tighter M5x23	PCS	2	5707292135933
11	0380060	Flexible v-brake pipe.	PCS	1	5707292425195
12	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
13	0371853	Strained wire, 1,5 mm	PCS	1	5707292135889
14	127592	Cable fix	PCS	2	5707292810540
15	0609737	Hexalobular (6 Lobe) flat head	PCS	4	5707292828743
16	136393	Kabelbøjle	PCS	1	5707292871169
17*	134307-075LTFSF	Solo stel m 75 mm wheel	PCS	1	5707292875082
	134307-100CBLTF	Solo stel m 100 mm wheel+CB	PCS	1	5707292813329
18*	139214-2	Footrest Simple frame size 2	PCS	1	5707292869623
	139214-1	Footrest Simple frame size 1	PCS	1	5707292869616
19	135020-1	Footrest size 1 Wombat Solo	PCS	1	5707292813619
20	0609826	Stainless insex M6x22	PCS	4	5707292394477
*Choose size, color, type					


Wombat Solo 1-2 lift frame tilt on brace



Item numbers for Wombat Solo 1-2 lift frame tilt on brace

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	121058	Push handle complete str 1	PCS	1	5707292623331
2	3503534	Gas pump lever	PCS	1	5707292845634
3	0380060	Flexible v-brake pipe.	PCS	2	5707292425195
4	0371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
5	5002998	Cable for new rabbit, 39 cm	PCS	1	5707292436016
6	0235075	Cable tape, black	PCS	1	5707292135308
7	0388006	Screwnipple 5x7,5mm	PCS	2	5707292136114
8	0387000	End nipple for outer casing	PCS	2	5707292136091
9	0380060	Flexible v-brake pipe.	PCS	1	5707292425195
10	121092-21	Wombat back size 2	PCS	1	5707292622822
11	139245-2	Wombat Solo seat size 2	PCS	1	5707292868909
12	134302-LFTSW	Cloth house complete str 1+2	PCS	1	5707292837332
13	0371853	Strained wire, 1,5 mm	PCS	1	5707292135889
14	1760005	Compensating nut 5mm	PCS	1	5707292141491
15	0380000	Cable tighter M5x23	PCS	1	5707292135933
16	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
17	127592	Cable fix	PCS	1	5707292810540
18	0609737	Hexalobular (6 Lobe) flat head	PCS	2	5707292828743
19	134307-100LFSF	Solo stel m 100 mm wheel	PCS	1	5707292875105
20	136393	Kabelbøjle	PCS	1	5707292871169
22	135020-2	Footrest size 2 Wombat Solo	PCS	1	5707292813671
23	0609826	Stainless inox M6x22	PCS	4	5707292394477


Wombat Solo 1-2 lift and tilt on brace



Item numbers for Wombat Solo 1-2 lift and tilt on brace

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	3503534	Gas pump lever	PCS	2	5707292845634
2	0371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
3	0380060	Flexible v-brake pipe.	PCS	2	5707292425195
4	5001600	Cable for back new Rabbit	PCS	1	5707292154491
5	0235075	Cable tape, black	PCS	4	5707292135308
6	0371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
7	5002998	Cable for new rabbit, 39 cm	PCS	1	5707292436016
8	0387000	End nipple for outer casing	PCS	2	5707292136091
9	0388006	Screwnipple 5x7,5mm	PCS	2	5707292136114
10	0380060	Flexible v-brake pipe.	PCS	2	5707292425195
11	121058	Push handle complete str 1	PCS	1	5707292623331
12	121092-21	Wombat back size 2	PCS	1	5707292622822
13	139245-2	Wombat Solo seat size 2	PCS	1	5707292868909
14	134302-LFBTFBW	Cloth house compl str1+2, SWIV	PCS	1	5707292837349
15*	134307-075SF	Solo Frame m. 75 mm wheel	PCS	1	5707292875068
	134307-100CBS	Solo Frame std m.100mm wheel+CB	PCS	1	5707292813060
16	139214-2	Footrest Simple frame size 2	PCS	1	5707292869623
17	135020-2	Footrest size 2 Wombat Solo	PCS	1	5707292813671
18	0609826	Stainless inox M6x22	PCS	4	5707292394477
*Choose size, color, type					

Wombat Solo 1-2, no tilt




C:\EPDM\Projects\1060\140927

Item numbers for Wombat Solo 1-2, no tilt

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-100CBSF	Solo Frame std m.100mm wheel+CB	PCS	1	5707292875129
	134307-075SF	Solo Frame m. 75 mm wheel	PCS	1	5707292875068
	134307-075S	Solo Frame std	PCS	1	5707292811134
2*	134302-SW	Cloth house complete str 1+2	PCS	1	5707292869197
	134302-S	Cloth house complete str 1+2	PCS	1	5707292869173
3	139245-2	Wombat Solo seat size 2	PCS	1	5707292868909
4	121092-21	Wombat back size 2	PCS	1	5707292622822
5	139214-2	Footrest Simple frame size 2	PCS	1	5707292869623
6	135020-2	Footrest size 2 Wombat Solo	PCS	1	5707292813671
7	0609826	Stainless insex M6x22	PCS	4	5707292394477

*Choose size, color, type

Wombat Solo back size 1-2




Item numbers for Wombat Solo back size 1-2

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	4913040	30mm adhesive hook, black	MTR	0,53	5707292786517
2*	82943	Back size 1, Wombat	PCS	1	5707292125989
	127158	Back 1 Wombat facelift	PCS	1	5707292789792
3	2501500	Ferrule d12x0,5-1,5 int.black	PCS	2	5707292143952
4	0445025	25mm press buckle, black	PCS	4	5707292465160
5	0658012	Button head M5x12 FZB	PCS	4	5707292139382
6	3503001	R82 handle M8x15mm screw	PCS	1	5707292413925
7	3502200	Tightening handle 42mm M6 bush	PCS	1	5707292146861
8	1109510	Starlock washer Ø8	PCS	2	5707292140456
9	1125618	Screenwasher fz6x18x3mm	PCS	1	5707292140494
10	82915	U-fitting for back	PCS	1	5707292125729
11	80920	Stick '8x54	PCS	1	5707292111364
12	0641200	Button head screw 6x12	PCS	4	5707292139252
13	82947-60	Back tilting device, black	PCS	1	5707292619822
14	82951	Inner back column	PCS	1	5707292530349
20	1905545	Bolt M6x60, FZB	PCS	1	5707292142306
21	1105080	Stainless washer 8mm	PCS	1	5707292140364
22	127164	Wombat size 2 back	PCS	1	5707292789808
MODUL 15	5082917	Back column for Wombat	PCS	1	5707292532268
MODUL-15	5082917	Back column for Wombat	PCS	1	5707292532268
MODUL-16	127819	Outer back column Long assm.	PCS	1	5707292786364

*Choose size, color, type


Wombat Solo seat size 1-2



Item numbers for Wombat Solo seat size 1-2

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	0641600	Button head screw 6x16	PCS	8	5707292139269
2	120701	Lock washer	PCS	1	5707292623416
7	1905710	Bolt M8x70 FZB wo/nut	PCS	1	5707292142368
8	0603665	Insex screw 8x65	PCS	1	5707292138286
9	0641620	Button head screw 6x20	PCS	2	5707292139276
10	0442070	Fixlock clamp 40 mm w/hole	PCS	2	5707292136398
11	126884	Seat size 1, Wombat	PCS	1	5707292793638
12	0506509	Rivet 4x10, stainless	PCS	4	5707292472441
13	4935300	100mm black velcro hook	MTR	0,15	5707292152107
14	0445025	25mm press buckle, black	PCS	2	5707292465160
15	0640920	Button head M5x20 FZB	PCS	2	5707292381033
16	82929	Distance piece	PCS	4	5707292125842
17	1103200	Fazet washer 8mm	PCS	1	5707292140296
18	1607000	Lock nut 8mm, large pack	PCS	1	5707292141378
19	1103500	Fazet washer 12mm	PCS	1	5707292140326
20	119445	Baseplate welded new	PCS	1	5707292621405
21	3502300	Tightening handle 65mm M8 bush	PCS	1	5707292146885
22	127065	Seat size 2, Wombat	PCS	1	5707292793669
23	4935300	100mm black velcro hook	MTR	0,21	5707292152107
24	0506509	Rivet 4x10, stainless	PCS	6	5707292472441
25	0603530	Insex low head M8x16 stainless	PCS	2	5707292412492
26	2502510	Plug 20x2,0-2,5 inside black	PCS	4	5707292557759


Clutch house complete Wombat Solo size 1+2



Item numbers for Clutch house complete Wombat Solo size 1+2

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	1606300	Lock nut 5 mm	PCS	3	5707292141354
2	136945	Insert cap small	PCS	1	5707292834287
3	128390	Cover	PCS	1	5707292810557
4	121245	Locking pin wombat solo	PCS	1	5707292810786
5	136440	Column assembly	PCS	1	5707292829825
6	0615022	M5x20 stainless hexagon screw	PCS	4	5707292837424
7	0609726	Insex 5x30 uh. Fzb	PCS	1	5707292433381
7	0615025	Insex screw 5x25	PCS	1	5707292138873
8	150677	Attaching plate	PCS	1	5707292918239
9	127725	Bearing disc	PCS	1	5707292794222
10	135987	Cover	PCS	1	5707292829870
11	2200545	Bloc-o-lift 200 N	PCS	1	5707292874528
12	1761000	Counter nut 10x1,0 fine thread	PCS	1	5707292141538
13	1103400	Fazet washer 10 mm	PCS	1	5707292140319
14	127836	Sliding dish	PCS	1	5707292794239
15	0640935	Button head M5x35 stainless	PCS	1	5707292485748
16	134717	Release handle for lift	PCS	1	5707292810779
17	128675-1	Handle w. seat lift pictogram	PCS	1	5707292830128
18	139080-1	clutch house sz 1-2, no tilt	PCS	1	5707292920911
18	139080	clutch house sz 1-2, no tilt	PCS	1	5707292868916
20	0380000	Cable tighter M5x23	PCS	1	5707292135933
21	127996	Bracket for tilt	PCS	1	5707290000000
22	658012	Button head M5x12 FZB	PCS	2	5710000000000
23	151104	Cloth house Wombat Solo str1+2	PCS	1	5707292917751
24	156511	Fork bracket	PCS	1	5707292968470
25	128261	Spring for tilt	PCS	1	5707292810809
26	1130500	Stainless washer M5, din9021	PCS	1	5707292516343
27	0640906	Button head M5x6 FZB	PCS	1	5707292401908
28	134225-1	Handle w seat tilt pictogram	PCS	1	5707292830142
29-30-31	SWIVEL FUNKTION 1-2	Wombat Solo str 1-2	PCS	1	5707292830487
32	150680	Bush f. attaching plate	PCS	4	5707292918246
33	160106	Bush for internal tube	PCS	1	5707292933720


Wombat solo size 3 lift/tilt on frame



Item numbers for Wombat solo size 3 lift/tilt on frame

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-100CBLTF	Solo stel m 100 mm wheel+CB	PCS	1	5707292813329
	134307-100LTFSF	Solo stel m 100 mm wheel	PCS	1	5707292875112
	134307-075LTFS	Solo stel m 75 mm wheel	PCS	1	5707292813015
2	139038-LFBTFBW	Cloth house compl str3, SWIVEL	PCS	1	5707292837370
3	1606500	Lock nut 6 mm	PCS	8	5707292141361
4	1103000	Fazet washer 6mm	PCS	8	5707292140289
5	135352	Seat Wombat Solo str 3	PCS	1	5707292825407
6	1101950	Washer Ø6,1-16x1 mm, FZB	PCS	8	5707292140210
7	0602401	Insex screw 6x30	PCS	8	5707292138125
8*	121391-7	back size 3 with high back	PCS	1	5707292784841
	121391	Wombat Living back size 3	PCS	1	5707292624161
9	0387000	End nipple for outer casing	PCS	1	5707292136091
10	0388006	Screwnipple 5x7,5mm	PCS	1	5707292136114
11	0380060	Flexible v-brake pipe.	PCS	1	5707292425195
12	0371853	Strained wire, 1,5 mm	PCS	1	5707292135889
13	1760005	Compensating nut 5mm	PCS	1	5707292141491
14	0380000	Cable tighter M5x23	PCS	1	5707292135933
15	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
16	127592	Cable fix	PCS	1	5707292810540
17	0609737	Hexalobular (6 Lobe) flat head	PCS	2	5707292828743
18	0235075	Cable tape, black	PCS	2	5707292135308
19*	139214-3	Footrest Simple frame size 3	PCS	1	5707292869647
	5002750	Cable for brake, size 1+2	PCS	1	5707292413635
20	0371853	Strained wire, 1,5 mm	PCS	1	5707292135889
22	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
*Choose size, color, type					


Wombat Solo size 3 lift on frame



Item numbers for Wombat Solo size 3 lift on frame

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-075LFS-1	Soloframe 75 mm hjul	PCS	1	5707292920652
	134307-100LFS-1	Solo stel m 100 mm hjul	PCS	1	5707292920799
	134307-100CBLF-1	Solo stel m 100 mm hjul+CB	PCS	1	5707292920713
	134307-100CBLFF-1	Solo stel m 100 mm hjul+CB	PCS	1	5707292920720
	134307-075LFS-1	Soloframe 75 mm hjul	PCS	1	5707292920652
	134307-075LFSF-1	Soloframe 75 mm hjul	PCS	1	5707292920669
	134307-100LFS-1	Solo stel m 100 mm hjul	PCS	1	5707292920799
	134307-100CBLF-1	Solo stel m 100 mm hjul+CB	PCS	1	5707292920713
	134307-100CBLFF-1	Solo stel m 100 mm hjul+CB	PCS	1	5707292920720
2*	139038-LFTSW1	Cloth house complete str 3	PCS	1	5707292920997
	139038-LFTS1	Cloth house complete str 3	PCS	1	5707292920980
3	1606500	Lock nut 6 mm	PCS	8	5707292141361
4	1103000	Fazet washer 6mm	PCS	8	5707292140289
5	135352	Seat Wombat Solo 3 complete	PCS	1	5707292825407
6	1101950	Washer Ø6,1-16*1 mm, FZB	PCS	8	5707292140210
7	602401	Insex screw 6*30	PCS	8	5707292138125
8*	121391	Wombat Living back size 3	PCS	1	5707292624161
	121391-7	back size 3 with high back	PCS	1	5707292784841
9	387000	End nipple for outer casing	PCS	1	5707292136091
10	388006	Screwnipple 5*7,5mm	PCS	1	5707292136114
11	380060	Flexible v-brake pipe.	PCS	1	5707292425195
12	371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
13	1760005	Compensating nut 5mm	PCS	1	5707292141491
14	380000	Cable tighter M5*18	PCS	1	5707292135933
15	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
16	127592	Cable fix	PCS	1	5707292810540
17	609737	Hexalobular (6 Lobe) flat head	PCS	2	5707292828743
18*	139214-3	Footrest Simple frame size 3	PCS	1	5707292869647
	235075	Cable tie, black, 3,6*200mm	PCS	1	5707292135308
19*	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
	5002750	Cable for brake, size 1+2	PCS	1	5707292413635
20	371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
22	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
*Choose size, color, type					


Wombat Solo size 3 lift and tilt on brace



Item numbers for Wombat Solo size 3 lift and tilt on brace

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-100CBSF	Solo Frame std m.100mm wheel+CB	PCS	1	5707292875129
	134307-100CBS	Solo Frame std m.100mm wheel+CB	PCS	1	5707292813060
2*	139038-LFBTFBW	Cloth house compl str3, SWIVEL	PCS	1	5707292837370
	139038-LFBTFB	Cloth house complete str 3	PCS	1	5707292811516
3	1606500	Lock nut 6 mm	PCS	8	5707292141361
4	1103000	Fazet washer 6mm	PCS	8	5707292140289
5	135352	Seat Wombat Solo str 3	PCS	1	5707292825407
6	1101950	Washer Ø6,1-16x1 mm, FZB	PCS	8	5707292140210
7	0602401	Insex screw 6x30	PCS	8	5707292138125
8*	121391-7	back size 3 with high back	PCS	1	5707292784841
	121391	Wombat Living back size 3	PCS	1	5707292624161
9	3503534	Gas pump lever	PCS	2	5707292845634
10	0371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
11	0380060	Flexible v-brake pipe.	PCS	2	5707292425195
12	5002756	Cable for recline, outer 1-4	PCS	1	5707292413673
13	0235075	Cable tape, black	PCS	2	5707292135308
14	0371853	Strained wires 1,5mm 190 cm.	PCS	1	5707292135889
15	5002750	Cable for brake, size 1+2	PCS	1	5707292413635
16	121218	Håndbøjle str 2 complete	PCS	1	5707292624178
17	0387000	End nipple for outer casing	PCS	2	5707292136091
18	0388006	Screwnipple 5x7,5mm	PCS	2	5707292136114
19	139214-3	Footrest Simple frame size 3	PCS	1	5707292869647
20	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
*Choose size, color, type					


Wombat Solo size 3



Item numbers for Wombat Solo size 3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-100CBS	Solo Frame std m.100mm wheel+CB	PCS	1	5707292813060
	134307-075SF	Solo Frame m. 75 mm wheel	PCS	1	5707292875068
	134307-100CBSF	Solo Frame std m.100mm wheel+CB	PCS	1	5707292875129
	134307-100SF	Solo Frame std m. 100 mm wheel	PCS	1	5707292875099
	134307-075S	Solo Frame std	PCS	1	5707292811134
2*	139038-TLSW	Cloth house complete size 3	PCS	1	5707292837356
	139038-TLS	Cloth house complete size 3	PCS	1	5707292811493
3	1606500	Lock nut 6 mm	PCS	8	5707292141361
4	1103000	Fazet washer 6mm	PCS	8	5707292140289
5	135352	Seat Wombat Solo str 3	PCS	1	5707292825407
6	1101950	Washer Ø6,1-16x1 mm, FZB	PCS	8	5707292140210
7	0602401	Insex screw 6x30	PCS	8	5707292138125
8*	121391-7	back size 3 with high back	PCS	1	5707292784841
	121391	Wombat Living back size 3	PCS	1	5707292624161
9	139214-3	Footrest Simple frame size 3	PCS	1	5707292869647
10	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
*Choose size, color, type					


Wombat Solo size 3, no tilt



Item numbers for Wombat Solo size 3, no tilt

Position	No.	Description	Unit of Measure	Quantity	Ean code
1*	134307-100CBS	Solo Frame std m.100mm wheel+CB	PCS	1	5707292813060
	134307-100CBSF	Solo Frame std m.100mm wheel+CB	PCS	1	5707292875129
	134307-100SF	Solo Frame std m. 100 mm wheel	PCS	1	5707292875099
	134307-075SF	Solo Frame m. 75 mm wheel	PCS	1	5707292875068
2*	139038-SW	std. løft på sæde m. swivel	PCS	1	5707292869326
	139038-S	Cloth house complete str 3	PCS	1	5707292869319
3	1606500	Lock nut 6 mm	PCS	8	5707292141361
4	1103000	Fazet washer 6mm	PCS	8	5707292140289
5	135352	Seat Wombat Solo str 3	PCS	1	5707292825407
6	1101950	Washer Ø6,1-16x1 mm, FZB	PCS	8	5707292140210
7	0602401	Insex screw 6x30	PCS	8	5707292138125
8*	121391-7	back size 3 with high back	PCS	1	5707292784841
	121391	Wombat Living back size 3	PCS	1	5707292624161
9	139214-3	Footrest Simple frame size 3	PCS	1	5707292869647
10	135020-3	Footrest size 3 Wombat Solo	PCS	1	5707292813688
*Choose size, color, type					


Wombat Living back size 3



Item numbers for Wombat Living back size 3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	1940010	Plastite 60, 2,9x6,5	PCS	6	5707292142573
2	4913040	30mm adhesive hook, black	MTR	0,62	5707292786517
3	82794	Backplate f. Wombat size 3	PCS	1	5707292400116
4	1940032	Eco-sight z M3,5x12 rustproof	PCS	8	5707292397584
5	82781	Mounting plate	PCS	2	5707292398840
6	82751	Back mounted with screws	PCS	1	5707292400123
7	2501500	Ferrule d12x0,5-1,5 int.black	PCS	1	5707292143952
8	3503001	R82 handle M8x15mm screw	PCS	1	5707292413925
9	0445025	25mm press buckle, black	PCS	4	5707292465160
10	0640920	Button head M5x20 FZB	PCS	4	5707292381033
11	3502300	Tightening handle 65mm M8 bush	PCS	1	5707292146885
12.1	127042	Backpile	PCS	1	5707292784858
12	172507	Backpile	PCS	1	5707292990365
13	82747	Shaft	PCS	1	5707292400130
14	1103200	Fazet washer 8mm	PCS	2	5707292140296
15	1907100	Steel bolt 8x70mm	PCS	1	5707292142467
16	0310163	Compression spring 1,3x12rustpr	PCS	1	5707292402790
17	82796-60	Lockcover, black	PCS	1	5707292619846
18	82795-60	Lockgrip, black	PCS	1	5707292619839
19	0530005	Cyl.pin Ø8x66 mm	PCS	2	5707292397805
20	0603975	Insex screw 8x110	PCS	1	5707292138330
21	82797	Assembl fitting seat/back	PCS	1	5707292400048
22	1607000	Lock nut 8mm, large pack	PCS	1	5707292141378
23	1606500	Lock nut 6 mm	PCS	1	5707292141361
24	0602188	Insex M6x8, FZB	PCS	1	5707292138071
25	0683925	Washer 8x14x1	PCS	2	5707292467980
26	1100495	Spring washer M8	PCS	1	5707292458674
*Choose size, color, type					


Wombat seat 3 for solo



Item numbers for Wombat seat 3 for solo

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	1606500	Lock nut 6 mm	PCS	4	5707292141361
2	0641600	Button head screw 6x16	PCS	2	5707292139269
3	0442070	Fixlock clamp 40 mm w/hole	PCS	2	5707292136398
12	82739-2	Seat sz3 Wombat, bushings,nuts	PCS	1	5707292804754
13	4935300	100mm black velcro hook	MTR	0,27	5707292152107
14	0506509	Rivet 4x10, stainless	PCS	12	5707292472441
15	4935300	100mm black velcro hook	MTR	0,15	5707292152107
16	1106000	Washer 6mm	PCS	4	5707292140371
20	2203322	Bloc-o-lift 684074 200n solid	PCS	1	5707292398147
21	1760800	Compensating nut 8x1,0 finep.t	PCS	1	5707292141521
22	82741	Shaft	PCS	1	5707292400109
23	2521840	Pipe holder f/16mmsquare,black	PCS	2	5707292399236
24	0602435	Insex screw 6x35	PCS	4	5707292138132
25	138381	Back tip	PCS	1	5707292868787
26	0530003	Cyl.pin Ø6x32 mm	PCS	1	5707292397638
27	1940025	Pan head screw 4x10 PZ FZB	PCS	2	5707292545374
28	112307	Lockplate	PCS	2	5707292591722
31	0445025	25mm press buckle, black	PCS	2	5707292465160
32		Button head M5x30 FZB	PCS	2	5707292139221
41	0602200	Insex screw 6x16	PCS	2	5707292138095
42	2502510	Plug 20x2,0-2,5 inside black	PCS	4	5707292557759


Clutch house complete Wombat size 3



Item numbers for Clutch house complete Wombat size 3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	1606300	Lock nut 5 mm	PCS	3	5707292141354
2	136903	Insert Cap	PCS	1	5707292834270
3	136547	Cover size 3	PCS	1	5707292831354
4	121245	Locking pin wombat solo	PCS	1	5707292810786
5	136440	Column assembly	PCS	1	5707292829825
6*	0631811	Insex uh M5x16 TufLok® Patch	PCS	4	5707292884114
	0615022	M5x20 stainless hexagon screw	PCS	4	5707292837424
7*	0609726	Insex 5x30 uh. Fzb	PCS	1	5707292433381
	0615025	Insex screw 5x25	PCS	1	5707292138873
8*	127728	Attaching plate	PCS	1	5707292787057
	150677	Attaching plate	PCS	1	5707292918239
9	127725	Bearing disc	PCS	1	5707292794222
10	135987	Cover	PCS	1	5707292829870
11	2200550	Bloc-o-lift 350 N	PCS	1	5707292830166
12	1761000	Counter nut 10x1,0 fine thread	PCS	1	5707292141538
13	1103400	Fazet washer 10 mm	PCS	1	5707292140319
14	127836	Sliding dish	PCS	1	5707292794239
15	0640935	Button head M5x35 stainless	PCS	1	5707292485748
16	138495	Release handle for lift	PCS	1	5707292842138
17	128675-1	Handle w. seat lift pictogram	PCS	1	5707292830128
18*	139095-1	Cloth house Wombat Solo str 3	PCS	1	5707292920928
	139095	Cloth house Wombat Solo str 3	PCS	1	5707292868947
20	0380000	Cable tighter M5x23	PCS	1	5707292135933
21	127996	Bracket for tilt	PCS	1	5707290000000
22	658012	Button head M5x12 FZB	PCS	2	5710000000000
23*	138348	Cloth house Wombat Solo str 3	PCS	1	5707292824707
	151079	Cloth house Wombat Solo str 3	PCS	1	5707292917737
24	138499	Fork Bracket	PCS	1	5707292842213
25	128261	Spring for tilt	PCS	1	5707292810809
26	1130500	Stainless washer M5, din9021	PCS	1	5707292516343
27	0640906	Button head M5x6 FZB	PCS	1	5707292401908
28	134225-1	Handle w seat tilt pictogram	PCS	1	5707292830142
29-30-31	SWIVEL FUNKTION 3	SWIVEL FUNKTION 3	PCS	1	5707292878878
32	150680	Bush f. attaching plate	PCS	4	5707292918246
33	160106	Bush for internal tube	PCB	1	5707292933720
*Choose size, color, type					


Frame assembly



Item numbers for Frame assembly

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	3006820	75mm black steinco wheel	PCS	4	5707292368553
2	1103400	Fazet washer 10 mm	PCS	4	5707292140319
3	127726-1	Release pedal	PCS	1	5707292810649
4	136102	Frame with cover plugs	PCS	1	5707292829122
5	1606300	Lock nut 5 mm	PCS	2	5707292141354
6	1102900	Fazet washer 5 mm	PCS	6	5707292140272
6.1	1100485	Spring washer M5, stainless	PCS	4	5707292478603
7	136009-1	Release pedal	PCS	1	5707292829139
8	0631955	Stainless button head M5x12	PCS	2	5707292139207
9	128090	Insert cover	PCS	1	5707292810533
11	0602454	Insex M5x40 stainless	PCS	4	5707292915214
12	127262	Taper ring	PCS	1	5707292810748
13	128070	Lock ring	PCS	1	5707292810731
14	3006938	100mm double wheel with brake	PCS	2	5707292574930
16	0602454-1	Insex M5x45 stainless	PCS	4	5707292927606


Footrest frame, size 1-2-3



Item numbers for Footrest frame, size 1-2-3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	0661800	Button head w/flange 8x20	PCS	2	5707292408075
2	138365	Footplate sz 1+2	PCS	1	5707292842299
3	135670	Foot Support profile	PCS	1	5707292813664
4	2504500	Ferrule e20x1,0-2,0 int.black	PCS	1	5707292144126
5	3502516	Tightening handle 45mm M8x16	PCS	1	5707292146915
6	1940017	Sheet metal screw uhx4,2x13fzb	PCS	4	5707292142597
7	82950	Hollow plate 16x32, Wombat	PCS	2	5707292126030
8	166392	Footrest Low, rotatable	PCS	1	5707292982056
9	1604000	Steel nut 10 mm	PCS	1	5707292141293
10	3502780	Handle 28mm,M10x160mm thread	PCS	1	5707292363145
11	82944	Footplate size 1, Wombat	PCS	1	5707292125996
12	3502790	Hand grip 28mm, M10x196mm	PCS	1	5707292363152
13	82940	Footplate size.2, Wombat	PCS	1	5707292125958
14	3502795	Knob 28mm,M10x??mm thread	PCS	1	5707292402110
15	82752	Footplate size 3, Wombat	PCS	1	5707292401373
16	138349	Footplate sz 3	PCS	1	5707292842282
18	166635	Footrest Long, rotatable	PCS	1	5707292982063


Footrest seat, size 1-2-3



Item numbers for Footrest seat, size 1-2-3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	3502516	Tightening handle 45mm M8x16	PCS	1	5707292146915
2	1940017	Sheet metal screw uhx4,2x13fzb	PCS	4	5707292142597
3	82950	Hollow plate 16x32, Wombat	PCS	2	5707292126030
4	1604000	Steel nut 10 mm	PCS	1	5707292141293
5	166392	Footrest Low, rotatable	PCS	1	5707292982056
6	3502780	Handle 28mm,M10x160mm thread	PCS	1	5707292363145
7	82944	Footplate size 1, Wombat	PCS	1	5707292125996
8	2504500	Ferrule e20x1,0-2,0 int.black	PCS	1	5707292144126
9	82927	Pipe for footplate, Wombat	PCS	1	5707292125828
10	3502300	Tightening handle 65mm M8 bush	PCS	1	5707292146885
11	120701	Lock washer	PCS	1	5707292623416
12	82934	Distance piece	PCS	2	5707292125897
13	1607000	Lock nut 8mm, large pack	PCS	1	5707292141378
14	1103200	Fazet washer 8mm	PCS	1	5707292140296
15	0701000	Stop screw with crater 8x10	PCS	1	5707292139870
16	1350022	Distance piece 12x15x8,2mm	PCS	2	5707292140906
17	82926	Turning fitting	PCS	1	5707292125811
18	1103500	Fazet washer 12mm	PCS	1	5707292140326
19	0603665	Insex screw 8x65	PCS	1	5707292138286
20	1905710	Bolt M8x70 FZB wo/nut	PCS	1	5707292142368
21	82734	Plate complete	PCS	1	5707292398932
22	74061	Tooth washer Ø30	PCS	2	5707292373816
23	82735	Fix for footrest	PCS	1	5707292399335
24	1101000	Center spring washer 8mm	PCS	4	5707292140173
25	0661145	Button head 8x45 FZB	PCS	4	5707292139535
26	82787	Rod for Wombat footrest 280 mm	PCS	1	5707292401472
27	0603660	Insex screw 8x60	PCS	1	5707292138279
28	82752	Footplate size 3, Wombat	PCS	1	5707292401373
29	3502795	Knob 28mm,M10x??mm thread	PCS	1	5707292402110
30	82940	Footplate size.2, Wombat	PCS	1	5707292125958
31	3502790	Hand grip 28mm, M10x196mm	PCS	1	5707292363152
32	166635	Footrest Long, rotatable	PCS	1	5707292982063

Fixed side supports, with covers (95121, 95121-S & 95121-23)



Item numbers for Fixed side supports, with covers (95121)

Position	No.	Description	Quantity	Unit of Measure	Ean code
1+2	121396	betræk til sidesupport	1	SET	5707292632050
3	4913040	30mm adhesive hook, black	0,2	MTR	5707292786517
4	82875	Sidesupports for Wombat	2	PCS	5707292364302
5+6	5001021	Screws etc. For side supports	1	SET	5707292369475

Item numbers for Fixed side supports, with covers (95121-S)

Position	No.	Description	Quantity	Unit of Measure	Ean code
1+2	121396	betræk til sidesupport	1	SET	5707292632050
3	4913040	30mm adhesive hook, black	0,2	MTR	5707292786517
4	82875	Sidesupports for Wombat	2	PCS	5707292364302
5+6	5001049	Screws etc. For side support	1	PCS	5707292969712

Item numbers for Fixed side supports, with covers (95121-23)

Position	No.	Description	Quantity	Unit of Measure	Ean code
1+2	171071	Inkobetræk t.faste sidestøtter	1	SET	5707292988881
3	4913040	30mm adhesive hook, black	0,2	MTR	5707292786517
4	82875	Sidesupports for Wombat	2	PCS	5707292364302
5+6	5001021	Screws etc. For side supports	1	SET	5707292369475