


x:panda Accessories

Exploded views / Part lists

Rev. date 2024.05


On the following pages you will find the parts used to assemble this product. Please note that some spare parts cannot be purchased individually but are part of a complete kit. Please also note that for some spare parts extra delivery time may occur.



Contents

Armrests size 1-3	3
Item numbers for Armrests size 1-3	3
Armrests size 4 dual rod	4
Item numbers for Armrests size 4 dual rod	4
Armrests size 4 single rod	5
Item numbers for Armrests size 4 single rod	6
Backrest extension.....	7
Item numbers for Backrest extension.....	7
Upper arm supports	8
Item numbers for Upper arm supports	9
Swing away shoulder supports.....	10
Item numbers for Swing away shoulder supports.....	10
Side supports	11
Item numbers for Side supports	12
Hip support - single rod - size 4.....	13
Item numbers for Hip support - single rod - size 4.....	14
Hip support - dual rod - size 4	15
Item numbers for Hip support - dual rod - size 4	16
Pommel removable	17
Item numbers for Pommel removable	17
Knee abduction pommel	18
Item numbers for Knee abduction pommel	18
Swing-away tray, Two rod system tray, armrests, toybar and grip handle.....	19
Item numbers for Swing-away tray, Two rod system, armrests, toybar and grip handle.....	20
Two rod system tray - Size 4	22
Item numbers for Two rod system tray - Size 4	22
Tray plast.....	23
Item numbers for Tray plast.....	23
Floor Frame	24
Item numbers for Floor Frame.....	25
Floor Frame accessories.....	26
Item numbers for Floor Frame accessories.....	27


Armrests size 1-3



Item numbers for Armrests size 1-3

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	116980	Release housing, left	PCS	1	5707292786395
2	120258	Release housing, Right	PCS	1	5707292786401
3	118979	Armrest size 1 left	PCS	1	5707292785572
3	128458	Armrest size 1 Tooth bar s.2 L	PCS	1	5707292810267
3	117344	Armrest size 2, L tooth bar s3	PCS	1	5707292786074
4	118992	Armrest size 1 right	PCS	1	5707292786067
4	128457	Armrest size 1 Tooth bar s.2 r	PCS	1	5707292810250
4	116760	Armrest Size 2, right	PCS	1	5707292785893


Armrests size 4 dual rod



Item numbers for Armrests size 4 dual rod

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	6175163	Armrest pu. Wombat sz 4. Right	PCS	1	5707292443151
2	82705	Armrest right size 3	PCS	1	5707292442734
3	75443	Mounting for table, right	PCS	1	5707292517173
4	82666-1	Armrest, wombat, extended	PCS	4	5707292546906
4	82680	Armrest rod ext.	PCS	4	5707292531148
5	75438	Armrest bracket, right	PCS	1	5707292517227
6	5082923	Bush for arm rest x:panda 4	PCS	4	5707292570505
7	75439	Armrest bracket, left	PCS	1	5707292517265
8	0661800	Button head w/flange 8*20	PCS	4	5707292408075
12	9500001	Lock Wombat	PCS	4	5707292404527
14	3652501	R82 handle M6*16mm screw	PCS	2	5707292413963
15	0611000	Insex screw counters. 6*20	PCS	4	5707292138682
16	75442	Mounting for table, left	PCS	1	5707292517203
17	0641620	Button head screw 6*20	PCS	4	5707292139276
18	82718	Armrest bracket left size 3	PCS	1	5707292442727
19	1310640	Washer, nature Ø16*10*0,5mm	PCS	4	5707292140760
20	1606500	Lock nut 6 mm	PCS	4	5707292141361
21	6175164	Armrest pu. Wombat sz 4. Left	PCS	1	5707292443168
22	1607000	Lock nut 8mm, large pack	PCS	2	5707292141378
23	75306	Snap lock R82 25	PCS	2	5707292437624
24	0640911	Button head M5x10 FZB tuflok	PCS	2	5707292428479


Armrests size 4 single rod



Item numbers for Armrests size 4 single rod

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	103147	Holder for armrest, right	PCS	1	5707292566089
2	103190	Tube for elevation right	PCS	1	5707292566096
3	75392	Holder for armrest, right	PCS	1	5707292454775
4	6175163	Armrest pu. Wombat sz 4. Right	PCS	1	5707292443151
5	6175164	Armrest pu. Wombat sz 4. Left	PCS	1	5707292443168
6	1310539	Distance piece 12/8,4*5mm	PCS	4	5707292423474
7	75391	Holder for armrest, left	PCS	1	5707292454768
8	0659580	Button head screw 6*30 mm	PCS	2	5707292139481
9	75395	Lock, armrest & table	PCS	2	5707292454799
10	103191	Tube for elevation left	PCS	1	5707292566102
11	0310113	Stainless spring 1,2*16indv.*20	PCS	2	5707292453143
12	3501635	R82 handle dia.,M8*37	PCS	2	5707292451910
13	103126	Holder for armrest, left	PCS	1	5707292566072
14	0602425	Insex M6*25, stainless	PCS	4	5707292404596
15	1605001	M8-1,25 stainless steel	PCS	4	5707292412454
16	103016	Bracket for side support	PCS	2	5707292566065
17	1105060	Stainless washer 6mm	PCS	4	5707292140357
18	0602207	Insex M6*18, stainless	PCS	4	5707292430380
19	1310645	Washer, nature Ø16*8*0,5mm	PCS	4	5707292446015
20	74063	Toothed disk ø30*9	PCS	4	5707292377968
21	0662000	Button head screw 8*20 stainl.	PCS	4	5707292139573
22	74784	Release f. Armrest	PCS	2	5707292419903
23	0532006	Stainless spring pin M5*22	PCS	2	5707292570215
24	0588014	Rustl.tightening spring 3*14	PCS	2	5707292137753
25	74681	Locking pin, armrest	PCS	2	5707292419132
26	5020101	Padded right armrest,x:panda 4	PCS	1	5707292566188
27	5020106	Padded left armrest,x:panda 4	PCS	1	5707292566195
28	0310100	Pressure spring 0,9*8*23*10,5v	PCS	2	5707292135513


Backrest extension



Item numbers for Backrest extension

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	8910652-160	Cover for back rest extension	STK	1	5707292542038
1	8910652-260	Cover for back rest extension	STK	1	5707292542045
1	8910652-360	Cover for back rest extension	STK	1	5707292542052
2	5031930	Padded back extens x-panda s 1	STK	1	5707292542410
2	5031931	Padded back exten. x-panda s.2	STK	1	5707292542458
2	5031932	Padded back exten. x-panda s.3	STK	1	5707292542465
3	76368	back extension bracket	STK	1	5707292542175
4	0609728	Insex, stainless 6*16 uh	STK	2	5707292138491
5	0603530	Insex low head M8x16 stainless	STK	1	5707292412492
6	75306	Snap lock R82 25	STK	2	5707292437624
7	0658010	Button head M5 x 10, stainless	STK	2	5707292139375
8	81155	Rod for backextension	STK	1	5707292112781
9	2504000	Ferrule e15*1,0-2,0 int.black	STK	1	5707292144102


Upper arm supports



Item numbers for Upper arm supports

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	117067	Joint f. upper arm sup. right	PCS	1	5707292613813
2	117066	Joint f. upper arm sup. right	PCS	1	5707292613806
3	0602200	Insex screw 6*16	PCS	2	5707292138095
4	114962	Joint f. upper arm sup. left	PCS	1	5707292613776
5	1103000	Fazet washer 6mm	PCS	4	5707292140289
6	0602190	Insex screw 6*12	PCS	4	5707292138088
7	0662500	Button head screw 8*25 stainl.	PCS	4	5707292139580
7.1	0661200	Button head screw 8*12 stainl.	PCS	4	5707292139542
8	1606500	Lock nut 6 mm	PCS	2	5707292141361
9	1105080	Stainless washer 8mm	PCS	4	5707292140364
10	6118920	M6 Nut	PCS	4	5707292433237
11	114970	Joint f. upper arm sup. left	PCS	1	5707292613820
12	117850	Padding upper arm sup. size2 L	PCS	1	5707292614018
12	117861	Padding upper arm sup. size1 L	PCS	1	5707292613981
13	8910407-20260L	Cover for adj. upper arm supp.	PCS	1	5707292788115
13	8910407-10160L	Cover for adj. upper arm supp.	PCS	1	5707292788078
14	8910407-20260R	Cover for adj. upper arm supp.	PCS	1	5707292788122
14	8910407-10160R	Cover for adj. upper arm supp.	PCS	1	5707292788085
15	117849	Padding upper arm sup. size2 R	PCS	1	5707292614056
15	117857	Padding upper arm sup. size1 R	PCS	1	5707292613998


Swing away shoulder supports



Item numbers for Swing away shoulder supports

Position	No.	Description	Unit of Measure	Quantity	Ean code
202	0602900	Insex screw 8*16	PCS	1	5707292138200
203	3706256	Knob 25, M6 thread	PCS	2	5707292147219
204	1100200	External snap ring	PCS	2	5707292140128
205	6108224	Sringhouse shoulder support	PCS	2	5707292185754
206	0310145	Pressure spring 0,9*8,6mm*15,5	PCS	2	5707292135612
207	6118800	Pin f/springhouse impala+panda	PCS	2	5707292186683
208	6109020	Holder f.shoulder support	PCS	2	5707292185938
209	0808008	Stainless stop screw 8*8	PCS	2	5707292140036
210	5700108	Glueing shoulder support-panda	PCS	2	5707292179517
211	82361	Mounting for shoulder support	PCS	1	5707292123657
211	82362	Mounting for shoulder support	PCS	1	5707292123664



Side supports



Item numbers for Side supports

Position	No.	Description	Unit of Measure	Quantity	Ean code
201+202	5700522-2	Complete backside support - 2	SET	1	5707292810472
201	5700522-2	Complete backside support - 2	SET	1	5707292810472
202	5700522-2	Complete backside support - 2	SET	1	5707292810472
203L	127527	Cover for left side supp. T.2	PCS	1	5707292807021
203R	127529	Cover for right side supp. T.2	PCS	1	5707292807038
204+205	5700520-2	Complete backside support - 1	SET	1	5707292810496
204	5700520-2	Complete backside support - 1	SET	1	5707292810496
205	5700520-2	Complete backside support - 1	SET	1	5707292810496
206L	127515	Cover for left side supp. T.1	PCS	1	5707292807007
206R	127517	Cover for right side supp. T.1	PCS	1	5707292807014
213+214	5700524-2	Complete backside support - 3	SET	1	5707292810489
213	5700524-2	Complete backside support - 3	SET	1	5707292810489
214	5700524-2	Complete backside support - 3	SET	1	5707292810489
215L	127531	Cover for left side supp. T.3	PCS	1	5707292807052
215R	127533	Cover for right side supp. T.3	PCS	1	5707292807069
216+217	5700526-2	Complete backside support - 4	SET	1	5707292810427
216	5700526-2	Complete backside support - 4	SET	1	5707292810427
217	5700526-2	Complete backside support - 4	SET	1	5707292810427
218L	127519	Cover for left side supp. T.4	PCS	1	5707292807076
218R	127521	Cover for right side supp. T.4	PCS	1	5707292807083
219	75203-1RAS	M8 slide with threads	PCS	4	5707292809841
220	0662501	Button head w/flange 8*25	PCS	4	5707292403759
230L	127535	Cover for left side supp. T.5	PCS	1	5707292807106
230R	127537	Cover for right side supp. T.5	PCS	1	5707292807113
231	5700528	cover back side size 5 left	PCS	1	5707292514004
232	5700529	Upholst back side supp 5 right	PCS	1	5707292514011
295	5700566	backside mounted size 6 R	PCS	1	5707292538413
295	5700565	backside mounted, size 6 L	PCS	1	5707292538406
296	5700565	backside mounted, size 6 L	PCS	1	5707292538406
296	5700566	backside mounted size 6 R	PCS	1	5707292538413
297L	127523	Cover for left side supp. T.6	PCS	1	5707292807137
297R	127525	Cover for right side supp. T.6	PCS	1	5707292807151
300	5027198	Sewn cover back s.5sp. r+l	PCS	2	5707292551337
301	5700519	Oplime rygsidest.str. 5 høj sp	PCS	1	5707292551320
302	5700518	Oplime rygsidest.str. 5 ven sp	PCS	1	5707292551245
304	5027196	Sewn cover back s.3sp. r+l	PCS	2	5707292572547
305+307	5700508	Oplime rygsidest.str. 3 spec.	SET	1	5707292572523
308	5700568	Oplime rygsidest.str. 6 h.spec	PCS	1	5707292574336
309	5027209	Sewn cover backsides s.6 spec	PCS	2	5707292574343
310	5700567	Oplime rygsidest.str. 6 v.spec	PCS	1	5707292574329


Hip support - single rod - size 4



Item numbers for Hip support - single rod - size 4

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	75490	Adjustable toothdisc Left	PCS	1	5707292541369
2	76171	Crank for side support	PCS	2	5707292491602
3	5082923	Bush for arm rest x:panda 4	PCS	2	5707292570505
5	9500001	Lock Wombat	PCS	2	5707292404527
6	75491	Adjustable toothdisc Right	PCS	1	5707292541376
10	76183	Handle M6 with bush	PCS	4	5707292505859
11	76201	Stay with M6	PCS	4	5707292503992
12	76188	Handle eccentric	PCS	4	5707292503350
13	0529045	Pin 4*15mm, stainless 303	PCS	4	5707292509833
14	76189	Saddle for excentrik handle	PCS	4	5707292503374
15	74905	Toothdisc bracket	PCS	2	5707292499080
16	82974	Threaded plate for hip support	PCS	2	5707292126184
17	6175171	Cushion size 1	PCS	2	5707292189301
18	82953	Plate for hip support	PCS	2	5707292126047
19	0609728	Insex, stainless 6*16 uh	PCS	4	5707292138491
20	0658010	Button head M5 x 10, stainless	PCS	4	5707292139375
21	1100485	Spring washer M5, stainless	PCS	4	5707292478603


Hip support - dual rod - size 4



Item numbers for Hip support - dual rod - size 4

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	0529045	Pin 4*15mm, stainless 303	PCS	4	5707292509833
2	76188	Handle eccentric	PCS	4	5707292503350
3	76189	Saddle for excentrik handle	PCS	4	5707292503374
4	76201	Stay with M6	PCS	4	5707292503992
5	103218	Mounting for hipsupport	PCS	2	5707292565914
6	76183	Handle M6 with bush	PCS	4	5707292505859
7	103223	Crank for side support	PCS	2	5707292565907
8	1100490	Spring washer M6 stainless	PCS	4	5707292411297
9	0659200	Button head screw 6*12 stainl.	PCS	4	5707292139436
10	0658010	Button head M5 x 10, stainless	PCS	4	5707292139375
11	1100485	Spring washer M5, stainless	PCS	4	5707292478603
12	74905	Toothdisc bracket	PCS	2	5707292499080
13	82953	Plate for hip support	PCS	2	5707292126047
14	0609728	Insex, stainless 6*16 uh	PCS	4	5707292138491
15	82974	Threaded plate for hip support	PCS	2	5707292126184
16	6175172	Cushion size 2	PCS	2	5707292189318


Pommel removable



Item numbers for Pommel removable

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	0602486	Insex M8*16	PCS	1	5707292434272
2	75206	Bracket for abduction block	PCS	1	5707292442284
3	0603530	Insex low head M8x16 stainless	PCS	1	5707292412492
4	101054	Seat piece	PCS	1	5707292554758
5	0602195	Insex M6*12 stainless	PCS	2	5707292404701
6	1105060	Stainless washer 6mm	PCS	2	5707292140357
7	101055	Pommel piece	PCS	1	5707292554765
8	6175110	Pommel turtle, black pur	PCS	1	5707292466624
8	6175111	Pommel size 1, black pur	PCS	1	5707292189103
8	6175112	Pommel size 2, black pur	PCS	1	5707292189110
8	6175113	Pommel sz. 3, black pur	PCS	1	5707292189127
8	6175114	Pommel size 4, black pur	PCS	1	5707292189134
9	1811275	Cutting screw 4*30mm	PCS	2	5707292141743


Knee abduction pommel



Item numbers for Knee abduction pommel

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	82978	Pommel brace Wombat	PCS	1	5707292126221
2	1811275	Cutting screw 4*30mm	PCS	2	5707292141743
3	3502516	Tightening handle 45mm M8*16	PCS	1	5707292146915
4	119879	Mount for pommel system	PCS	1	5707292621320
5	1100500	Spring washer 8mm, stainless	PCS	2	5707292140142
6	0661500	Button head screw 8*16 stainl.	PCS	2	5707292139559
7	113633	Pommel with foam	PCS	1	5707292600882
8	95430-10160	Cover f. knee abduction pommel	PCS	1	5707292600646

Swing-away tray, Two rod system tray, armrests, toybar and grip handle




R82 nothing compares to a smile R82 A/S Parallelvej 3 8751 Gedved Denmark	Created by JYBA	Revision 3	Revision date 04-12-2020
	Project X-panda	Memo 79313	Table

Item numbers for Swing-away tray, Two rod system, armrests, toybar and grip handle

Position	No.	Description	Unit of Measure	Quantity	Ean code
120	75361	Table size 1	PCS	1	5707292455994
120	75362	Table size 2	PCS	1	5707292456298
123	75358	Cover plate for table 1, right	PCS	1	5707292456236
124	75359	Cover plate for table 1, left	PCS	1	5707292456243
126	1940038	Eco-syn uh Ø4*14	PCS	16	5707292465597
127	1631060	Stainless lock nut 6mm	PCS	2	5707292141446
128	0602435	Insex screw 6*35	PCS	2	5707292138132
135	5037981	Upholstered armrest, right	PCS	1	5707292456441
136	5037980	Upholstered armrest, left	PCS	1	5707292456434
137	75342	Bracket for armrest	PCS	2	5707292456052
138	0602195	Insex M6*12 stainless	PCS	4	5707292404701
139	1105060	Stainless washer 6mm	PCS	4	5707292140357
141	75330	Bracket	PCS	1	5707292456083
142	148835	Rocker shaft f. Table	PCS	1	5707292911629
143	148836	Rocker shaft	PCS	1	5707292911971
144	75334	Bush for rocker table	PCS	1	5707292456151
145	0659550	Button head M6*20, stainless	PCS	2	5707292139467
146	0609816	Stainless insex screw 8*16	PCS	2	5707292138569
147	2502002	Ferrule d18*0,8-2,5 black	PCS	2	5707292433947
148	0661200	Button head screw 8*12 stainl.	PCS	4	5707292139542
149	75338	Rocker shaft	PCS	1	5707292456168
150	0102500	O-ring 12,5*2 nbr70	PCS	2	5707292467362
151	1100500	Spring washer 8mm, stainless	PCS	4	5707292140142
152	1100490	Spring washer M6 stainless	PCS	2	5707292411297
154	1100495	Spring washer M8	PCS	2	5707292458674
155	75345	Aktiviti pipe	PCS	1	5707292456175
156	75363	Holder for activity pipe	PCS	2	5707292456328
157	73819	Safety pin	PCS	2	5707292441447
158	5505015	Strap for toy bar	PCS	2	5707292470171
159	3503001	R82 handle M8x15mm screw	PCS	2	5707292413925
160	0603534	Insex low head M8x20	PCS	2	5707292575302
182	6153700	Knob spring f.activity centre	PCS	4	5707292188601
186	82472	Mounting f.activity centre	PCS	2	5707292124494
187	82471	Mounting f.activity centre	PCS	2	5707292124487
188	1102000	Planing washer 8mm	PCS	2	5707292140227
189	3504000	Tightening handle 42mm M8*20	PCS	2	5707292147097
191	1603000	Steel nut 6 mm	PCS	2	5707292141279
192	0615500	Eye screw 6*16	PCS	2	5707292138897
193	82460	Brace, right, activity centre	PCS	1	5707292124401
194	82463	Intermidate, activity centre	PCS	1	5707292124432
195	2521500	Pipeholder Ø14/15, white	PCS	6	5707292144379
196	82450	Tableplate f. Activity centre	PCS	1	5707292124302
199	0621600	Insex screw counters. 6*16	PCS	4	5707292138965
201	82461	Brace, left, activity centre	PCS	1	5707292124418
300	75861	Joint part-table-right	PCS	1	5707292480583
301	75860	Joint part-table-left	PCS	1	5707292480576
302	1631080	Stainless lock nut 8mm	PCS	2	5707292141453
303	0603535	Insex low head M8*25	PCS	2	5707292575319
304	75857	Pipe for brace on tray	PCS	1	5707292480521

305	75859	Joint-rail-right/left	PCS	2	5707292480552
306	0603530	Insex low head M8x16 stainless	PCS	2	5707292412492
307	2500631	fex tubing CPE/EPDM	PCS	1	5707292574053
308	1100495	Spring washer M8	PCS	2	5707292458674


Two rod system tray - Size 4



Item numbers for Two rod system tray - Size 4

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	80221	Table for x:panda 4	PCS	1	5707292517326
2	0589500	Insex keys n5	PCS	1	5707292137869
3	77557	Slide bar (M8x2)	PCS	2	5707292544995
4	76351	Mounting plate	PCS	2	5707292541741
5	0631990	Button head with flange M6*8	PCS	4	5707292517333
5	1940011	Plastite 60, 4*8	PCS	8	5707292545978
6	1125618	Screenwasher fz6*18*3mm	PCS	4	5707292140494
8	1105080	Stainless washer 8mm	PCS	4	5707292140364
9	1100500	Spring washer 8mm, stainless	PCS	4	5707292140142
10	0661500	Button head screw 8*16 stainl.	PCS	4	5707292139559


Tray plast



Item numbers for Tray plast

Position	No.	Description	Unit of Measure	Quantity	Ean code
1	76350	Table for Stingray, Grey	PCS	1	5707292541727
2	0589500	Insex keys n5	PCS	1	5707292137869
3	77557	Slide bar (M8x2)	PCS	2	5707292544995
4	76351	Mounting plate	PCS	2	5707292541741
5	1940011	Plastite 60, 4*8	PCS	8	5707292545978
8	1105080	Stainless washer 8mm	PCS	4	5707292140364
9	1100500	Spring washer 8mm, stainless	PCS	4	5707292140142
10	0661500	Button head screw 8*16 stainl.	PCS	4	5707292139559
6+7+11	5075364	Brackets for table sz 4	SET	1	5707292543608
9	1100500	Spring washer 8mm, stainless	PCS	4	5707292140142


Floor Frame



Item numbers for Floor Frame

Position	No.	Description	Quantity	Unit of Measure	Ean code
1	821111	Floor Frame	1	PCS	5707292990730
2		Incl.:			
3		Short legs Long legs			

Floor Frame accessories



Item numbers for Floor Frame accessories

Position	No.	Description	Quantity	Unit of Measure	Ean code
2	821001-01	Short leg, 1 pcs	4	PCS	5707292990761
3	82002-02	Long legs, 2 pcs	2	PCS	5707292990754
4	95145-00	Foot plate sz 0 w/o mount rod Model Wombat	1	PCS	5707292995650
	95145-01	Foot plate sz 1 w/o mount rod Model Wombat	1	PCS	5707292995667
5	8910820-0	Rod for footrest, sz 1:110	1	PCS	5707292990747
	8910820-1	Rod for footrest, sz 1,155 mm	1	PCS	5707292457936
6	8910809	Daily adj. f. foot sup	1	PCS	5707292992949