



Fixed-Tilt Cruiser

Standard Specifications

MEASUREMENTS:

INCHES

HCPCS:E1236

- * - Plus or minus .5 inches due to fabric
- Δ - Weight of chair is without front rigging
- - Adjustable push handle is not available on 18 inch models
- ◇ - Models with Textilene Fabric, Shoulder Strap Height is 2 inches shorter
- ^ - Seat back height
- (C) - Cordura upholstery
- (T) - Textilene upholstery

	CX 10	CX 12	CX 14	CX 16	CX18
Seat Width* 	10	12	14	16	18
Seat Depth 	6 – 11	8 – 13	10 – 15	11 – 16	14 – 21
Seat Back Height^ 	21C/20T	22C/20T	25C/23T	27C/24T	25C/25T
Seat to Footplate* 	5 – 21	5 – 25	4 – 22	4 – 22	8 – 23
Seat to Floor* 	23	25	21	23	26
Seat-to-Back Angle 	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°
Seat Angle 	30°	30°	30°	30°	30°
Overall Height 	37	38	37	40	38
Overall Length 	48	44	49	52	55
Overall Width 	22	23.5	25.5	27.5	28.5
Folded Length 	47	48	47	53	50
Folded Width 	15	15	16	16	15
Folded Height 	17	17	17	17	17
Headrest Extension 	6	8	8	8	10
Shoulder Strap Height ◇ 	11 – 19	12 – 20	15 – 23	17 – 25	16 – 22
Push Handle Adjustments ○ 	28 – 43	29 – 46	28 – 45	31 – 48	N/A
Weight Capacity Regular/Transit 	75/66 lbs	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair Δ 	27 lbs	27 lbs	28 lbs	32 lbs	30 lbs



Fixed-Tilt Cruiser

Standard Specifications

MEASUREMENTS:

CENTIMETERS

- * - Plus or minus 1.25 cm due to fabric
- Δ - Weight of chair is without front rigging
- - Adjustable push handle is not available on 18 inch models
- ◇ - Models with Textilene Fabric, Shoulder Strap Height is approximately 5 cm shorter
- ^ - Seat back height
- (C) - Cordura upholstery
- (T) - Textilene upholstery

	CX 10	CX 12	CX 14	CX 16	CX18
Seat Width* 	25.5	30.5	35.5	40.5	45.5
Seat Depth 	15.2 - 27.9	20 - 33	25.4 - 38	27.9 - 40.6	35.6 - 53.3
Seat Back Height^ 	53C /51T	56C/51T	64C/58T	71C/61T	64C/64T
Seat to Footplate* 	13 - 53.4	13 - 63.5	10 - 56	10 - 57	20.3 - 58
Seat to Floor* 	58	63.5	53.3	58.4	66
Seat-to-Back Angle 	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°
Seat Angle 	30°	30°	30°	30°	30°
Overall Height 	94	97	94	102	97
Overall Length 	122	112	124.5	132	140
Overall Width 	56	60	65	70	72.4
Folded Length 	119	122	119	135	127
Folded Width 	38	38	41	41	38
Folded Height 	43	43	43	43	43
Headrest Extension 	15	20	20	20	25.4
Shoulder Strap Height ◇ 	28 - 48	30.5 - 51	38 - 58	43 - 63.5	41 - 56
Push Handle Adjustments ○ 	71 - 109	74 - 117	71 - 114	78 - 122	N/A
Weight Capacity Regular/Transit 	34/30 kg	34/30 kg	45.5/45.5 kg	77/77 kg	144/77 kg
Weight of Chair Δ 	12 kg	12 kg	13 kg	14 kg	14 kg