

Choose the frame and seat combination to find the max. load



Product		Panda Futura Active Size 2	Panda Futura Active Size 2½	Panda Futura Active Size 3	Panda Futura Active Size 4	Panda Futura Active Size 5	Panda Futura Standard Size 1	Panda Futura Standard Size 2	Panda Futura Standard Size 2½	Panda Futura Standard Size 3	Panda Futura Standard Size 4	Panda Futura Standard Size 5	x:panda Size 1S	x:panda Size 1M	x:panda Size 1L	x:panda Size 2S	x:panda Size 2M	x:panda Size 2L	x:panda Size 3S	x:panda Size 3M	x:panda Size 3L	x:panda Size 4S	x:panda Size 4M	x:panda Size 4L	Serval seat Size 1	Serval seat Size 2	Nandu seat Size 1	Nandu seat Size 2	Nandu seat Size 3	Nandu seat Size 4	Stingray seat Size 1	Stingray seat Size 2
Serval Frame Size 1	Ordinary use	30	35	-	-	-	30	30	50	-	-	-	30	30	-	30	-	-	-	-	-	-	-	-	35	-	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	50	-	-	-	30	30	-	30	-	-	-	-	-	-	-	-	35	-	-	-	-	-	-	-
Serval Frame Size 2	Ordinary use	30	50	-	-	-	30	30	50	-	-	-	30	30	30	30	30	-	-	-	-	-	-	-	-	50	-	-	-	-	-	-
	Transportation in vehicles	-	-	-	-	-	30	30	50	-	-	-	30	30	30	30	30	-	-	-	-	-	-	-	-	43	-	-	-	-	-	-