
Pre-sale information (EN 12182)

A:

The latest versions of all instructions are always available on the R82 website R82.org/support and can be printed in larger sizes

B:

The Nandu is a small seat for active children who need little support while seated. The seat is padded and adjustable for a perfect fit. The Nandu is available in four sizes, each with up to 15 cm continuous height adjustment and up to 8 degrees forward tilt. The back rest can be reversed to provide abdominal support when in the forward leaning position. The Nandu is constructed for one user at a time and is suitable for indoor use.

C:

The Nandu is suitable for users needing a simple indoor seat, e.g. children with CP, GMFCS level 2-3.

D:

See measurement chart on R82.org

E:

See measurement chart on R82.org

F:

N/A

G:

See available standard options on R82.org

H:

N/A

I:

See User guide on R82.org/support

J:

N/A

K:

See User guide on R82.org/support

L:

N/A

M:

N/A

N:

N/A

O:

N/A

P:

N/A