



Etac presents:
The Clean family
A cleaner Clean





Discover the Clean family

For almost 20 years, Etac Clean has been a best-selling mobile shower commode chair. There have been copies and imitators, but the original has maintained its position. But we have no intention of resting on our laurels. Our commitment and our collective knowledge mean that we are always striving for continuous improvement – and today, Etac Clean is better than ever.

Etac Clean is a prize-winning range comprising four shower chairs. Each model has unique features which enable a positive, safe bathroom experience for its users. Over the years, the original chair has been improved in small ways that together make a big difference. Clean is now available in a new design with soft, smooth lines which make the chair even easier to clean, and to keep clean. What's more, we are offering Clean in the new colours of lagoon green and white. It's still Clean, but better and more hygienic. A cleaner Clean.



reddot design award

A close-up photograph of the Etac Clean shower commode chair. The image shows a light teal-colored plastic armrest and backrest. The backrest has a grid of oval-shaped perforations. The armrest is ergonomically shaped with a curved top and a flat base. A white vertical support post is visible behind the chair. The overall design is clean and functional.

Etac Clean

shower commode chair

Etac Clean is secure, safe and simple for both users and carers. With practical functionality and design, every detail of the chair has been developed to give the user a positive and safe bathroom experience.

Easy to fit and manouevre

Etac Clean is extremely easy to operate, thanks to the width of the push handle. The clever design can be seen throughout the entire chair. The upright handle ensures that Clean is a perfect fit, even over wall-mounted toilets, and the chair is open at the back to facilitate access.

Comfortable and practical foot support

The foot support is unique and can be pushed under the seat when not in use. The soft curved surface provides good comfort for the feet, which are also protected by the foot support's side protection.

All other variants of Clean also have these features, but they also have additional functionality to meet other needs.



Features Clean

Slim but roomy

The smart design gives Etac Clean a large and comfortable seating area even though the chair is neat and slim.

Easy access

The chair is open at the back to facilitate access for hygiene purposes.

Foot support designed for feet

The practical foot support has a soft, rounded shape which gives the feet extra comfort and protection.

Fits over the toilet

Thanks to the upright operating handle, the chair is a perfect fit, even over wall-mounted toilets.





Etac **Clean Comfort** shower commode chair

All the features of Etac Clean can also be found in this chair, but with even greater comfort. The seat, arm support and back support are padded and guarantee soft and comfortable seating.

Soft and warmth for those who need it

Etac Clean Comfort enhances the feeling of well-being for those who need a little extra comfort and support in the bathroom. The water-repellent and heat-reflective cover of the back support and arm supports helps the user to stay warm. The contoured seat offers comfortable support for the thighs and envelops the user for a safer and more comfortable feeling.

Etac Clean Comfort has the same basic benefits as other members of the Clean family – but now with even greater comfort.



Features Clean Comfort

Comfort and support

The soft enveloping seat and back support provide stability, warmth and softness.

Integrity in exposed situations

Clean Comfort is designed to increase the feeling of integrity and dignity. The comfort padding protects and hides the user without reducing accessibility for the carers.





Etac Clean
Height Adjustable
shower commode chair

In addition to all the features of the basic model, this variant allows you to adjust the height of the seat.

Adjust the height for a perfect fit

The possibility of adjusting the seat height means that Etac Clean Height Adjustable can be tailor-made to suit the user's needs – for example to make transfers easier or to fit perfectly over the toilet. Adjustment is simple and no tools are required. Etac Clean Height Adjustable has the same basic benefits as the other members of the Clean family – but is also height-adjustable up to 60 cm, and has a greater range of foot supports and calf supports.



Features Clean Height Adjustable

The same chair for all heights

The smart design means that the chair is a perfect fit above all wall-mounted toilets, regardless of height.

Additional support for the legs

In the event of additional need for leg and foot support, various accessories are available, e.g. split foot supports and calf supports.



Etac Clean 24" self propelled shower commode chair

This variant of Etac Clean is suitable for those who can operate the chair independently and want to maximise their mobility.



Easy to manoeuvre, for the independent user

The Etac Clean 24-inch wheels are positioned to make it easy to manoeuvre the chair without compromising user safety or comfort. It is easy to transfer to and from the chair thanks to the large, unrestricted seat in front of the wheels. The brake handles are also positioned below the seat so that they do not get in the way during transfer, while remaining easy to reach.

Etac Clean 24" has the same basic benefits as the other members of the Clean family – but allows the user to operate the chair independently.




Features Clean 24"

Self propelled

Easy to operate for the independent user.

Safe and simple transfer

The flat seat and the brake position facilitate lateral transfer.





A cleaner Clean.

Visit Etac.com for more information.

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 – 58 73 30
Fax +46 371 – 58 73 90
info@etac.se www.etac.com


Creating Possibilities